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**English version** on next page

**Users with visual impairments can find the PDF files together with further information on our website:**

**<[www.viteacare.com](http://www.viteacare.com)>.**

**☞ Contact your specialist dealer when required.**

**Alternatively users with visual impairments can have the documentation read out by a helper.**

# Table of contents

<b>Introduction .....</b>	<b>35</b>
<b>Product description.....</b>	<b>36</b>
Indications.....	36
Contraindications .....	36
<b>User safety.....</b>	<b>37</b>
Electromagnetic interference – EMI.....	39
Safety check.....	40
Wheelchair description.....	40
Technical specifications.....	42
<b>Practical advice on use .....</b>	<b>43</b>
Wheelchair preparation.....	43
Before use.....	43
Anti-tip wheels .....	45
Control panel operation.....	46
Joystick - control lever .....	46
Reclining seat and backrest .....	47
Backrest reclining .....	47
Footrest angle adjustment.....	47
Seat lift - lift mechanism.....	47
Getting into the wheelchair .....	48
Wheelchair driving training .....	48
Driving on uneven terrain .....	48
Seat belt .....	49
Fuse.....	50
Drive disengagement clutch .....	50
Getting out of the wheelchair.....	51
Wheelchair transport.....	51
Batteries .....	52
Safe battery charging.....	52
Charging process .....	53
Charging issues .....	54
Batteries – storage.....	54
<b>Maintenance and repair .....</b>	<b>54</b>
Modifications.....	54
Inflating the tyres.....	55
<b>Care and maintenance .....</b>	<b>55</b>
Maintenance.....	55
Cleaning.....	56
Disinfection.....	56
Product life cycle.....	56
Disposal of the electric wheelchair .....	56
Maintenance plan.....	57
<b>Symbols .....</b>	<b>58</b>
<b>Contact details .....</b>	<b>60</b>

## Introduction

DE LUXE LIFT is a wheelchair created with the latest technologies. It has been designed to enable maximum adaptation to the needs of the user. It has a strong, robust structure and modern electronic controls and safety systems. To ensure safety and extend lifespan, the wheelchair should be regularly inspected by a service centre (at least once a year). In case of any doubts or damage, please contact the distributor. Make sure that the enclosed “warranty certificate” has been filled out and submitted to the distributor to ensure the best possible after-sales service.

# ATTENTION

Always observe basic precautions according to the information provided below regarding the safe use of the product. Carefully read the entire instructions before use. Please keep the instructions for future reference.

### WARNING!

Improper use can lead to death, serious injury, unsafe conditions or product malfunction.

### ATTENTION!

Improper use may lead to personal injury and/or damage to the product.

### SUGGESTIONS!

Follow the instructions below to keep your product in good working order.

If the product is used by another person, you must provide them with copies of these instructions for use.

These instructions for use relate to the design and specifications of the product at the time of publication. Due to design changes, some illustrations and photos in the manual may not correspond to the product you have purchased. We reserve the right to make design changes.

The manufacturer disclaims all liability for personal injury or property damage that may result from improper or unsafe use of this product.

## Product description

The electric wheelchair has been created for people with significant motor disabilities. The wheelchair is controlled with a joystick or another control device used to control movement direction. It is a solution for people who cannot use a conventional wheelchair.

The convenient controller ensures excellent manoeuvrability and enables the user to fully manoeuvre the electric wheelchair. Electric wheelchairs can easily revolve around their axis.

## Indications

The product is intended for disabled people with musculoskeletal conditions to increase their mobility.

## Contraindications

Contraindications to sitting.

### SUGGESTION!

To preserve their life, the batteries will be disconnected for the duration of transport. Before using the wheelchair for the first time, make sure to connect the red cable with the input marked with the “+” symbol.

### SUGGESTION!

Avoid full discharge of the batteries. Charge the battery after a full day of use and at least once every 2 weeks.

## User safety

### WARNING!

Every “serious incident” that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State where the user is established.

### ATTENTION!

The product is intended for disabled persons with movement disorders. If you are not sure if your health permits you to use the product, consult a health professional.

### WARNING!

Do not modify the product in any way; this will void the warranty.

### WARNING!

Exceeding the maximum user weight will void the warranty. The manufacturer is not responsible for injury and damage if the maximum user weight is exceeded.

### WARNING!

Before using the wheelchair, read these instructions for use and keep them for future reference.

- a. Read and follow the information included in the instructions for use. It is prohibited to use the product and available optional accessories without carefully reading and understanding these instructions. If the warnings or instructions are incomprehensible, contact the supplier to avoid injuries or accidents. Never try to use the wheelchair in a way that exceeds the limits described in these instructions.
- b. Never switch off the control panel while driving because this may cause the wheelchair to stop suddenly, potentially leading to an injury.
- c. Consult the doctor to check if the drugs you are taking do not affect your ability to drive mechanical vehicles.
- d. To use the wheelchair safely, drive on flat surfaces where all four wheels have adequate contact with the ground.
- e. If the battery level indicator drops to the red level, recharge the wheelchair as soon as possible. If the batteries are fully discharged, they may be damaged, or their lifespan may be shortened.
- f. Keep metallic objects away from battery connections. Electrocution risk.
- g. When driving on curves, always slow down and maintain a stable centre of gravity. When driving on sharp curves, always slow down.
- h. When reversing, make sure that there are no obstacles or objects behind the wheelchair. For safety reasons, reversing speed should be as low as possible.

- i. Children are not allowed to play near the wheelchair or operate it.
- j. Some parts of the electric wheelchair are sensitive to temperature changes.
- k. Remember that the temperature of the components and frame parts may increase when exposed to the sun or another heat source, and at low temperatures, there is a risk that the metallic parts of the wheelchair will be very cold. This may cause burns or frostbite.
- l. Keep the wheelchair in clean and dry conditions.
- m. Do not drive at night without adequate lighting.
- n. Do not remove anti-tip wheels or drive wheels.
- o. Keep your hands away from moving parts when driving. If you wear loose-fitting clothing, make sure that it does not get caught in the drive wheels.
- p. Do not connect the battery charger to an extension cord.
- q. It is prohibited for persons other than authorised service personnel to dismantle the controller, motor or charger – this will void all applicable warranty terms.
- r. Do not use the wheelchair to move on public roads and motorways.
- s. Do not use the wheelchair under the influence of alcohol, medications or other substances that affect driving ability.
- t. Take caution when driving in busy areas and shopping centres.
- u. Never use the wheelchair as a chair for transport in motor vehicles.
- v. Do not try to lift the wheelchair by holding onto parts such as the seat, armrest or parts of the housing.
- w. When driving, do not use a mobile phone or other personal communication devices (e.g., CB radio).
- x. Do not use the wheelchair outdoors during rainfall or snowfall or in high (over 30°C) or low (below 0°C) temperatures.
- y. Do not drive the wheelchair onto port quays, river banks or lake shores that are not protected.
- z. The wheelchair is not suitable for driving in difficult terrain conditions (on unpaved roads).
- aa. Do not stop, turn or reverse suddenly to avoid falls and injuries.
- ab. When driving, keep your forearms on the armrests and your feet on the footrest platform.
- ac. Do not stand up when driving the wheelchair. Do not lean to the side or lift objects from the ground.
- ad. Do not drive onto raised objects (steps, kerbs) higher than 50 mm.
- ae. The wheelchair is not suitable for moving on stairs. Always use ramps or lifts.
- af. The wheelchair is not suitable for passing over gaps in the ground wider than 100 mm.

**WARNING!**

Drive the wheelchair in a straight line up or down the incline, do not reverse. Irregular forward or backward movement combined with changes in inclination may increase the risk of wheelchair tipping.

The maximum user weight is 160 kg

## **Electromagnetic interference – EMI**

Tests of electromagnetic interference (EMI) showed that the interference has an adverse impact on the operation and control of mobile devices with an electric drive.

EMI may come from many sources, such as mobile phones, two-way radios, radio stations, TV stations, amateur radio transmitters (HAM), wireless computer connections, microwave signals, paging transmitters and transmitters used by the vehicles of emergency services. EMI waves may cause unintended movement of the wheelchair or damage the controller. Every mobile device with an electric drive has a specific immunity to electromagnetic interference. The higher the immunity, the better the protection. The intensity of interference can be measured in volts per metre (V/m). The current technology enables the immunity of 20 V/m, which ensures protection against the most common sources of electromagnetic interference. Without any modifications, the wheelchair has an immunity of 20 V/m. As it gets closer to the source, electromagnetic interference becomes stronger. Particular caution should be taken with electromagnetic interference generated by manual radio equipment. A high EMI level may accidentally appear close to the wheelchair control system, potentially affecting its movement and braking. To prevent potential disruptions to the control system, follow the warnings below.

**WARNING!**

Do not use radio transmitters such as a CB radio or a mobile phone near the wheelchair when it is switched on.

## Safety check

### ATTENTION!

Before using the wheelchair, check the following components:

- Electrical connections and wire harnesses – protect from and remove corrosion.
- Batteries – check if they are fully charged and not corroded.
- Brakes – check the functioning of the hand brake.
- Tyre condition – make sure the tyres are not damaged or unevenly worn out – in case of a problem, contact the distributor.

## Wheelchair description


Illustration 1. Side view of the wheelchair.


Illustration 2. Control panel.


Illustration 3. Charger.

## Technical specifications

	DE LUXE LIFT	
Model	PCBL 1620	PCBL 1820
Maximum user weight	160 kg	
Wheel axle distance	450 - 500 mm	
Seat width	430 mm	480 mm
	16 inches	18 inches
Seat hight from ground	560 - 860 mm	
Overall width	640 mm	
Maximum height with headrest	1100 - 1400 mm	
Height of backrest	570 mm	
Overall length	1190 - 1240 mm	
Backrest width	430 mm	480 mm
Permissible tyre pressure	Information shown on the tyre	
Frame material	Aluminum	
Materials	Nylon, polyurethane, rubber	
Battery type	2 x 75AH / 12V	
Motor	450 W	
Energy range	25 km	
Maximum speed	10 km/h	
Turning radius	1400 mm	
Front wheels	230 mm / 9" (pneumatic)	
Rear wheels	355 mm / 14" (pneumatic)	
Charger	DC 24V / 8A	
Functions	Electric backrest and seat tilt adjustment, electric footrest angle adjustment, anti-tip wheels	
Maximum safe inclination	5% (2,86°)	
Description of the operating environment	  	
Storage and operating temperature	From +5°C to +40°C	
Humidity for storage and use	From 30% to 70 %	

Speed and range depends on the user's weight, type of terrain, charge and battery condition. Information is current at the time of publication, the manufacturer reserves the right to change specifications.

## Practical advice on use

### Wheelchair preparation

The wheelchair is supplied in cardboard packaging, and delicate parts are additionally protected against damage during transport. Unpack it carefully to avoid damage.

### Before use

- Check if the batteries are fully charged.
- Check if the driving mode is enabled (drive engagement clutch lever).
- Always plan the route in such a way as to avoid dangerous areas (boggy, rough terrain) and obstacles that cannot be crossed in a wheelchair.
- Check if the seat is correctly set.
- Check if the speed matches the abilities and skills of the user.
- Check if the horn works.

#### ATTENTION!

- When getting into the wheelchair, do not support yourself on the armrests to prevent the wheelchair from tipping and causing an injury.
- Do not stand on the footrests to prevent the wheelchair from tipping (except in the vertical position).
- Remember to check if the wheelchair is off because if the joystick is touched accidentally when you are getting into the wheelchair, the wheelchair may move, causing you to fall and suffer an injury.
- During the transfer, try to be as close to the wheelchair as possible to avoid slipping and falling.

#### WARNING!

The following activities may have a negative impact on the control and stabilisation of the wheelchair:

- Holding onto a pet and taking it out on a leash.
- Carrying passengers (including pets).
- Hanging an object on the controller.
- Towing or being pushed by another vehicle.

Assembly and adjustment of the main parts as indicated below:

- To remove the batteries, loosen the screw shown in illustration 4 and pull out the rail on which the batteries are mounted.
- To install the batteries, perform the operations in reverse order - insert the rail and then lock by tightening the screw shown in illustration 4.

Illustration 4. Battery installation.

- To adjust the height of the armrest, loosen the butterfly screw, place the armrest in the correct position and then tighten the screw (Illustration 5.). Check that the armrest does not move.


Illustration 5. Armrest height adjustment.

- To move the armrests to the side and so increase the available seat width, unscrew the screw with an Allen key (6 mm) located under the seat plate (Illustration 6.).


Illustration 6. Armrest width adjustment.

- To adjust the extension of the control panel, loosen the clamp **A** and adjust the distance of the panel from the armrest. Then tighten the screw firmly and check that the rail with the panel does not move.
- The control panel can be mounted on any armrest (right or left) and also at the back in a dedicated mount located behind the backrest. It can then be used by the carer.


Illustration 7. Control panel adjustment.

- The control panel can also be tilted outwards as shown in illustration 8.

Illustration 8. Tilting the control panel control.

## Anti-tip wheels

To unfold the anti-tip wheels, tilt them as shown in illustration 9.

Illustration 9. Unfolding anti-tip wheels.

## Control panel operation

The electronic control panel is used to drive the electric wheelchair. It supplies all operating systems. It can be used to control the movement of the wheelchair and monitor the battery level and the condition of the entire electrical system. A typical control panel is placed on one of the armrests and is connected to the battery and the motor.

- To start the electrical system of the wheelchair, press the power button (Illustration 3.).
- To change wheelchair speed, press the left button to reduce speed or press the right button to increase the maximum speed (Illustration 3.).
- To change the driving mode to seat and backrest angle control, press the function selection button (Illustration 3.).

## Joystick - control lever

The joystick can be used to control the direction and speed of the wheelchair. When you move the joystick from the neutral position (central setting), the electromagnetic brake of the wheelchair will be released, enabling its movement. The more you move the joystick away from the neutral position, the faster the wheelchair will go. Returning the joystick to the neutral position activates the brake system, slowing the wheelchair down until it stops completely. Moving the joystick from the neutral position backwards causes the wheelchair to move backwards. Moving the joystick to the side causes the wheels to turn left or right.

### WARNING!

The wheelchair has electromagnetic brakes. Take caution when driving because if the joystick suddenly returns to the neutral position, the wheelchair will stop as quickly as possible.

### WARNING!

Take special care that no parts of your clothing or body are in danger of being caught when lowering the seat or reclining the backrest. Special care must be taken.

### WARNING!

When raising the wheelchair seat to its maximum, it is recommended not to use (not to apply) the permitted maximum speed. This may result in a loss of stability and a fall.

## Reclining seat and backrest

- The wheelchair has a tilt function for the seat including the backrest.
- To tilt the seat and back at the same time, press the function selector button (illustration 10.)
- A sideways movement of the joystick highlights the respective elements on the wheelchair icon that you wish to control.
- When the backrest and seat are illuminated, it is possible to tilt these elements by moving the joystick forwards or backwards.


Illustration 10. Seat inclination.

## Backrest reclining

- To tilt the backrest, select the backrest highlight on the wheelchair icon after pressing the function select button.
- By moving the joystick forwards or backwards, adjust the position of the backrest as desired.


Illustration 11. Backrest reclining.

## Footrest angle adjustment

- To tilt the footrest, after pressing the function selection button and tilting the joystick sideways, select the highlight of the required footrest on the wheelchair icon.
- Set the desired angle of the footrest by moving the joystick forwards or backwards..

Illustration 12. Footrest tilt adjustment.

## Seat lift - lift mechanism

- To raise the seat with backrest, after pressing the function selector button and tilting the joystick to the side, select the options to highlight all elements on the wheelchair icon.
- Use the joystick to move the joystick forwards or backwards to set the desired height..

Illustration 13. Adjustable seat lift.

## Getting into the wheelchair

### WARNING!

Getting in or out of the electric wheelchair can be performed only when the wheelchair is off and when the drive disengagement clutch lever is set in the driving mode (lever in the bottom position). Unintended movement of the joystick (driving and turning lever) can result in the uncontrolled movement of the electric wheelchair – Accident risk!

## Wheelchair driving training

1. Make sure to learn how to handle the wheelchair and perform all the required manoeuvres to move safely. First, learn how to move forward. Make sure to set the lowest possible speed.
2. Once you learn how to move forward, try to trace an “S” on the floor in the wheelchair.
3. Then, you can train driving backwards. Regardless of the programmed forward speed, the wheelchair will always move slower when reversing.

## Driving on uneven terrain

Moving on uneven terrain presents more risks. It is important to observe the following safety rules. Failure to observe them might lead to severe injuries.

- Do not drive onto inclines greater than 100. Remember that unpaved or wet surfaces further reduce traction and brake power when driving up or down the incline.
- In the vertical position, do not drive up inclines greater than 30. Speed is automatically reduced adequately to inclination and position. In extreme conditions, the system will prevent further movement.
- Do not reverse down inclines.
- Do not drive over uneven, unpaved surfaces (sand, gravel, grass).
- When driving on a decline, move at minimal speed.
- Do not get off the wheelchair or get into the wheelchair when it is parked on an inclined surface (always place the wheelchair on an even, horizontal surface first).
- When driving on inclines, always drive in a straight line along the slope.

## Seat belt

### WARNING!

Remember to always fasten the seat belt before driving.

Make sure that no objects are compressed by the belt. This way, you will avoid painful pinch points. Subsequent installation of a lap belt can be performed only by a specialised service centre!

The lap belt is not a part of the fastening system of the electric wheelchair and/or its user during transport in a mechanical vehicle.

The lap belt is used to stabilise the users in the sitting position and prevent them from falling forward out of the electric wheelchair.

To fasten the lap belt, pull both of its ends forward and close the latch so that you hear a click. To open the lap belt, press the button and slide the ends of the belt apart.

Illustration 14. Correctly fastened lap belt.

### ATTENTION!

The length of the lap belt can be adjusted but remember not to fasten it too tight.

## Fuse

When the wheelchair's power system overheats, press the fuse to break the circuit and protect it from damage. The fuse is located on the battery housing at the rear of the wheelchair (Illustration 15.).

Illustration 15. Fuse location.

### WARNING!

Do not touch the motor to avoid the risk of burns.

## Drive disengagement clutch

On the wheelchair, next to the motors on each side, there are clutch levers (yellow switch colour) to engage or disengage the electric drive and push the wheelchair manually.

- To disengage the drive for manual pushing, turn the yellow switch to the horizontal position relative to the wheelchair frame.
- Re-engage the drive (drive mode) by turning the yellow switch located on the motor to the vertical position.

### WARNING!

Do not disengage the electric drive if there is no one to push the wheelchair.  
Do not disengage the electric drive if the wheelchair is located on a slope because it might roll down, potentially leading to injury.

**Remember that when the electric drive is disengaged, the electromagnetic brake system is disengaged as well.**

Always stand behind the wheelchair when disengaging the electric drive. Do not do this while in the wheelchair. After moving the wheelchair to the desired location, remember to re-engage the electric drive.

## Getting out of the wheelchair

### WARNING!

Never get into or out of the wheelchair without switching it off first. This prevents the wheelchair from moving if the joystick is touched by accident.

1. Stop the wheelchair completely.
2. Switch the wheelchair off using the control panel.
3. Dismantle the armrest on the side on which you can safely leave the wheelchair.
4. Carefully and safely get off the seat and sit in a safe position next to the wheelchair.
5. To make it easier to get into the wheelchair, you can leave the armrest dismantled.

## Wheelchair transport

The wheelchair should only be transported using suitable vehicles. During transport, it should be fastened with straps fastened to the metal chassis frame. The wheels should be locked. It is prohibited to stay in the wheelchair during transport because this presents a risk of falling and severe injuries. The wheelchair cannot be transported in unsuitable vehicles – without a special fastening system – because this may damage the wheelchair. Always make sure that the wheelchair is appropriately fastened in the vehicle; otherwise, it may cause injuries to the persons in the vehicle and severely damage the vehicle itself.

When moving up a ramp into the vehicle, the wheelchair should be secured from behind to prevent it from tipping.

All additional components of the wheelchair should be removed from the wheelchair during transport, or they should be attached to the wheelchair in such a way that they do not harm people or damage the transport vehicle when it brakes suddenly.

## Batteries

1. The wheelchair has long-life batteries. The batteries are sealed and do not require maintenance (so-called maintenance-free batteries).
2. The following situations may increase energy consumption and reduce the range of the wheelchair: driving up slopes, overloads, winter.
3. Complying with the following rules will extend the life of the batteries:
  - Avoid full discharge of the batteries. Charge the battery after a full day of use.
  - During the first charging, keep the battery plugged in for 40 hours until it is fully charged.
  - If you do not use the wheelchair regularly, charge the batteries every three days using a full charging cycle to preserve their life.
4. Do not remove the cover from maintenance-free batteries. Do not add water to the battery.
5. Leaks from an acid battery may cause burns or eyesight loss. If the acid from the battery comes into contact with the skin, eyes, clothing or any other object, rinse it immediately with plenty of running water and immediately report for a medical consultation.
6. Do not use batteries other than the original batteries dedicated to the wheelchair to avoid damage to the wheelchair.
7. Make sure that both the wheelchair and the batteries are correctly protected during transport.
8. Batteries wear out, which is why the range of the wheelchair decreases over time. Replace the batteries when the range of the wheelchair decreases by half even when the battery is fully charged.
9. Replace both batteries at the same time to ensure that they can function effectively for a long time.

## Safe battery charging

- Use a standard earthed 3-pin plug and an earthed outlet.
- Do not use an extension cord or adapters to avoid the risk of fire, electric shock or injuries.
- Do not use the same outlet for activities other than charging the battery.
- Only use the original charger dedicated to the wheelchair to avoid damage.
- Do not leave the charger in the port if the batteries are not charged.
- Keep the charger and the batteries away from water and moisture, do not handle the charger with wet hands.
- Charge the wheelchair in a well-ventilated area.
- Do not cover the fan of the charger to enable heat dissipation.
- Do not use the equipment until the batteries are fully charged. Do not sit in the wheelchair while the batteries are charging.
- Do not stop battery charging until the charge level indicator shows that the batteries are fully charged.
- Do not place the charger directly on the battery.

- Check if the charger specifications match the parameters of the batteries and the mains voltage. Connecting a battery with different specifications may cause it to be damaged and result in severe injuries to the user.
- In particular, do not expose the wheelchair and charger to high temperatures or smoke during charging to avoid ignition and explosion.
- To avoid electrocution, do not touch the charger or the battery with wet hands or other objects and do not expose them to moisture.
- Check if the switch of the charger is placed in the “O” position before the charger is connected to the mains and before the batteries are connected to the charger.
- Do not expose the charger to rain, water and moisture.
- Charging should be performed at ambient temperatures from 00C to 400C; charging in different conditions may be unsuccessful or damage the batteries.
- If the charging process indicator is not lit up, disconnect the cables for 10 minutes and re-connect them again. If this does not help, contact the service centre.
- Contact the service centre if the charge level indicator remains yellow after charging for more than 12 hours.

## Charging process

1. Park the wheelchair as close to an earthed outlet as possible.
2. Switch the wheelchair off.
3. Make sure that the electric power supply of the wheelchair is switched on (the clutch is in the correct position – driving mode).
4. Connect the charger to the charging port.
5. Connect the charger to an earthed outlet, remember to check if the outlet has the right voltage, i.e., 115 V / 230 V.
6. Switch on the charger (red – low battery level, yellow – medium battery level, green – fully charged batteries).
7. If the LEDs are blinking in an irregular fashion, switch the charger off. Check the connection of the charger and switch it back on.
8. If the charging procedure is correct, the LEDs indicating the battery level will light up (as described above).
9. The optimum charging time is 12 h (it is slightly longer in winter).
10. When charging ends, disconnect the charger from the outlet first and then remove the plug from the charging port.

## Charging issues

1. No indication of voltage after the charger is connected to the power supply.
  - Check the connection of the cables.
  - Check if the charger port is working.
  - Check if the fan is working.
2. Red and yellow charging level indicators are lit, but the fan does not work despite the fact that the power is on.
  - Check if the charger is correctly connected to the batteries.
3. The red LED of the charging level is lit up steady, and the green LED is blinking, there is no charging signal despite the fact that the power is on.
  - Check the reverse polarity fuse.
  - Check the connection to the batteries.
  - Check the polarity of the batteries.
  - Check if the batteries are placed in the correct position.

## Batteries – storage

- If the wheelchair is not going to be used for a longer period, remove the batteries and keep them in a dry, cool place.
- Avoid full discharge of the batteries.
- If you use the wheelchair every day, charge the batteries every day.
- If you use the wheelchair once a week, charge the batteries once a week.
- Keep the batteries fully charged.
- If you do not use the wheelchair regularly, charge the batteries at least every two weeks using a full charging cycle to preserve their life.
- Never replace the batteries with batteries other than the original batteries supplied with the wheelchair.
- If the batteries are faulty or have to be replaced, contact the distributor.
- Do not discard the batteries into the waste bin – hand them over to the distributor.
- To prevent short-circuits and damage to the battery, never touch their contacts with metallic objects.

## Maintenance and repair

Perform the inspection regularly, it should be preceded by a thorough cleaning of the product (to identify any hidden damage). If you need to perform repairs or have questions concerning the product, contact the distributor.

## Modifications

Do not modify, remove, switch off or add any parts, components or functions of the wheelchair.

## Inflating the tyres

- Unscrew and remove the plastic cap from the valve and put it in a safe place. Do not lose it because the valve will get dirty, making the tyres harder to inflate and accelerating air loss.
- Put the nozzle of the pump against the valve and check if it adheres correctly. If necessary, press the pump against the valve to prevent the air from escaping to the sides.
- If the pump has a lever with a clamp on the nozzle, clamp it.
- If the nozzle of the pump and the valve have a threaded end, screw them together. This way, the air will be injected directly into the tube, and the nozzle will not detach itself from the valve.
- When inflating, keep checking the pressure in the tyres every now and then.
- After inflating the tyres, remove or screw out the nozzle from the valve while pressing your thumb against the nozzle outlet. Make sure you do this skilfully enough to prevent too much air from escaping.
- Screw the plastic cap onto the valve.

Illustration 16. "Schrader" valve  
(automotive)

### ATTENTION!

Inflate the wheels according to the correct tyre pressure specified by the manufacturer on the sides of the tyre.

Use pumps with a pressure gauge or a manual bike pump.

## Care and maintenance

### ATTENTION!

Insufficient or negligible care and maintenance of the electric wheelchair limit the liability of the manufacturer.

## Maintenance

To ensure long and efficient operation, the wheelchair should be subject to basic maintenance activities and inspections carried out by the service centre. The inspection should be carried out once a year. In case of any doubts or damage, contact a service centre that has all the required qualifications to solve the issue.

## Cleaning

### WARNING!

Follow the safety instructions provided by the manufacturer of the cleaning and disinfecting products. Failure to observe this requirement may cause skin irritation or early damage to the upholstery or the electric wheelchair.

### Do not use:

- solvents,
- toilet cleaning products,
- chlorine-based products,
- sharp brushes and hard objects,
- corrosive detergents.

*Do not clean the wheelchair using pressure or steam cleaners!*

## Disinfection

If the product is used by many persons (e.g., at a nursing home), use commercial disinfectants. Before disinfecting the wheelchair, clean the outer surface. Only use proven and approved disinfectants.

To obtain information about proven and approved disinfectants and disinfection methods, contact the local sanitary competent authority.

### WARNING!

Disinfectants may damage the outer surface, reducing the long-term functionality of the components. Always observe the instructions concerning the use of the products specified by the manufacturer.

## Product life cycle

The expected life of the medical device is 4 years if it is used as intended and if all maintenance and repair activities are performed as required. The life cycle of the medical device also depends on the frequency of use, the operating environment and care of the device. The life cycle of the device can be extended using spare parts. Spare parts are normally available for 5 years after the model is withdrawn from production. The expected life cycle of the medical device is not equivalent to an additional guarantee.

## Disposal of the electric wheelchair

The electric wheelchair has to be disposed of according to the applicable local and national codes. Information about the correct disposal of the packaging, metallic parts of the frame, plastic parts, electronic systems, batteries, neoprene, silicone and polyurethane materials can be obtained from the local company responsible for waste disposal or from the authorised supplier.

## Maintenance plan

When	What	Note
<b>Before driving</b>	Perform a general functional check.	Perform the check yourself or with an assistant.
	<b>Inspection of the magnetic brake.</b> Check the functioning of the lever/clutch switching the driving mode.	Perform the check yourself or with an assistant. If the wheelchair can be pushed, immediately request brake repair from the service centre – accident hazard!
Especially <b>before driving after dusk</b>	<b>Lighting.</b> Check the lighting system and headlight to verify if they function correctly.	Perform the check yourself or with an assistant.
<b>Every 2 weeks</b> (depending on the travelled distance)	<b>Tyre pressure check.</b> Tyre pressure: 3.5 bar = 50 psi	Perform the check yourself or with an assistant. Use a pressure gauge to do this.
	<b>Screws.</b> Check the tightness of bolts and nuts.	Perform the check yourself or with an assistant. Securely tighten any loose screws. If necessary, contact the service centre.
<b>Every 6–8 weeks</b> (depending on the travelled distance)	<b>Affixing the screws.</b> Check the nuts and screws to determine if they are tightened correctly.	Perform these activities yourself or with another person. Securely tighten any loose wheel nuts/screws and re-tighten them after 10 hours of operation or 50 km. Where necessary, contact the service centre.

<b>Every 2 months</b> (depending on the travelled distance)	<b>Tyre profile check</b> The minimum tyre tread depth is 1 mm	Perform the check yourself or with another person. If the profile is worn out or the tyres are damaged, request a professional service centre to perform the repair.
<b>Every 6 months</b> (depending on the frequency of use)	<b>Check</b> - Cleanliness - General condition	Perform this action yourself or with another person.
Recommendations of the manufacturer: <b>Every 12 months</b> (depending on the frequency of use)	<b>Inspection work</b> - Vehicle - Charger	To be performed by a professional service centre.

## Symbols

The symbols used below on the electric wheelchair are used to identify warnings, prohibitions and mandatory actions. Read them all and understand their meaning.

	Medical device
	Date of manufacture
	Manufacturer
	Distributor
	Conformité Européene (compliance with the European Community)
	Catalog number (to re-order)
	Serial number
	Lot code

	Unique Device Identifier
	Refer to the instructions for use
	Follow the instructions for use
	Warning
	Protect from humidity
	Possibility of jamming fingers
	Permitted temperature
	For indoor/outdoor use
	Maximum user weight
	Maximum load capacity
	Maximum safe inclination
	The product is designed to be transported by car
	The product is not intended to be used as a seat for transport in a motor vehicle
	Indicates attachment points on the scooter

#### SEPARATE COLLECTION FOR WASTE OF ELECTRICAL AND ELECTRONIC EQUIPMENT.

Do not dispose of battery in municipal waste. The symbol indicates separate collection for battery is required. Such products should be handed over to collection sites for this type of waste for proper disposal, renovation or recycling. Alternatively, in some countries of the European Union or other European countries you may return your products to your local retailer when purchasing a similar new product. By correct disposal of this product, you will help to conserve valuable natural resources and help prevent potential negative effects on the environment and human health, which could be caused by incorrect disposal of waste. Further information is available from your local authority or waste collection point. Penalties may apply in accordance with local regulations in the event of improper disposal of this product. For entities in the European Union countries. If you want to dispose of an electrical or electronic device, ask your dealer or supplier for the necessary information.

Liquidation in countries outside the European Union: This symbol is valid in the European Union. If you want to dispose of this device, obtain the necessary information on proper disposal at local authorities or from your dealer.

## Contact details

### Manufacturer

mdh Sp. z o.o.

Address: ul. Maratońska 104, 94-007 Łódź, Poland

tel. +48 42 674 83 84, fax. +48 42 636 52 21

[www.mdh.pl](http://www.mdh.pl)

[www.viteacare.com](http://www.viteacare.com)

We reserve the right to make technical changes due to technical progress.