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**English version** on next page

**Users with visual impairments can find the PDF files together with further information on our website:**

**<[www.viteacare.com](http://www.viteacare.com)>.**

**☞ Contact your specialist dealer when required.**

**Alternatively users with visual impairments can have the documentation read out by a helper.**

# Table of contents

<b>Introduction from mdh Sp. z o.o.....</b>	<b>42</b>
Indications.....	43
Contraindications .....	43
<b>Safety instructions .....</b>	<b>43</b>
<b>Safety instructions .....</b>	<b>45</b>
General safety rules .....	45
Modifications .....	45
Security check.....	45
Capacity .....	45
Transport the scooter .....	46
Inclination guidance.....	46
Driving outdoors and precautions in bad weather conditions .....	46
Physical limitations .....	47
Information on electromagnetic interference - EMI .....	47
<b>Product information.....</b>	<b>48</b>
Clutch lever for disengaging the drive .....	49
Battery charge indicator.....	49
Drive control lever .....	50
Power switch .....	50
Speed control dial .....	50
Horn.....	50
<b>Operate.....</b>	<b>50</b>
Before you start using .....	50
Getting on the scooter.....	51
Pre-ride inspection and adjustment.....	51
Operating a scooter .....	51
Getting off the scooter .....	52
Power off function after a certain time.....	52
<b>Regulation and adjustment of scooter components .....</b>	<b>53</b>
Delta steering column and rear view mirrors adjustment (Option).....	53
Steering column angle adjustment.....	53
Seat rotation adjustment.....	54
Armrest spacing adjustment .....	54
Seat height adjustment.....	55
<b>Disassembly and assembly the scooter .....</b>	<b>55</b>
Disassembly .....	55
Disassembly procedures .....	55
Disassembly the frame (MINI II - D only) .....	56
Assembly .....	56

<b>Basic troubleshooting .....</b>	<b>57</b>
<b>Batteries and charging .....</b>	<b>58</b>
First charge.....	58
Charging the scooter.....	59
Operating temperatures .....	59
Battery operation and charging.....	59
Charging the batteries.....	59
<b>Care and maintenance .....</b>	<b>61</b>
Tires.....	61
Wheel Replacement.....	61
External surfaces .....	61
Cleaning and disinfection.....	61
Battery terminal connections.....	62
Electric wiring harnesses .....	62
Daily checks .....	62
Weekly checks.....	62
Monthly checks .....	62
Yearly checks .....	62
ABS plastic covers.....	63
Axle bearings and motor/transaxle .....	63
Console, Charger, and Rear Electronics .....	63
Nylon Lock Nut Replacement.....	63
Storage of the scooter .....	63
Disposal of an electric scooter .....	64
<b>Technical parameters .....</b>	<b>64</b>
<b>Warranty.....</b>	<b>65</b>
Two (2) years limited warranty.....	65
One (1) year limited warranty .....	65
Warranty does not include.....	66
Batteries.....	66
Inspection and warranty service.....	66
Replacement parts.....	66
<b>Product safety symbols.....</b>	<b>67</b>
Symbols.....	67
<b>Contact details .....</b>	<b>72</b>

## Introduction from mdh Sp. z o.o.

This scooter has been manufactured by skilled and dedicated personnel. It has been designed and manufactured to the highest quality standards. Before using your scooter for the first time, please read these operating instructions. Improper use of the scooter may result in personal injury or traffic accidents. If you still have questions after reading this manual, please contact your specialist dealer.

# ATTENTION

Always follow the basic precautions according to the information given below regarding safety when using the product. Read the entire manual carefully before use. Keep this manual for future reference.

### WARNING!

Improper use can lead to death, serious injury, unsafe conditions or product malfunction.

### ATTENTION!

Improper use may lead to personal injury and/or damage to the scooter.

### SUGGESTION!

Follow the instructions below to keep your scooter in good working order.

If the electric scooter is used by another person, you must provide them with copies of these instructions for use.

These instructions for use relate to the design and specifications of the product at the time of publication. Due to design changes, some illustrations and photos in the manual may not correspond to the scooter you have purchased. We reserve the right to make design changes. The product shipped in its original packaging requires partial assembly - follow the information in the "Disassembly and assembly the scooter" section.

The manufacturer disclaims all liability for personal injury or property damage that may result from improper or unsafe use of this product.

## **Indications**

The product is intended for disabled people with musculoskeletal conditions to increase their mobility.

## **Contraindications**

Contraindications to upright standing. Injury or paresis of the upper limbs preventing the use of a power wheelchair.

## **Safety instructions**

- Please note that the device is a small electric scooter and is intended mainly for indoor use. Special care should be taken caution when driving around bends, obstacles and inclines.
- Always ensure that the electric scooter is fully unfolded before getting on. Do not fold the scooter until you have fully dismounted the scooter, folded the scooter seat back and lowered the armrests to their lowest position.
- The folding mechanism of the electric scooter is equipped with an anti-trapping function. This is primarily to prevent damage to the armrests and backrests. Folding will be locked when the user is seated in the scooter. Assistance may be required to unfold the scooter. Folding the scooter while you are in it is strictly prohibited.
- Do not operate the folding mechanism while the user is on the seat.
- Read and follow the information in the instructions for use. It is forbidden to use this product and the available optional equipment without first carefully reading and understanding this manual. If you do not understand the warnings, cautions or instructions, contact your supplier, otherwise personal injury or accident may result. Never attempt to use the electric scooter beyond its limitations as described in this manual.
- Place your feet on the footrest when driving. Do not stand on the footrest.
- To operate an electric scooter safely, it should be ridden on flat surfaces where all four wheels have sufficient contact with the ground.

- If the battery level indicator drops to red, recharge the scooter as soon as possible. Completely discharging the batteries may damage them or shorten their life.
- Keep metal objects away from the battery connections. Risk of electric shock.
- When cornering, always reduce speed and keep the center of gravity stable. Always reduce speed when negotiating sharp bends.
- When reversing your electric scooter, make sure there are no obstacles behind you.
- It is prohibited for children to play near or operate the scooter.
- Some parts of the electric scooter are susceptible to temperature changes. The controller can only operate at temperatures between -10°C ~ 50°C (-14°F to 122°F).
- In extremely cold temperatures, the batteries may freeze and the electric scooter may not operate. In extremely high temperatures, the scooter may run at a lower speed due to the controller's heat protection function, which prevents damage to the motors and other electrical components.
- Store in clean and dry conditions.
- It is forbidden to drive at night without proper lighting.
- Do not remove the anti-tip wheels or transport wheels.
- Keep hands and feet away from moving parts when driving. Take care that loose clothing does not get caught in the drive wheels.
- Do not connect the battery charger to an extension cord.
- Disassembly of the controller, motor or charger by anyone other than an authorized service is prohibited and will invalidate any applicable warranty conditions.
- It is prohibited to drive the scooter on public roads and motorways.
- It is forbidden to use an electric scooter under the influence of alcohol, medication or any other substance affecting the ability to drive.
- Be careful when driving in busy areas and shopping centers.
- Under no circumstances should the scooter be used as a seat when being transported by motor vehicles.
- Do not attempt to lift the scooter by gripping anything other than the frame (for example, the seat, armrest or chassis components).
- Please note that the temperature of the components and frame parts may increase when exposed to the sun or other heat source, while at low temperatures there is a danger of very high cooling of the metal parts of the wheelchair. This can cause burns or frostbite on parts of the body.

**WARNING!**

The product is intended for use by people with disabilities and movement disorders. If you are not sure whether your health condition allows you to use the product, consult a healthcare professional.

**WARNING!**

Any "serious incident" related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides.

# Safety instructions

## General safety rules

Read the instructions for use carefully with understanding before using the electric scooter.

Users may encounter difficulties with maneuvering when passing through narrow doorways, ascending and descending ramps, driving around corners and navigating uneven terrain. Reduce speed, take your time and maneuver the electric scooter carefully.

## Modifications

Do not modify, remove, disable or add any parts, components or features to your electric scooter. The manufacturer offers a wide range of electric scooter accessories to meet your needs.

### WARNING!

Do not modify the electric scooter in any way; this will void the warranty.

## Security check

ATTENTION - Check the following before use:

- Electrical connections and wiring harnesses - protect and remove corrosion.
- Batteries - fully charged and not corroded.
- Brakes - check that the electromagnetic brake is functioning.
- Folding mechanism - make sure the scooter is fully unfolded.
- Condition of tires- make sure they are not damaged or excessively worn  
- if there is a problem, contact your dealer.

## Capacity

This electric scooter is designed and built to transport one person with a maximum weight of 136 kg. The scooter is not designed to carry goods or objects other than those previously described.

### Warning!

Exceeding the maximum user weight will void the warranty. The manufacturer is not responsible for injury or damage caused by failure to comply with the maximum user weight.

## Transport the scooter

When transporting the electric scooter by vehicle, store it safely in the back of a van, truck or in the boot of a car. Moving parts should be removed or properly secured for transport.

Do not sit in the vehicle when is moving.

## Inclination guidance

Remain cautious when going up and down hills. If necessary, reduce your speed when going downhill. If you are moving up a hill try to keep the scooter moving forward. If necessary, stop. Be careful and slowly accelerate your scooter forward. If you are descending, reduce the speed of the electric scooter to the lowest setting and then proceed with caution. If the scooter is descending faster than you anticipated, slowly release the drive control lever until it stops and then push the lever slightly to continue descending safely.

### WARNING!

Drive the electric scooter forward down or up a slope, never reverse. Irregular movements or making back and forth movements may increase the probability of tipping over. Changes in slope can reduce stability. Under no circumstances should you travel up a slope with a gradient greater than the maximum stability angle of this electric scooter.

## Driving outdoors and precautions in bad weather conditions

Although the electric scooter is designed to be used indoors to help with daily living activities, the manufacturer has designed it to perform perfectly outdoors on level surfaces.

### ATTENTION!

When driving outdoors, avoid uneven terrain, soft surfaces, tall grass, loose gravel, loose sand, fresh and salt water, stream banks, lakes and oceans. If you are unsure about a surface, avoid it.

Exposure to adverse weather conditions should always be avoided. If you find yourself in bad weather conditions while using your electric scooter, immediately go to the nearest shelter. Dry the electric scooter completely before use, before charging or before storing.

## Physical limitations

Use awareness, caution, prudence and common sense when operating the scooter. Always be aware of your own limitations and rules of use when operating an electric scooter.

### WARNING!

It is forbidden to use the scooter while under the influence of alcohol or psychotropic drugs. Never smoke while sitting in the scooter. Damaged or worn upholstery must be replaced as it increases the risk of fire.

### ATTENTION!

When operating an electric scooter, please be aware of precautions, warnings and safety issues if you are taking medications prescribed by a doctor or available without a prescription.

## Information on electromagnetic interference - EMI

Studies of electromagnetic interference (EMI) have shown that it adversely affects the operation and control of electrically powered mobile devices.

EMI can come from a variety of sources, such as cell phones, two-way radios, radio stations, TV stations, amateur radio (HAM) transmitters, wireless computer links, microwave signals, paging transmitters and transmitters used by emergency vehicles. EMI waves can cause unintended movement of the scooter or damage to the controller. Every electrically powered mobile device has a certain immunity to electromagnetic interference. The higher the immunity level, the greater the protection. The intensity of interference can be measured in volts per meter (V/m). Current technology allows for an immunity level of 20 V/m, which provides protection against the most common sources of electromagnetic interference. Your scooter, without any modifications, has an immunity level of 20 V/m. As you get closer to the source, electromagnetic interference becomes more intense. Special attention should be paid to electromagnetic interference from hand-held radio devices. It is possible to unintentionally bring high levels of EMI close to the scooter's control system, which can affect movement and braking. To prevent possible interference with the scooter's control system, follow the warnings below.

### WARNING!

Do not use hand-held radio transmitters such as a CB-radio or cell phone near a powered-up scooter.

**WARNING!**

Drive the electric scooter in a straight line up or down a slope, never reverse. Irregular forward and backward movements with changes in slope can increase the risk of scooter rollover.

**WARNING!**

Pay attention and stay away from nearby transmitters, such as radio or TV stations.

## Product information

The MINI II scooter is easy to use and designed primarily for indoor use. It is foldable, making it very compact for storage and transportation. The following illustration shows all parts relevant to the user/driver. The section on operation and maintenance describes the parts used. The detachable MINI II-D model has the drive disconnect clutch lever on top of the rear cover. The MINI II-F model is one-piece, and the clutch lever that disconnects the drive is on the rear of the wheelchair.

Figure 1. Important parts

- |                      |                                |
|----------------------|--------------------------------|
| 1. Operating console | 6. Footplate                   |
| 2. Handlebars        | 7. Seat with backrest          |
| 3. Steering column   | 8. Lever disengaging the drive |
| 4. Battery pack      | 9. Rear wheel with drive       |
| 5. Front wheel       |                                |

## Clutch lever for disengaging the drive

Whenever you need or want to push the scooter for a short distance, you can switch the drive disconnect clutch lever.

Locate the lever at the top of the rear frame of the wheelchair, on the right side (for MINI II-D model only).

Push the clutch lever forward to disconnect the drive and brake. Now you can push the scooter.

Pull the lever toward the rear of the scooter to engage the drive and brake system. Push mode disconnected.

### WARNING!

Remove the key from the power switch before disconnecting or engaging the drive in the scooter. Never sit in a scooter that is on an incline while switching the drive clutch lever. If the scooter has the drive disconnected then the brakes are also not active.

- Disconnect the motor drive only on a level ground.
- Make sure the key is removed from the power switch.
- Stand next to the scooter when engaging or disengaging the drive clutch lever.
- Never sit in the scooter while engaging or disengaging the drive.
- When you have finished pushing, engage the drive to lock the brakes.


Figure 2. Control panel

- |                             |                       |
|-----------------------------|-----------------------|
| 1. Battery charge indicator | 4. Speed control dial |
| 2. Drive control lever      | 5. Horn               |
| 3. Power switch             |                       |

## Battery charge indicator

When the key is fully inserted into the power switch, the charge indicator shows the power level of the batteries. For more information about the battery charge indicator, look at chapter "Batteries and Charging."

## Drive control lever

The lever allows you to control forward and while reversing, to get:

- Maximum speed set with the speed control dial.
- Place your right hand on the right handle and your left hand on the left handle.
- Use your right thumb while pressing the right side of the lever to disengage the scooter's brakes and move forward.
- Release the lever and allow the scooter to brake completely before pressing the left side of the lever allowing you to move backward.
- When the lever is fully released, it automatically returns to the center "stop" position, the brakes are automatically applied.

## Power switch

- Insert the key completely into the power switch to start the scooter.
- Remove the key from the power switch to turn off the scooter.

When you have finished riding / using the product, remove the key to prevent unauthorised use of the scooter.

## Speed control dial

The dial allows you to set the maximum speed limit.

- The bird symbol represents the slowest setting.
- The airplane symbol represents achieving maximum speed.

## Horn

This button activates the warning signal. The scooter must be powered on for the horn to be active. Do not hesitate to use the horn if it can prevent an accident or injury.

## Operate

### Before you start using

- Have you fully charged the batteries? See the section "Batteries and charging"
- Have the manual drive disconnect clutch lever in the rear position.
- Never leave the lever in the forward position if you do not intend to push the wheel-chair manually.

## Getting on the scooter

### WARNING!

Never attempt to get on or off the scooter without first removing the key from the power switch. This prevents the scooter from being moved by accidentally switching the drive control lever.

1. Make sure the key is out of the power switch.
2. Stand next to the scooter.
3. Unlock the seat rotation lever and rotate the seat so that it faces you.
4. Make sure the seat is locked in this position.
5. Assume a safe and comfortable position on the seat.
6. Unlock the seat rotation lever and position the seat facing forward.
7. Make sure the seat is locked in this position.
8. Make sure your legs are securely located on the footplate.

## Pre-ride inspection and adjustment

- Is the seat at the right height? See the section “Regulation and adjustment of scooter components”.
- Is the seat securely locked?
- Is the steering column angle properly adjusted and locked? See the section “Regulating and adjusting scooter components”.
- Is the key inserted into the power switch?
- Is the horn work correctly?
- Are there no people, animals or other obstacles along the route?
- Have you planned your route to avoid going through difficult terrain and, if possible, as few inclines as possible?

## Operating a scooter

### WARNING!

The following operations may have a negative effect on the steering and stability of the scooter:

- Keeping and walking an animal attached to a leash.
- Carrying passengers (including animals).
- Hanging an object from the steering column.
- Towing or pushing by another vehicle.

### WARNING!

At all times when using your scooter, keep both hands on the handlebars and your feet on the footplate. This position provides the best control of the vehicle.

- Set the speed control dial to the desired speed.
- Press the selected side of the drive control lever with your thumb.
- The automatic electromechanical brake will disengage and the scooter will smoothly accelerate to the set speed.
- Pull the left handlebar grip to turn left.
- Pull the right handlebar grip to turn right.
- Move the steering column to the center position to go straight ahead.
- To stop, slowly release the drive control lever. The electronic brakes come on automatically when the scooter stops.

Note: The scooter, with the same speed control dial setting, moves slower in reverse than in forward gear.

## Getting off the scooter

### WARNING!

Never attempt to get on or off the scooter without first removing the key from the power switch. This prevents the scooter from being moved by accidentally switching the drive control lever.

1. Stop the scooter completely.
2. Remove the key from the power switch.
3. Release the seat rotation lever and rotate the seat so that it faces the side of the scooter.
4. Make sure the seat is locked in this position.
5. Carefully and safely get off the seat and stand next to the scooter.
6. To make it easier to get on the scooter, you can leave the seat facing to the side of the scooter.

## Power off function after a certain time

To keep the life of the batteries, the scooter is equipped with an energy-saving function to automatically turn off the power after a certain period of time. If you accidentally leave the key inserted in the power switch, but do not use it within a 15-minute interval, the scooter controller will turn off the power automatically.

If the timed shutdown function has been activated, you will need to do the following to return to normal use:

1. Slide the key out of the power switch.
2. Slide the key back in to power up the scooter.

# Regulation and adjustment of scooter components

## Delta steering column and rear view mirrors adjustment (Option)


Figure 3. Delta steering column and rear-view mirrors

It is possible to adjust the rear-view mirrors:

1. Slide the cap upward ③.
2. Loosen the mirror nut.
3. Grasp gently ②
4. Rotate the tube ② and set the mirror ① in the desired position.

## Steering column angle adjustment


Figure 4. Steering column adjustment dial

Steering column - option 1: The steering column is fixed with a locking device. It supports two positions, vertical locking or unlocking . The tilt angle is not adjustable.

- To fold the scooter, slide the lock out of the opening of the steering column.
- Hold the steering column in its highest position and release the lock. The lock will automatically return to the hole and secure it in the steering column

Steering column - option 2: The column is locked with a knob. This supports the ability to adjust the angle. The height is adjustable according to the user's requirements

- Turn the knob counterclockwise which will release the mount and allow you to fold the scooter.
- Keep the steering column at the correct height and tighten the knob clockwise to secure the position.

**WARNING!**

Remove the key from the power switch before adjusting the steering column or seat. It is prohibited to adjust the steering column or seat if the scooter is in motion.

**ATTENTION!**

To completely lower the steering column, you must first remove the seat and battery pack. Remove the seat by pulling it upward off the scooter.

## Seat rotation adjustment

The seat rotation lever ② located on the right side of the base, secures the seat in several positions. See Figure 5.

1. Pull the seat rotation lever ② to unlock the seat.
2. Rotate the seat to the desired position.
3. Release the lever to lock the seat in place.


Figure 5. Seat adjustment

## Armrest spacing adjustment

The spacing of the armrests can be adjusted inward or outward.

1. Loosen the armrest adjustment knobs ① - Figure 5.
2. Move the armrests to the desired width.
3. Tighten the armrest adjustment knobs.

**ATTENTION!**

Raise the armrests to make it easier to get on and off the scooter.

## Seat height adjustment

The seat can be set on different height.  
See figure 6.

1. Remove the seat and battery pack from the scooter.
2. Use the attached ring to pull and remove the locking pin from the bottom of the seat post.
3. Raise or lower the upper seat post to the desired height.
4. Keep the upper post at the height where the adjustment holes align with the lower post.
5. Fully insert the locking pin.
6. Install the seat and battery pack.

### Identification key

1. Upper seat post
2. Locking pin
3. Lower seat post


Figure 6. Seat height adjustment

## Disassembly and assembly the scooter

### Disassembly

You can disassemble the scooter to several components: ① seat, ② rear frame, ③ front frame, ④ basket, ⑤ battery pack. Disassembly and assembly does not require any tools. Always carry out disassembly on a level, dry surface in a room big enough to move freely, at a distance of 1.5m in each direction from the scooter.

### Disassembly procedures

#### WARNING!

Do not lift weights beyond your physical capabilities. When disassembling and assembling the scooter, if necessary, ask for assistance.

### Identification key

1. Seat
2. Rear frame
3. Battery pack
4. Front frame
5. Basket


Figure 7. Disassembled scooter

1. Remove the seat by lifting it completely above the scooter. If you encounter resistance when removing the seat, lift the rotation adjustment lever and rotate the seat left and right while lifting up.
2. Remove the battery pack by lifting it off the scooter.

## Disassembly the frame (MINI II - D only)

### ATTENTION!

Before attempting to disassemble, bend the rear of the scooter over the anti-tip wheels, as shown in Figure 8.


Figure 8. Frame release lever

1. After removing the seat and battery pack (see the removal procedure on the previous page), lower the steering column and tighten the steering column adjustment knob. See Figure 4.
2. Lift the frame release lever and move the frame half forward until the hooks are no longer resting on the rear frame bar.
3. Slowly separate the two parts of the scooter frame. See Figure 8.

## Assembly

1. Use the front frame release lever (MINI II-D only) to raise and align the front frame hooks to the rear frame bar. See Figure 8.
2. When the front frame hooks are over the rear frame bar, lower the front frame while tilting the rear frame. This will automatically activate the locking mechanism. See Figure 8.
3. Raise the steering column and fully tighten the steering column adjustment knob. See Figure 4.
4. Install the battery pack.
5. Install the seat and rotate until it locks in place.

### WARNING!

Be especially careful when unfolding and folding the wheelchair, lack of attention may cause body parts to be jammed by moving and folding parts!

## Basic troubleshooting

### WARNING!

If the battery symbol is flashing on the battery level indicator, the scooter has detected a fault. If this happens, follow the procedures below:

1. Release the drive control lever, then switch the scooter off and on again. If this does not clear the fault, proceed below to step 2.
2. Switch off the scooter and fully charge the battery. If this does not clear the fault, go below to step 3.
3. Count the number of flashes. The battery symbol will flash several times before an extended pause. Count the number of flashes between breaks and then pass this on to your dealer. The following fault codes are listed below:

Number of flashes	Fault description	Meaning and procedure
1	Battery low	Batteries are discharged • Recharge the batteries
2	Low battery fault	Batteries are discharged. Recharge the batteries. Check the batteries and their associated wires and connections.
3	High battery fault	The battery voltage is too high. This could be the result of overcharging and/or moving down the slope. When moving down the slope, reduce speed to minimize regenerative charging.
4	Current Limit Time-out or Controller too hot	Motor exceeds its maximum rated current for too long. • The scooter may have stopped. Turn the power off, wait a few minutes, then turn the power back on. • The motor may be defective. Check the motor and the cables and their connections.
5	Failure of the parking brake	The parking brake release switch is active or the parking brake is faulty. • Check the parking brake cables and their connections. • Ensure any associated switches are in their correct positions.

6	Drive lock	<p>The stop function is active or a charger locked or OONAP conditions have occurred.</p> <ul style="list-style-type: none"> <li>• Release stop condition (e.g. seat raised)</li> <li>• Disconnect the battery charger.</li> <li>• Ensure that the drive control lever is in the neutral position when starting.</li> <li>• The drive control lever may need to be calibrated.</li> </ul>
7	Failure of the speed potentiometer	<p>The drive control lever, speed potentiometer, SRW or wires and their connections may have failed.</p> <ul style="list-style-type: none"> <li>• Check the drive control lever, speed potentiometer and wires and their connections.</li> </ul>
8	Motor voltage fault	<p>Motor or wiring and connections have failed.</p> <ul style="list-style-type: none"> <li>• Check the motor and wiring and its connections.</li> </ul>
9	other error	<p>The controller has internal failure.</p> <ul style="list-style-type: none"> <li>• Check all connections and wiring.</li> </ul>

## Batteries and charging

Before charging the batteries, read the instructions and prepare the charger. The battery status indicator on the steering wheel console indicates the approximate power of the batteries by means of a colour code. From right to left, where green indicates fully charged, yellow indicates that some power has been used and red indicates an immediate need to charge the batteries. To determine this exactly, check the battery charge indicator while driving on dry, level ground at top speed.

### First charge


The scooter requires two long, full charge cycles (12 volt), sealed and maintenance-free batteries.

They are charged by the charging system outside the vehicle. Before first use, first charge the scooter batteries for 8 to 14 hours. Keep the batteries fully charged so that the scooter can run smoothly.

## Charging the scooter

The scooter can be charged through the charging socket located on the battery pack. Connect the charger to the scooter or to the battery pack and then plug into an electrical outlet. When the batteries are fully charged, disconnect the charger from the scooter.

## Operating temperatures

The scooter is designed for use, ideally in temperatures between  $-10^{\circ}\text{C}$  ~  $50^{\circ}\text{C}$  ( $-14^{\circ}\text{F}$  to  $122^{\circ}\text{F}$ ). In extremely cold conditions, the batteries may freeze, depending on charge, use and composition. In extremely high temperatures, the scooter may run at a lower speed. This is due to a built-in mechanism to protect against damage to the gearbox and other electronic components when operating in hot conditions.

## Battery operation and charging

If you use your electric scooter on a daily basis, charge the batteries as soon as possible after you have finished using it for the day. The scooter will be ready every morning. When using the scooter on a daily basis, we recommend charging for 8 to 14 hours. We also recommend that you continue to charge for 4 hours after the charger has indicated the end of the cycle. If you use your scooter once a week or less frequently, charge the batteries at least once a week for 24 hours.

## Charging the batteries

### SUGGESTION!

Keep the batteries charged and avoid discharging them completely. To keep the batteries in good condition and to extend their life, we recommend charging a minimum for 48 hours, once a month.

The LED light on the charger indicates the different charging status.

Red light - means that, it is in the process of charging.

Green light - means that, bulk charging (95% charged) is complete and the charger continues to cycle charging.

#### ATTENTION!

There is a charger in the scooter which blocks its functions. The scooter will not start and the battery charge indicator will not be active when charging the battery pack fitted to the scooter.

**WARNING** - Removal of the ground pins may cause an electrical hazard. If necessary, when accessing a 2-pin socket, use a suitable, approved 3-pin adapter.

**WARNING** - Never use an extension cord to plug in the charger. Plug the charger directly into a properly installed standard electrical outlet.

**WARNING** - It is prohibited, while charging the batteries, for children to play unsupervised near the scooter. We recommend that you do not charge the scooter batteries if the scooter is occupied.

**WARNING** - Explosive gases may be generated when charging the batteries. Keep the scooter and charger away from sources that can cause ignition such as fire, sparks and ensure proper ventilation when charging the batteries.

**WARNING** - Charge the scooter batteries with the supplied external charger. Do not use any other type of charger.

**WARNING** - Before each use, check the charger, cables and connectors for damage. If you find any damage, contact your supplier.

**WARNING** - Do not open the battery charger case. If the charger does not work correctly, contact your authorized supplier.

**Warning** - Be aware that during charging, the battery charger case may be hot. Avoid contact with skin and other surfaces that may be damaged by heat.

**WARNING** - If the charger has not been tested and approved for outdoor use then it should not be exposed to extreme weather conditions. If the battery charger is exposed to bad or extreme weather conditions, allow the charger to adjust to differences in environmental conditions before using it indoors.

## Care and maintenance

Your scooter requires minimal care and maintenance. If you are not sure of your care and maintenance capabilities listed below, you can have the inspection and service carried out by an authorized supplier. The components listed below require periodic inspection and maintenance.

### Tires

Regularly check your electric scooter's tires for signs of damage or wear.

### Wheel Replacement

The scooter is equipped with polyurethane (PU) tires. Replace the entire wheel if the tire has become worn or damaged. For information on replacing wheels, contact your authorized electric scooter supplier.

**WARNING!**

The wheels on the scooter must be serviced/replaced by a qualified technician. During this procedure, make sure that, the key is removed from the power switch and the scooter does not have the clutch lever switched disengaging the drive.

### External surfaces

Bumpers, tires and rims can be occasionally cleaned with rubber and vinyl cleaner.

### Cleaning and disinfection

Use a damp cloth and a mild, non-abrasive cleaner to clean the plastic and metal parts of the electric scooter. Avoid using products that can scratch the scooter's surfaces. If necessary, use approved disinfectants. Before using, make sure the disinfectant is safe for the scooter's surface.

**WARNING!**

Follow the safety instructions provided by the manufacturer of cleaning and disinfecting products. Not complying may result in skin irritation or early damage to the upholstery or electric scooter.

## **Battery terminal connections**

Make sure the connections remain tight and non-corroded.

## **Electric wiring harnesses**

- Regularly check wire connections.
- Regularly check the insulation of the wires, including those from the charger, for wear or damage.
- Have an authorized supplier repair or replace any damaged connectors, connections or insulation before using your electric scooter again.

## **Daily checks**

- With the power off, check the operation of the drive control lever. Make sure that the lever is not bent or damaged and that when released, it returns to the neutral position. Self-repair is prohibited. In case of problems, contact an authorized supplier.
- Visually inspect the steering column wiring. Make sure it is not burnt, cut or wires sticking out. Contact your authorized supplier if there are problems.
- Check for worn (flat) spots on the tires. Worn tires can affect stability.
- Check armrests for loose, strained or damaged parts. If there are problems, contact your authorized supplier.
- Check brakes. Checks should be made on a level surface without obstacles within 3 feet (1 meter) of the scooter.

## **Weekly checks**

- Check the controller and charger connections for corrosion. If necessary, contact an authorized supplier.
- Check the pressure if there is a pneumatic tire on the equipment. If the tire does not hold pressure, contact an authorized supplier to replace the inner tube.

## **Monthly checks**

- Check that the anti-tip wheels do not rub against the ground when operating the electric scooter.
- Check the anti-tip wheels for wear and tear. Replace if necessary.
- Check for tire wear. Contact an authorized supplier for repairs.
- Keep the scooter clean and away from foreign objects such as mud, sand, hair, food, drinks, etc.

## **Yearly checks**

- Take your electric scooter to an authorized supplier for an yearly inspection, especially if you use your scooter daily. This will ensure the correct operation of the scooter and help avoid future breakdowns.

## **ABS plastic covers**

If your electric scooter has a shiny chassis surface then it has been treated with a clear fixative coating. You can apply a small coat of car wax to help restore the glossy surface. If your scooter's chassis has a matte finish, use only products designed for matte finish surfaces. Do not use wax or other products designed for glossy paintwork.

## **Axle bearings and motor/transaxle**

All these components are pre-lubricated, sealed and do not require further lubrication.

## **Console, Charger, and Rear Electronics**

Keep these areas free of moisture. Allow these areas to dry thoroughly if they have been exposed to moisture before operating your electric scooter again.

## **Nylon Lock Nut Replacement**

During periodic maintenance, assembly or disassembly, each securing nylon nut must be replaced with a new one. Nylon nuts must not be reused as this may damage the nylon threads and reduce the effectiveness of the fit. Replacement nylon nuts are available from local hardware stores or an authorized supplier.

## **Storage of the scooter**

If you plan not to use the scooter for a period of time then follow these rules:

- Fully charge batteries before storage.
- Remove the battery pack from the electric scooter.
- Store the scooter in a warm, dry room.
- Avoid storing in a place where the scooter will be exposed to temperature changes.
- Recommended storage temperatures: -10°C to 50°C (-14°F to 122°F).

Batteries that are regularly completely discharged, not regularly charged, stored at inappropriate temperatures, or stored without being charged may be permanently damaged resulting in unreliable operation and shortened battery life. To ensure proper operation, it is recommended that batteries be charged regularly, throughout the period of prolonged storage of the scooter.

For longer storage, several boards can be placed under the frame of the electric scooter to lift it above the ground. This will reduce the weight resting on the tires and reduce the probability of flat tires during contact with the ground.

## Disposal of an electric scooter

The electric scooter must be disposed of in accordance with applicable local and national regulations. For information on the proper disposal of packaging, metal frame parts, plastic parts, electronics, batteries, neoprene, silicone and polyurethane materials, please contact your local waste disposal company or authorized supplier.

## Technical parameters

Model	MINI II - D	MINI II - F
Type	Scooter, rear wheel drive, Class A	
Maximum user weight	136 kg	
Maximum speed (driving forward)	6 km/h 8 km/h	
Minimum braking distance at max. speed	6 km/h: 1000mm 8 km/h: 1500mm	
Range	24V 12Ah Gell batteries: 12 km	
Maximum obstacle height	45 mm	
Maximum safe inclination	6° (10,5%)	
Minimum turning radius	1150 mm	
Minimum width when reversing	1250 mm	
Dimensions unfolded (L x W x H)	1030mm x 490mm x 880mm	1050mm x 490mm x 880mm
Total mass with batteries	43.5 Kg (24V/12Ah Gel)	41.5 Kg (24V/12Ah Gel)
Battery weight	9 kg	9 kg
Maximum weight of single part	15 kg	24 kg
Effective seat depth	370 mm	
Effective seat width	410 mm	
Seat surface height at front edge	380 mm	
Footplate height	110 mm	
Seat surface angle	-8°	
Motor	270 Watt, electromagnetic brake	
Batteries	Gel Batteries 24V/12Ah	
Controller	Dynamic 50A	
Protection class	IPX5	

Battery charging	Input AC100-240VAC, Output DC29.4V, 2Amp
Protection class of battery charger	IPX1
Insulation class of battery charger	Class 2
Maximum motor sound level	68 dB
EMC compatible according to	ISO 7176-21
Ignitability of upholstery according to	EN 1021-2
Diameter rear wheels (number)	8x2.5 cala
Diameter front wheels (number)	8x2.5 cala
Suspension	N/A
Control (wig-wag)	Lever to start/stop movement
Operating console	Rotatable switch for speed control
Storage and use temperature	-10°C to +50°C
Humidity to store and use	30% - 70%
<p>We reserve the right to make technical changes.</p> <p>The tolerance in the measurements is <math>\pm 15 \text{ mm} / 1.5 \text{ kg} / 1.5^\circ</math>. The theoretical turning radius decreases if the scooter is often used on inclines, uneven surfaces or for driving up curbs. The maximum distance covered has been tested under ideal conditions in accordance with ISO7176-4.</p>	

Speed and range depends on the user's weight, type of terrain, charge and battery condition. Information is current at the time of publication, the manufacturer reserves the right to change specifications.

## Warranty

### Two (2) years limited warranty

Structural frame components, including: seat support tube and welded frames.

### One (1) year limited warranty

The MINI II scooter is warranted for a total of (12) twelve months from the date of purchase and against defects due to defects in design and materials. This warranty does not reduce your rights, but is in addition to them.

All electronic components, including controllers and battery chargers, have a one (1) year warranty. Powertrain, including: differential, motor and brakes. Servicing of the controller and battery charger must be performed by an authorized supplier. Any attempt to open or disassemble these components will void the component warranty.

## **Warranty does not include**

The warranty does not cover components that may need to be replaced due to normal wear and tear during operation (tires, belts, light bulbs, upholstery, plastic bumpers, engine brushes, fuses and batteries) or have been damaged during misuse or accident for which the manufacturer or supplier is not responsible. This warranty does not cover labor or service calls.

## **Batteries**

Batteries are covered by a six (6) month warranty from the manufacturer.

NOTE - The warranty does not cover deterioration of battery performance caused by being left in a fully discharged state, in cold conditions for an extended period of time, or by excessive use.

## **Inspection and warranty service**

Warranty service can be performed by an authorized supplier. Please contact your authorized supplier for the current costs associated with a service visit.

## **Replacement parts**

The availability of replacement components and/or spare parts is dependent on the capabilities of the supplier, not the manufacturer. For more information on replacement components and/or spare parts, please contact your authorized supplier.

	Unique Device Identifier
	Follow the instructions for use
	Refer to the instructions for use
	Indicates the maximum weight capacity. Do not exceed the specified weight capacity of the product
	Warning
	Indicates the danger of electric shock
	EMI-RFI - This product has been tested and passed at an immunity level of 20 V/m
	Indicates the risk of explosion
	Corrosive chemicals contained in batteries
	Indicates a mandatory action that should be performed as specified in order to avoid a hazardous condition/situation. Failure to perform mandatory actions can lead to personal injury and/or equipment damage
	Protect from humidity
	Possibility of jamming fingers
	Pinch/Crush danger points
	Permitted temperature

	For indoor/outdoor use
	Contains Lead
	Maximum user weight
	Maximum safe inclination
	Indicates attachment points on the scooter
	The product is designed to be transported by car
	The product is not intended to be used as a seat for transport in a motor vehicle
	Do not stand on the footrest
	Use only AGM or Gel-Cell batteries to reduce the risk of leakage or explosive conditions
	Wear eye protection
	N = Neutral (Brakes disengaged; lever disengaging the drive turned on)
	D = Drive (Brakes engaged; unit in Drive Mode)
	Indicates a prohibited action that should not be performed at any time or in any circumstances. Performing a prohibited action can cause personal injury and/or equipment damage
	Do not use batteries with different capacity (Ah). Do not mix old with new batteries. Always replace both batteries at the same time

	Do not mix AGM and Gel-Cell batteries. Always use batteries of the same type and chemistry
	Keep tools and other metal objects away from battery terminals. Contact with tools can cause electric shock
	It is forbidden for children to play, unsupervised, while charging the scooter's batteries
	Removal of the grounding pin may cause an electrical hazard. If necessary, properly install an approved 3-pin adapter in an electrical outlet with access to a 2-pin plug
	Do not connect an extension cord to the AC/DC converter or the battery charger
	Indicates flammable material. Do not expose to heat sources such as open flame or sparks. Do not transport batteries with flammable or combustible items
	Avoid exposure to rain, snow, ice, salt, or standing water whenever possible. Maintain and store in a clean and dry condition
	Do not remove anti-tip wheels
	Keep your hands away from the tires when driving. Be aware that loose fitting clothing can become caught in the drive tires
	Do not operate the drive disengagement clutch on an incline
	Do not operate the drive disengagement clutch on a decline
	Contact your local recycling center or authorized provider for information on proper recycling of product components

## DISPOSAL OF WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT

This symbol on the products or accompanying documents means that used electrical or electronic products should not be disposed of with the normal household waste. For proper disposal, refurbishment or recycling, please take these products to the collection points for this type of waste.

Alternatively, in some EU countries or other European countries, you may return your product to your local retailer when purchasing a similar new product. By disposing of this product correctly, you will help to conserve valuable natural resources and support the prevention of potential negative effects on the environment and human health, which could result from incorrect disposal. For further information, please contact your local municipality or collection point waste collection points. Penalties may be imposed in accordance with local regulations if this product is disposed of in an inappropriate manner. For parties within the European Union, If you wish to discard electrical or electronic equipment, please obtain the necessary information from your retailer or supplier.

Elimination in countries outside the European Union. This symbol is valid in the European Union. If you wish to decommission this appliance, please obtain the correct decommissioning information from your local authorities or from your retailer.

## Contact details

### Manufacturer

mdh Sp. z o.o.

Address: ul. Maratońska 104, 94-007 Łódź, Poland

tel. +48 42 674 83 84, fax. +48 42 636 52 21

[www.mdh.pl](http://www.mdh.pl)

[www.viteacare.com](http://www.viteacare.com)

We reserve the right to make technical changes due to technical progress.