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
**English version** on next page

**Users with visual impairments can find the PDF files together with further information on our website:**

**< [www.viteacare.com](http://www.viteacare.com) >.**

**👉 Contact your specialist dealer when required.**

**Alternatively users with visual impairments can have the documentation read out by a helper.**

# Table of contents

<b>Introduction .....</b>	<b>44</b>
Indications.....	45
Contraindications .....	45
<b>Safety of use .....</b>	<b>45</b>
Information on electromagnetic compatibility (EMC) .....	48
Safety control .....	48
<b>Product description .....</b>	<b>49</b>
Wheelchair description.....	49
Technical specifications.....	51
<b>Practical advice on use.....</b>	<b>52</b>
Preparing the wheelchair.....	52
Before use.....	52
Folding the backrest.....	53
Assembly of elements.....	53
Footrest length adjustment .....	54
Control Panel operation .....	55
Joystick - control handle .....	55
OPTION - Rear caregiver control panel .....	56
OPTION - Wireless remote control.....	58
Position lights.....	62
Armrests .....	63
Drive disengagement clutch .....	63
<b>Safety .....</b>	<b>64</b>
Safety belt.....	64
Anti-tip wheels.....	65
Transporting the wheelchair .....	65
Inflating tires.....	66
<b>Use .....</b>	<b>66</b>
Getting in the wheelchair .....	66
Wheelchair driving training .....	67
Driving on uneven terrain .....	67
Getting out of the wheelchair – end of product use .....	67
<b>Batteries.....</b>	<b>68</b>
Safety of battery charging.....	68
Charging process .....	69
Batteries - Storage .....	70
<b>Service and repair .....</b>	<b>70</b>
Modification.....	70
<b>Care and maintenance .....</b>	<b>70</b>
Cleaning .....	70
Disinfection.....	71
Reuse (reusable product) .....	71
Product life cycle.....	71
Disposal of the electric wheelchair .....	72
Maintenance plan.....	72
Replacing motor brushes.....	74
Fault codes and troubleshooting.....	74
<b>Symbols.....</b>	<b>77</b>
<b>Contact Information.....</b>	<b>80</b>

## Introduction

Cross Smart is a wheelchair created in accordance with the latest technologies. Its design allows for maximum adaptation to the user's needs. It combines a robust and durable construction, modern electronic control solutions and safety systems. In order to maintain safety and maintain the full efficiency of the product as long as possible, it should be regularly inspected at the service point (at least once a year). In case of any doubts or damage, please contact the distributor. Make sure that you have filled the attached warranty form and that you have submitted it to the distributor, as this will enable us to ensure the best possible service.

# ATTENTION

Always observe the basic precautions according to the information given below regarding safety when using the product. Read the entire manual carefully before use. Please keep these instructions for future use.

### WARNING!

Improper use can lead to death, serious injury, dangerous conditions or product malfunction.

### ATTENTION!

Improper use may lead to personal injury and/or damage to the product.

### SUGGESTIONS!

Follow the instructions below to keep the product in good working order.

If the device is used by another person, a copy of these instructions for use must be provided.

These instructions for use apply to the design and specification of the product at the time of publication. Due to design changes, some illustrations and photos in the manual may not correspond to the product purchased. We reserve the right to make design changes.

The manufacturer disclaims any liability for personal injury or property damage that may result from improper or unsafe use of this product.

## Indications

An electric wheelchair is a device that improves the quality of life, designed to provide mobility.

## Contraindications

Contraindications to the sitting position.

Suggestion!

To maintain the life of the battery, during transport, they will be disconnected. Remember to connect the red cable with the input marked with the "+" symbol before using the wheelchair for the first time.

Suggestion!

Avoid completely discharging the batteries. Charge the battery after all-day use and at least once every 2 weeks.

## Safety of use

**WARNING!**

Every "serious incident" that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State where the user resides.

**ATTENTION!**

The product is intended for disabled people with musculoskeletal conditions to increase their mobility. If you are not sure whether your health condition allows you to use the device, consult a health care professional.

**WARNING!**  
Do not modify the product in any way; this will void the warranty.

**WARNING!**  
Exceeding the maximum user weight will void the warranty. The manufacturer is not responsible for injury and damage if the maximum user weight is exceeded.

**WARNING!**  
Please read these instructions for use before using the wheelchair and keep them for future reference.

1. Read and observe the information contained in the instructions for use. It is forbidden to use this product and the available optional equipment without first carefully reading and understanding these instructions. If the warnings, cautions or instructions are incomprehensible, contact the supplier, otherwise personal injury or an accident may occur. Never attempt to use the wheelchair beyond its limitations described in this instructions for use.
2. Never turn off the control panel while driving, as this may cause the wheelchair to come to a sudden stop, which may result in injury.
3. Consult a doctor if the medications taken do not affect the ability to drive motor vehicles.
4. To operate the electric wheelchair safely, you must drive on flat surfaces where all four wheels are in contact with the ground.
5. If the battery level indicator drops to the red level, charge the wheelchair as soon as possible. Completely discharging the batteries can damage them or shorten their lifetime.
6. Keep metal objects away from battery connections. Danger of electric shock.
7. When driving on curves, always slow down and maintain a stable centre of gravity. When taking sharp turns, always reduce the speed.
8. Take extra caution when reversing. For safety, the reverse speed should be limited to the minimum speed.
9. Children are not allowed to play near the wheelchair or operate it.
10. Some parts of the electric wheelchair are sensitive to temperature changes.
11. Remember that the temperature of the components and frame parts may increase when exposed to the sun or another heat source, and at low temperatures, there is a risk that the metallic parts of the wheelchair will be very cold. This may cause burns or frostbite.
12. Store the wheelchair in clean and dry conditions.
13. It is forbidden to drive at night without proper lighting.
14. Do not remove anti-tip wheels or drive wheels.
15. Keep hands and feet away from moving parts when driving. Ensure that loose clothing does not become trapped in the drive wheels.

16. Do not connect the battery charger to the extension cord.
17. Disassembly of the controller, engine or charger by persons other than the authorized service is prohibited and invalidates all applicable warranty conditions.
18. It is forbidden to move the wheelchair on motorways. When moving on public roads, follow all current traffic rules for the country.
19. It is forbidden to use the electric wheelchair under the influence of alcohol, drugs or other substances affecting the ability to drive the vehicle.
20. Take caution when driving in busy areas and shopping centres.
21. Never use the wheelchair as a chair for transport in motor vehicles.
22. Do not try to lift the wheelchair by holding onto parts such as the seat, armrest or parts of the housing.
23. When driving, it is forbidden to use a mobile phone or other personal communication devices (e.g. CB-radio).
24. Do not use the wheelchair outside during rain or snowfall, high (above 30° C) and low (below 0° C) temperatures.
25. Do not drive the wheelchair onto unprotected wharves, river banks or lake shores.
26. The wheelchair is not suitable for use in difficult terrain (unpaved roads).
27. Do not stop, turn or reverse suddenly to avoid falls and injuries.
28. When driving, keep your forearms on the armrests and feet on the footrest platform.
29. It is forbidden to stand while driving the wheelchair. Do not lean to the side or lift objects from the ground.
30. Do not drive onto raised objects (steps, kerbs) higher than 50 mm.
31. The wheelchair is not suitable for moving up or down the stairs. Always use ramps or elevators.
32. The wheelchair is not designed to negotiate ground gaps larger than 100 mm.

**WARNING!**

Drive the electric wheelchair in a straight line up or down the slope, never reverse. Irregular forward and backward movements with slope changes can increase the risk of the wheelchair tipping over.

The maximum user weight is 120 kg

## Information on electromagnetic compatibility (EMC)

Despite complying with applicable EMC directives and standards, an electrically powered wheelchair may interfere with the operation of other electrical devices, such as mobile phones or alarm systems. In addition, this may affect the performance of the wheelchair itself. If you notice this behavior, it is a recommended idea to switch off the third-party device or power wheelchair and move it out of interference range. Remember that other devices emitting electromagnetic fields can also affect the operation of the wheelchair.

## Safety control

**ATTENTION!**

The following points should be checked before use:

- Electrical connections and wiring harnesses – protect and remove corrosion.
- Batteries – Fully charged and not corroded.
- Brakes - check that the handbrake is working.
- Folding lock mechanism - make sure wheelchair is fully unfolded.
- Condition of the tyres – make sure that they are not damaged and excessively worn - if there is a problem, contact the distributor.

## Product description

The electric wheelchair for people with disabilities has been designed for people who experience significant limitation of movement. The control is carried out using a joystick or other control device that allows you to set the direction of movement. It is a solution for people who cannot move in a traditional wheelchair.

The convenient controller allows you to fully maneuver the electric wheelchair, especially since it is very maneuverable. Electric wheelchairs can successfully rotate around their own axis.

## Wheelchair description

Illustration 1. Front of the wheelchair.


Illustration 2. The back of the wheelchair.


Illustration 3. Control panel (joystick)

## Technical specifications

Model	CROSS SMART
Seat width	430 mm
Maximum user weight	120 kg
Maximum load capacity	120 kg
Maximum speed forward	6 km/h
Maximum speed reverse	2,5 km/h
Energy range	20 km ± 10%
Battery / weight	20Ah/24V (gel) / 13 kg
Engine power	2 x 250W
Seat surface height from the floor	540 mm (with cushion)
Charger (output power)	24V/2A
Seat Depth	430 mm
Backrest height	440 mm
The total length of the wheelchair	1100 mm
Total wheelchair height	940 mm
Total width of the wheelchair	660 mm
Width when folded	440 mm
Armrest height from the ground	720 mm
Seat cushion thickness	50 mm
Dimensions of the folded wheelchair for transport	L710x W.500 x H.790 mm
Front wheels	200 x 50 mm (8" PU)
Rear wheels	320 x 60 mm (12" pneumatic)
Turning radius	1200 mm
Footrests	Swivel, adjustable
Weight of the wheelchair with batteries / without	55 kg / 42 kg
Safety	Adjustable safety belt, anti-tip wheels
Maximum safe inclination	8° (14%)
Other functions	Folding backrest for transport, reclining footrests, reclining armrests, cross-folding frame, front shock absorbers
Minimum braking distance at maximum speed	1000 mm
Frame material	Steel
Description of the environment in which the product will be used	  
Storage and usage temperature	Od +5°C do +40°C
Storage and usage humidity	From 30% to 70%
Tools necessary for adjustment	size 5 hex key

Speed, max distance of the vehicle depends on the user weight, battery load, battery capacity and the terrain on which the vehicle is used. The information is current at the time of publication, the manufacturer reserves the right to make changes to the specifications.

# Practical advice on use

## Preparing the wheelchair

The wheelchair is delivered in a cardboard package, and the fragile parts are additionally protected against damage during transport. Unpack it carefully to avoid any damage.

### Before use

- Check all electrical connections, make sure that they are stable and not corroded.
- Check that the batteries are fully charged.
- Switch on the lighting – particularly recommended at night, at dusk and in conditions of limited visibility (fog, rain, snow) – so that the wheelchair user is clearly visible to other road users.
- Check that the wheelchair is engaged in drive mode (clutch levers engaging the drive in the upper position).
- Always plan your route in such a way as to avoid dangerous (marshy, difficult areas) and obstacles that you will not be able to overcome with a wheelchair.
- Whether the speed is adapted to the user's capabilities and skills.
- Does the horn work.
- Fasten seat belt.

#### ATTENTION!

- When transferring to the wheelchair, do not rest your body weight on the armrests as this can lead to the wheelchair tipping over and causing injury.
- Do not stand on the footrests because this could tip the wheelchair over.
- Remember to check whether the wheelchair power is switched off because unintentionally touching the joystick can cause the wheelchair to move, fall and, consequently cause injury.
- During the transfer, try to be as close to the wheelchair as possible to avoid slipping and falling.

#### WARNING!

The following actions can influence the steering and stability of the wheelchair negatively:

- Holding and walking an animal strapped on a leash.
- Carrying passengers (including pets).
- Hanging the item on the controller.
- Towing or pushing by another vehicle.

## Folding the backrest

- Press the lock button located on the backrest shown in illustration 4.

Illustration 4. Folding the backrest.

- Then pull the handles down to fold the backrest.
- To fold the backrest, perform the steps in the reverse order by lifting the backrest, grabbing the handles and then lifting it up until the folding lock does not click, which will be indicated by a click.

Illustration 5. Folded backrest.

## Assembly of elements

Assembly of the main parts should be carried out in accordance with the instructions given below:

- Install the battery by applying it to the bracket under the wheelchair seat (Illustration 6.). Note that the hooks are applied to the bracket rail.

Illustration 6. Install battery.

- Then connect the battery plug to the wheelchair socket as shown in illustration 7.

Illustration 7. Connecting the battery power connector.

- To install the footrest, hang the footrest mounting brackets on the footrest hangers (Illustration 8.)

Illustration 8. Preparing to install the footrest.

- Then lower the platforms and rotate the footrest until the footrest tilt lock is closed.

Illustration 9. Installation of the footrest.

- To tilt the footrests again, pull the lock lever and tilt the footrests outside the wheelchair. (Illustration 10.)

Illustration 10. Footrest tilt.

## Footrest length adjustment

- To adjust the length of the footrests, unscrew the locking screw shown in Illustration 11 using a 5 mm hex key.
- Then adjust the extension of the footrest rail and tighten the locking screw.
- Check that the footrest is firmly installed.

Illustration 11. Adjusting the length of the footrest.

## Control Panel operation

The electronic control panel is an element that allows you to steer the power wheelchair. It powers all operating systems. It allows control and steering of the wheelchair's movements, monitors the battery charge level and the status of the entire electrical system. Control panel is located on one of the armrests and is connected to the battery and motor.

- To start the electrical system of the wheelchair, press the power button (Illustration 3.).
- To change the speed of the wheelchair, press the left button to reduce the speed or the right button to increase the maximum speed (Illustration 3.). A single LED lighting up indicates the lowest speed. The highest speed is indicated by the lighting of all LEDs.
- In dangerous situations, an warning signal can be used by pressing the horn button (Illustration 3.).

## Joystick - control handle

Thanks to the joystick, you can control the direction and speed at which the wheelchair moves. When you move the joystick from the neutral position (middle position), the electromagnetic brake of the wheelchair will be released allowing it to move. The greater the swing from the neutral position, the faster the wheelchair will move. When the joystick returns to the neutral position, the brake system is activated, which slows down until the wheelchair comes to a complete stop. Moving the joystick backwards from the neutral position causes the wheelchair to move backwards. Swinging the joystick to the side causes the wheels to turn left or right.

### WARNING!

The wheelchair is equipped with electromagnetic brakes. When moving, be careful because the immediate return of the joystick to the neutral position stops the wheelchair at the shortest possible distance.

Installation and adjustment of the control panel (joystick):

- Unscrew the screw locking the control panel handle (Illustration 12.).
- Set the panel in the correct position.
- Then lock by tightening the screw.

Changing side of the control panel:

- Contact the distributor.

Illustration 12. Control panel adjustment screw.

## OPTION - Rear caregiver control panel

The electric rear wheelchair power steering controller is an intelligent motor drive control system. The controller uses intelligent control that provides a reliable, stable and efficient speed control scheme. The intuitive LCD display makes it easy to control a wheelchair. As with other electronic products, proper maintenance and use can ensure the best reliability of the control system.


Illustration 13. Description of the rear controller.

Illustration 14. Rear controller.

The rear controller has the following functions:

- power switch,
- speed control,
- choice of direction of wheelchair assistance,
- display that informs about the above functions.



Illustration 15. Rear controller dimensions given in millimeters.

- 1 Switching on/off the power and activity of the controller. Pressing the button activates the controller and turns on the display screen. Pressing the button again turns off the power and turns off the screen.

### WARNING!

Unless an emergency situation occurs, do not use the switch button to stop the wheelchair while driving - if used regularly in the above-mentioned situation, it may shorten the life of the wheelchair parts or cause personal injury.


**2** Increase / decrease wheelchair speed button. The caregiver can adjust the speed to their own preference or to the conditions while driving the wheelchair. The speed is divided into three sections , which are 20%, 40% and 60% of the maximum speed. The maximum speed should not exceed 6 km/h).

- The speed of gear 1 is 20% of the highest speed (one of the speed indicator light is on).
- Gear 2 is 40% of the highest speed (two speed light indicators are on).
- Gear 3 is 60% of the highest speed (three speed light indicators are on).

Pressing the button  will increase the speed. The button works in cyclic mode - after increasing the speed to gear 3, when pressed again, it will return to gear 1.

**3** Direction selection handle/joystick. The direction of movement of the joystick can control the direction of movement of the electric wheelchair. Direction and speed depends on the range of motion of the joystick.

- Pressing the joystick forward, will move the electric wheelchair forward.
- Pressing the joystick backwards will cause it to move backwards.
- If the joystick is not pressed, the controller will automatically turn off the power after the time set by the manufacturer and enter sleep mode. Press the power button to wake the controller from sleep mode.

**4** LCD screen displaying information on set functions.

### Rear controller maintenance

All parts of the electric wheelchair controller should be checked regularly for loosening, damage, or improper connection of terminals, connectors, and cables. All cables must be protected against damage.

Damaged parts must be replaced, for this purpose, contact the distributor.

- All parts of the controller must not be exposed to dust, dirt and liquid. Wipe them with a cloth soaked in warm water if necessary. Do not use solvent-based cleaners or abrasives.
- Do not attempt to open the housing or repair it yourself, otherwise the product warranty will be voided. If you have any question, please contact the manufacturer.

#### WARNING!

It is the user's responsibility to keep all components in working order at any time. Prolonged exposure to direct sunlight or chemical solvents used in the household can damage the plastic parts, making the controller inoperable.

### Rear controller safety

- It is forbidden to repair independently.
- It is the user's responsibility to keep the wheelchair in good working order and to ensure that all parts, including cables, are protected as much as possible to prevent damage.
- When the wheelchair is not in use, turn off the controller power.
- If the controller malfunctions or shows abnormal reactions, such as overheating, smoking and sparking of electrical wires, immediately turn off the power system and contact the distributor.
- In case of malfunction caused by electromagnetic interference, immediately turn off the power to the wheelchair controller.
- If the LEDs on the display flash while driving - this indicates an error, turn off the power system and contact the distributor.

### OPTION - Wireless remote control

1 Switching on / off the power and activity of the remote control. Pressing and holding the button for 1 second turns on the power of the remote control. Pressing and holding the button again turns off the power.

2 Control handle (joystick). Controls the movement of the wheelchair by swinging it in a given direction.

3 Button to reset the connection to the wheelchair controller.

4 Decrease speed button. Press to reduce the speed of the wheelchair.

5 Increase speed button. Press to increase the speed of the wheelchair.

6 Horn button. Press the button, in dangerous situations, to activate the audible warning.

7 LED light. It informs about the function activity by means of light indications.

8 Battery cover.



Illustration 16. Remote control description.

## LED light information signals.

<b>Operating status of the remote control</b>	<b>LED indicator color and status</b>
Activated	The red and green lights flash alternately for 1 second.
Off	The red light flashes quickly for 1 second.
Remote control and wheelchair connection	The light flashes continuously alternating red and green.
Active remote controller	The green light flashes slowly.
Low battery in remote control	The red light flashes slowly.

### **Safety of use the remote control**

- To reduce the risk of fire and electric shock, do not allow the remote control to get wet or damp.
- Do not use the remote control near water, do not immerse it in any liquid, do not allow any liquid to drip or spray on it.
- Do not install or use the remote control near heat sources (such as heaters, heating appliances, furnaces or other heat generating devices).
- Do not open or remove the enclosure for maintenance.
- Regardless of the cause of the damage to the remote control, it must be repaired by professional service personnel of the manufacturer.
- Before wiping and cleaning the remote control with a dry cloth, turn off the product and make sure that it is disconnected from the wheelchair controller.
- Do not use chemicals such as strong acids, alkalis, gasoline, alcohol, etc. to wipe the surface of the remote control. If cleaning is required, use neutral solvents or clean water.
- The shortwave radio signal emitted by wireless products may interfere with the normal use of other electronic devices. Please turn off the remote control in places that are not allowed to use this type of product.
- Do not use the remote control within 20 cm of electronic pacemakers or other medical devices, as radio waves may interfere with the normal use of electronic pacemakers or other medical devices.
- In order to ensure the safety and extend the life of the remote control when it is not used for a long time or during a storm, the battery should be removed from the remote control.

### **Use of the remote control**

#### **Preparation before connection:**

- Make sure that the wireless remote control is turned on (the signal lamp lights up when the device is turned on).
- Turn off all other controllers and turn on only the one to be connected. It is necessary to prevent interference of "Bluetooth" signals that can affect the efficiency of the connection.
- Turn on the power to the controller in the wheelchair to be connected and make sure that the wireless remote control is in the active state.


Illustration 17. Connection of the wireless remote control to the wheelchair controller.

**Automatic wireless remote control connection process:**

- Make sure that the controller in the wheelchair to be connected is turned on and the other controllers are disconnected from the power supply.
- Press and hold the power button of the wireless remote control for about 1 second and place it near the wheelchair controller. During this time, the control lamp of the wireless remote control flashes alternating red and green, indicating that the connection is in progress. Wheelchair controller will emit three short signals and the speed status will change to gear 2 and the wireless remote control status indicator will slowly flash green to indicate the end of the automatic connection.
- After successful connection, all available functions of the wheelchair can be controlled remotely by remote control. If the connection fails, you can try again and again until it is successful. In case of unsuccessful connection, contact the distributor.

**To check that the wireless connection has been made correctly:**

- Gently press the lever (joystick) of the wireless remote control to check whether the wheelchair is moving.
- Switch off and on the power of the remote control and observe if the controller in the wheelchair informs about a successful connection.
- Use the "speed up", "speed down" buttons on the wireless remote control to check whether the speed indicator of the controller in the wheelchair changes the displayed value.
- Use the "signal" button on the wireless remote control to check the correct operation of the acoustic warning signal in the wheelchair.

**WARNING!**

To improve the safety and functionality of the remote control, if it is not active for 35 seconds, it will turn off automatically. After switching on the remote control, if the control panel in the wheelchair is not detected, wait 35 seconds before automatically switching off the remote control. If you want to control the wheelchair again, turn on the control panel in the wheelchair again and then the remote control. To ensure a certain level of safety of use, when using wireless wheelchair controller, the speed is automatically adjusted to the highest speed in the second gear.

### Deleting a saved connection (preparing for a new connection):

- Turn off the power of the wireless remote control.
- Press and hold the buttons **3** and **1** at the same time for more than 3 seconds as shown in illustration 18, and the red and green light of the indicator lights will illuminate 3 times (if the wheelchair controller is on, the connection to the remote control will be automatically removed).
- After deleting the connection, the status indicator light flashes alternating red and green, and the disconnection is successfully completed. The wireless remote control will no longer automatically connect to the wheelchair controller.
- Long press of the power button will disable the remote control.
- If deleting of the connection fails, try to repeat the operation several times. If the repeated attempts do not bring the correct result, contact the distributor.




Illustration 18. Deleting saved connections.


### Install the wireless remote control battery.

1. Hold the battery cover and pull it down to remove.

2. Install the battery in the battery compartment (pay attention to the correct alignment of the positive and negative battery poles).

3. Place the cover on the battery shelf and slide to close completely.

## Position lights

The wheelchair has manually operated front and rear lights.

The position lights define the 'outline' of the wheelchair – they show its presence, width and direction, allowing pedestrians and drivers to see the person in the wheelchair from a distance. This increases safety after dark, at dawn, in fog and in other conditions of limited visibility, even when the wheelchair is moving slowly.

Illustration 19. Front light.

Illustration 20. Rear light.

To turn on the lights:

- First, check whether the power plug is connected to the battery box (Illustration 21).

### SUGGESTION!

The socket used to power the lighting is also the socket for connecting the charger.

- Locate the power switch on the right side of the wheelchair frame, below the armrest (Illustration 22). After turning on the power, ensure that all lights are working properly.

Illustration 21. Lighting power connection.

### ATTENTION!

The lights are powered independently of the control unit (joystick). If the wheelchair is not in use, it is recommended to switch off the position lights.


Illustration 22. Light power switch.

## Armrests

- To tilt the armrest, pull the locking lever shown in illustration 23. and then tilt the armrest back (Illustration 24.).

Illustration 23. Reclining armrests.

Illustration 24. Reclining armrests.

### WARNING!

Special care should be taken when folding the wheelchair to ensure that no clothing or body parts are at risk of being trapped. Take special care.

## Drive disengagement clutch

There are two clutch levers in the wheelchair to enable or disengage the electric drive and manually push the wheelchair.

- To disconnect the drive for manual pushing, lower the lever as far down as possible - position **B** in Illustration 25.
- Restarting the drive (travel mode) should be done by lifting the lever up - position **A**.



Illustration 25. Disconnect the motor drive.

### WARNING!

Do not turn off the electric drive if there is no person who can operate the wheelchair.

Do not switch off the electric drive when the wheelchair is on an incline as it could roll away with the risk of injury.

**WARNING!**

**Remember that when the electric drive is switched off, the electromagnetic brake system is also switched off.**

Always stand behind the wheelchair when the electric drive is switched off. It is forbidden to do this while sitting in a wheelchair. After moving the wheelchair to the desired place, remember to switch the electric drive back on.

## Safety

### Safety belt

**WARNING!**

Remember to fasten seat belt before each ride.

Care should be taken to ensure that no objects are jammed under the belt. This way you will avoid painful pressure points.

The hip seat belt is not part of the fastening system of an electric wheelchair and/or its user when transported in a motor vehicle.

The hip seat belt is used to stabilize the sitting position and prevent the user from falling out of the electric wheelchair forward.

To install the hip seat belt, pull both ends forward and close the latch so that it locks audibly. To open the hip seat belt, it is necessary to press the buckle buttons and open the belt ends.

Illustration 26. Properly fastened safety seat belt.

**ATTENTION!**

The length of the hip belt can be adjusted, making sure not to tighten it too much.

## Anti-tip wheels

- To set the extension of the anti-tipping wheels, pull the lock **A**
- Then slide out the wheel mounting rail and adjust the holes **B** so that the lock can be reinserted (Illustration 27.).


Illustration 27. Anti-tip wheels.

## Transporting the wheelchair

The wheelchair should only be transported by vehicles adapted for this purpose. The wheels should be blocked. It is forbidden to be in a wheelchair during transport, this may cause a fall and serious injuries. Always make sure that the wheelchair is properly fixed in the transport vehicle, otherwise it may cause injuries to people moving in this vehicle as well as serious damage to the transport vehicle itself.

When driving up the ramp into the transport vehicle, protection should be provided at the rear, to protect the wheelchair from tipping over.

All additional elements of the wheelchair during transport should be removed from it or attached to the wheelchair so that they do not cause damage to people or transport equipment during sudden braking.

For transport, disconnect the battery power supply in reverse order as shown in Figure 7. in the “Assembly of elements” section.

To fold the wheelchair, follow these steps:

1. Disassemble the battery in the reverse order shown in illustration 28.
2. Remove the seat cushion.
3. Then carefully lift the scissor hinge blocking the cross frame upwards, taking care not to crush your hands and fingers and pull the middle part of the seat upholstery upwards.
4. Fold the backrest as shown in illustrations 4 and 5.
5. Remove the footrests.
6. Remove the control panel by unscrewing the adjustment screw as shown in illustration 12.

Illustration 28. Folding the wheelchair.

## Inflating tires

- Unscrew and remove the plastic cap from the valve and put it in a safe place. It is important not to lose it, because the valves will get dirty, which will lead to difficulties in pumping and faster air loss.
- Apply the pump nozzle to the valve and check for proper adhesion. If necessary, it is advisable to press the pump against the valve so that air does not escape to the sides.
- If the pump has a lever with a clamp at the nozzle, clamp it. If the nozzle of the pump and the valve have a threaded tip, they should be screwed together. As a result, air will flow straight into the inner tube and the nozzle will not detach from the valve.
- When inflating, check the tyre pressure periodically.
- Once pumping is complete, remove or unscrew the nozzle from the valve by applying pressure to its outlet with your thumb. It is important to do this efficiently so that not too much air escapes.
- Screw the plastic nut onto the valve.

Illustration 29. "Schrader" valve  
(car valve)

### ATTENTION!

Observe the correct wheel inflation in accordance with the values given by the manufacturer on the tyre sides.

Use a pump with a pressure gauge or a bicycle hand pump.

## Use

### Getting in the wheelchair

#### WARNING!

Entering or exiting the electric wheelchair can only be performed with the wheelchair switched off and the drive disengagement clutch lever set in the driving mode (lever in the upper position).

Unintentional movement of the joystick (driving and steering lever) can lead to uncontrolled movement of the electric wheelchair! – Danger of accident!

## Wheelchair driving training

1. Handling the wheelchair and performing all manoeuvres must be well mastered in order to move safely. First, you need to master moving forward. Ensure that the minimum speed is set.
2. When you have mastered moving straight forward, try to drive by circling the shape of the letter S with a wheelchair.
3. Then you can start backwards training. Regardless of the programmed forward speed, the wheelchair will always move backwards more slowly.

## Driving on uneven terrain

Navigating rough terrain poses more risks. It is important to follow the safety rules. Failure to do so, may result in serious injuries.

- Do not climb hills above 8°degrees of inclination. Note that uneven or wet surface will further reduce grip and weaken brake performance when ascending and descending a hill.
- Do not go backwards from the hills.
- Do not move on uneven, unpaved surfaces (sand, gravel, grass).
- When descending a hill, move at a minimum speed.
- Do not get off the wheelchair or sit on the wheelchair standing on an inclined surface (always place the wheelchair on a level, horizontal surface first).
- When overcoming hills, always approach in a straight line in accordance with the slope of the terrain.

## Getting out of the wheelchair – end of product use

### WARNING!

Never get in or out of the wheelchair without first turning off the power. This prevents the wheelchair from moving by accidentally switching the joystick of the controller.

1. Stop the wheelchair completely.
2. Turn off the power supply using the control panel (Illustration 3.).
3. Tilt the armrest on the side that will allow the wheelchair to be abandoned safely.
4. Carefully and safely dismount from the seat and take a safe position next to it.
5. To make it easier to get in the wheelchair, you can leave the armrest tilted.

## Batteries

The wheelchair is equipped with long-life batteries. The batteries are protected and require basic maintenance and regular charging.

- Perform a full charge cycle of at least 12 hours before using the wheelchair for the first time.
- To maintain the best performance of the battery always perform a full charge cycle.

**ATTENTION!**


Failure to fully charge the battery shortens its life.

### Safety of battery charging

- Use a standard 3-pin grounded plug and a grounded wall socket.
- Do not use extension cords or adapters to avoid the risk of fire, electric shock and injury.
- Do not use the same wall socket for other activities when charging the battery.
- Only use the factory charger designed for this wheelchair to avoid damage.
- Do not leave the charger in the socket if the batteries are not charging.
- Keep the charger and batteries away from water and moisture.
- Charge the wheelchair in a well-ventilated area.
- Do not cover the charger fan to allow heat to escape.
- Leave the unit idle until the batteries are fully charged. Do not sit on the wheelchair while charging the batteries.
- Do not stop charging the batteries until the charge indicates that the batteries are fully charged.
- Do not place the charger directly on the battery.
- Check that the parameters of the charger correspond to the parameters of the batteries and the voltage in the electric socket. Connecting a battery with other parameters can lead to its damage.
- Especially during charging, avoid exposing the wheelchair and charger to high temperatures and smoke to avoid ignition and explosion.
- To avoid electric shock, it is forbidden to touch the charger or battery with wet hands or other objects and do not expose to moisture.
- Check that the charger switch is set to the "O" position before connecting the charger to the electric socket and the batteries connected to it.
- Do not expose the charger to rain, water or moisture.
- The charging process should take place at an ambient temperature of 0°C to 40°C, in other conditions charging may fail or lead to damage the batteries.
- If the charging process indicator does not light up, disconnect the cables for 10 minutes and try connecting them again after this time. If this does not help, contact the distributor.
- Contact the distributor if the charge level indicator remains yellow after charging for more than 12 hours.

## Charging process

1. Park the wheelchair as close to the grounded socket as possible.
2. Switch off the wheelchair.
3. Make sure that the wheelchair is powered on (the clutch is in the right position – travel mode).
4. Connect the charger to the charging socket in the wheelchair.
5. Connect the charger to a grounded socket, remember to check that the socket has the correct voltage of 230V.
6. Switch on the charger.
7. The red LED light on the charger indicates an active power supply.
8. The yellow LED light on the charger indicates an active charging process.
9. The green LED light on the charger indicates that the charging process is complete.
10. The optimal charging time is 12 hours.
11. When charging is complete, disconnect the charger from the wall socket first, and then remove the plug from the charging socket.

The battery can be charged outside the wheelchair.

- remove the battery and charge it in a place free of moisture by connecting the charger to the socket in the battery (Illustration 27.)


Illustration 30. Socket in the battery.

### WARNING!

Do not leave the charger connected to the charging port if the charger is not actively charging the wheelchair.

The charger is for indoor use only.

- Do not expose it to rain or water.
- Do not cover the charging terminal of the battery.
- Always keep the battery dry, never use it high full power and never charge it for more than 12 hours.

## Batteries - Storage

- If the wheelchair is to be unused for a long time, remove the batteries and store them in a dry and cool place.
- Avoid completely discharging the batteries.
- If you use the wheelchair every day, charge the batteries every day.
- If you use the wheelchair once a week, charge the batteries once a week.
- Store the batteries fully charged.
- If the wheelchair is not used regularly, charge it at least once every two weeks in a full charge cycle to maintain battery life.
- Never replace the batteries with batteries other than those supplied with the wheelchair.
- If the batteries are damaged or need to be replaced, contact the distributor.
- Do not dispose of batteries in the waste bin. Dispose of in accordance with local waste management regulations.
- To prevent short-circuiting and damaging the battery, never touch their terminals with metal objects.

## Service and repair

Inspection should be performed regularly, and should be preceded by a thorough cleaning of the product (so as to reveal any hidden damage). If repairs are needed or if you have any questions about the product, contact your distributor.

## Modification

Do not modify, remove, disable or add any parts, components or functions of the electric wheelchair that are not supplied by the manufacturer.

## Care and maintenance

### ATTENTION!

Insufficient or abandoned care and maintenance of the electric wheelchair leads to limitation of the manufacturer's liability.

## Cleaning

In order to ensure long-term and efficient operation of the wheelchair, it should be subjected to basic maintenance and inspections carried out by the service point. Inspections should take place once a year. If you have any doubts or damage, contact the service center, which has all the competences to solve the problem.

**WARNING!**

Observe the safety instructions given by the manufacturer of cleaning and disinfecting products. Failure to do so may result in skin irritation or premature damage to the upholstery of the electric wheelchair.

**It is not permitted to use:**

- solvents,
- toilet cleaners,
- agents with chlorine content,
- sharp brushes and hard objects,
- corrosive cleaners.

*Do not clean the wheelchair with pressure or steam units!*

**Disinfection**

If the product is used by many people (e.g. in a social care home), you should use commercially available disinfectants. Before proceeding with disinfection, clean the surface. Only biocidal products authorised in Poland by the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products or by another competent authority of a given EU Member State may be used for disinfection.

**WARNING!**

Disinfectants may damage the outer surface, reducing the long-term functionality of the components. Always observe the instructions concerning the use of the products specified by the manufacturer.

**Reuse (reusable product)**

The product is suitable for reuse. The number of times it can be reused depends on how often and in which way the product is used. Before reuse, thoroughly clean and disinfect the product and check its safety features. If necessary, take repair action. If necessary, take repair action.

**Product life cycle**

The expected life of the medical device is 4 years if it is used as intended and if all maintenance and repair activities are performed as required. The life cycle of the medical device also depends on the frequency of use, the operating environment and care of the device. The life cycle of the device can be extended using spare parts. Spare parts are normally available for 5 years after the model is withdrawn from production. The expected life cycle of the medical device is not equivalent to an additional guarantee.

## Disposal of the electric wheelchair

The electric wheelchair must be disposed of in accordance with current local and national regulations. For proper disposal of packaging, metal components, plastic components, neoprene, silicone and polyurethane materials, contact your local waste disposal company or authorised distributor.

## Maintenance plan

When	Scope of control
<b>Before the ride</b>	<p>In general, check for impeccable performance. After turning off the power, check the joystick. Ensure that it is not bent or damaged and that it returns to the neutral position when released. Make sure that the rubber cover around the base of the joystick is not damaged. Check the cover visually. Check the controller cable visually. Ensure that it is not ragged, cut or exposed wires. Check that there are no flat places on solid tires that may have a negative impact on the stability of the wheelchair.</p> <p>In case of problems and questions, contact the distributor.</p>
<b>Weekly</b>	<p>Make sure that all parts of the control panel (joystick) are well attached to the handle.</p> <p>Check the tightening of the screws and tighten with moderate force if necessary.</p> <p>Check the correct tyre pressure. If the tire does not hold air, contact the distributor.</p> <p>Check the operation of the electromagnetic brake. This test must be carried out on a level surface with at least a meter of free space around the electric wheelchair. Try to push the wheelchair with the motor clutches engaged into driving mode. In case of doubt, always contact the distributor.</p>

<p><b>Monthly</b> Checking the brakes</p>	<p>Turn on the controller and reduce the speed level of the electric wheelchair. After one second, check the battery level indications - make sure that the LED light remains on. Slowly push the joystick forward until you hear the electromagnetic brakes click - release the joystick immediately. The operation of each electromagnetic brake must be audible within a few seconds of the movement of the joystick. Repeat this test three times by pushing the joystick backwards, then left and then right.</p>
<p><b>Monthly</b> (depending on the length of the traveled sections)</p>	<p><b>Inspection of tire profile</b> <b>Check drive tyres for wear:</b> Rear wheels - minimum tread depth = 1mm. Replacement of tires should be commissioned to an authorized service. <b>Front wheels:</b> Check front, self-aligning wheels for wear. If necessary, contact your distributor. <b>Anti-tipping wheels:</b> Check that the anti-tipping wheels do not rub against the ground when using the electric wheelchair. If necessary, adjust them. Check that the anti-tipping wheels are not excessively worn. In case of doubt, always contact the distributor. Keep the electric wheelchair clean and protected from foreign bodies such as mud, dirt, hair, food, beverages, etc.</p>
<p><b>Every 6 months</b> (depending on the frequency of use)</p>	<p><b>Check</b></p> <ul style="list-style-type: none"> <li>- Cleanliness</li> <li>- Overall condition</li> </ul>
<p>Manufacturer's recommendations <b>Every 12 months</b> (depending on the frequency of use)</p>	<p><b>Maintenance work in the manufacturer's service department</b></p> <ul style="list-style-type: none"> <li>- Vehicle</li> <li>- Charger</li> </ul>

## Replacing motor brushes

### ATTENTION!

If the following symptoms are observed when using the wheelchair:

- loss of power when starting or driving under load
- sparking visible through the engine vents
- unusual noises such as crackling, knocking or vibrations
- burning smell or smoke from the engine compartment

Stop using the device immediately and contact an authorised distributor if you notice any of the above symptoms. Do not attempt to repair the device yourself to avoid the risk of electric shock, fire or loss of control of the wheelchair.

## Fault codes and troubleshooting

The table below describes the error indication and some possible actions that can be taken to correct the problem. The listed activities are not ordered in any specific order and are only suggestions. The intention is that one of the suggestions can help to solve the problem. In case of doubt, always consult an authorized supplier.

No.	Error description	Error analysis	Possible solution
1	The wheelchair does not work.	<ol style="list-style-type: none"><li>1. The vehicle is not powered on.</li><li>2. The connection between the battery and the controller is loose.</li><li>3. The battery is empty.</li></ol>	<ol style="list-style-type: none"><li>4. Switch on the unit.</li><li>5. Check the connection elements and reconnect them.</li><li>6. Recharge battery</li></ol>


2	Does not work when the battery is charged.	<p>7. The charger is defective</p> <p>8. The battery is damaged</p> <p>9. The controller is damaged</p>	<p>Disconnect the battery from the controller and connect the charger to the controller - if it can be switched on, the charger and the controller are in good condition. Battery may be damaged</p> <p>Contact the distributor.</p>
3	The acoustic and light signals flicker once.	Low battery voltage.	<p>Charge battery immediately. If the red light is still on when the battery is charged, disconnect the battery and connect the controller to the charger to check that the charger is working properly. The battery may be damaged - always contact the distributor.</p>
4	The acoustic and light signals flicker twice.	Engine disconnected.	Reconnect the plug.
5	The acoustic and light signals flicker 3 times.	Left motor brake fault	Reconnect the plug.
6	The acoustic and light signals flicker 4 times.	Right engine disconnected.	Reconnect the plug.
7	The acoustic and light signals flicker 5 times.	Right motor brake fault.	Reconnect the plug.
8	The acoustic and light signals flicker 6 times.	Possible overcurrent of the controller.	<p>1. Restart the control panel.</p> <p>2. Disconnect the control panel and reconnect - contact the distributor.</p>
9	The acoustic and light signals flicker 7 times.	Control panel malfunction possible.	Place the control panel (joystick) in the neutral position. Check the connection of the joystick.


10	The acoustic and light signals flicker 8 times.	Possible malfunction of the controller.	Contact the distributor.
11	The acoustic and light signals flicker 9 times.	Possible malfunction of the controller.	Contact the distributor.
12	When charging, the controller is very hot.	Possible malfunction of the controller.	Contact the distributor.
13	The green LED light of the charger is active but the charger is not charging the battery.	Possible malfunction of charger or battery.	Check output voltage of charger and battery. The normal voltage of the charger should be in the range of 26.4~ 28.4V, and the minimum battery voltage 21V~25V. Contact the distributor.
14	The charger plug is loose or damaged.	Charger fault.	Contact the distributor.
15	Unable to charge battery.	Possible malfunction of charger or battery.	Check output voltage of charger and battery. The normal voltage of the charger should be in the range of 26.4V~28.4V, and the minimum battery voltage 21V~25V. Contact the distributor.
16	The motor is running loudly.	Motor damage.	Contact the distributor.
17	The drive disengagement clutch is locked.	Clutch friction is too high.	Contact the distributor.
18	It is difficult to push the wheelchair manually.	Motor drive clutch engaged.	Disconnect the clutch to manual pushing mode. See page 63.

19	The motor works, but the wheelchair does not move.	Motor drive clutch disengaged.	Switch the clutch to driving mode. See page 63.
20	Unable to disengage drive clutch.	Drive clutch fault.	Contact the distributor. Motor drive clutch must be replaced.
21	The footrest does not lock in the set position.	Footrest handle fault.	Contact the distributor.

## Symbols

The symbols used below are used to identify warnings, mandatory actions and prohibited actions. You should read everything and understand their meaning.

	Medical device
	Date of manufacture
	Manufacturer
	Authorized representative in the European Union
	Distributor
	Conformité Européene (compliance with the European Union)
	Catalog number (to re-order)
	Serial number
	Batch code

	Unique Device Identifier
	Please refer to the instructions for use
	Please refer to the instructions for use
	Warning
	Protect from moisture
	Possibility of jamming fingers
	Permitted temperature
	For indoor/outdoor use
	Maximum user weight
	Total permissible mass
	Maximum safe inclination
	The product is designed to be transported by car
	The product is not intended to be used as a seat for transport in a motor vehicle

#### SELECTIVE COLLECTION OF ELECTRICAL AND ELECTRONIC EQUIPMENT.

This symbol on the products or accompanying documents indicates that used electrical or electronic products and batteries must not be disposed of in municipal waste. For proper disposal, restoration or recycling, such products should be handed over in collective places for this type of waste. Alternatively, in some countries of the European Union or other European countries, you can return your products to a local seller at the time of purchase of a similar new product. By correctly disposing of this product, you will help preserve valuable natural resources and support the prevention of potential negative impacts on the environment and human health, which could result from incorrect waste disposal. More information is available from the communal offices or waste collection sites. In case of incorrect disposal of this product, penalties may be imposed in accordance with local regulations. If you want to dispose of this device, please obtain the necessary information on the correct disposal at your local authorities or from your dealer.

**Liquidation in countries outside the European Union:** This symbol is valid in the European Union. If you want to dispose of this device, please obtain the necessary information on the correct disposal at your local authorities or from your dealer.

## Contact Information

### Manufacturer

mdh Sp. z o.o.

Address: Maratońska 104, 94-007 Łódź, Poland

tel. +48 42 674 83 84, fax. +48 42 636 52 21

[www.mdh.pl](http://www.mdh.pl)

[www.viteacare.com](http://www.viteacare.com)

We reserve the right to make technical changes due to technical progress.