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English version on next page


Table of contents

Introduction from mdh Sp. z o.o.....	37
Indications and contraindications	38
Indications.....	38
Contraindications	38
Usage	38
Safety instructions	38
Information on electromagnetic interference - EMI	41
Safety in handling the battery	42
Description and installation for models Q1-10 i Q2-16.....	43
Description of the components for model Q1-10	43
Description of the model components Q2-16.....	44
Specifications of the drives Q1-10 i Q2-16	45
Assembly and adjustment.....	46
Battery and charging.....	51
Description and installation - Q5 model.....	53
Description of the Q5 model components	53
Technical parameters - model Q5.....	54
Assembly and regulation.....	54
Battery and charging.....	56
Technical parameters of the manual wheelchair	58
Care and maintenance	59
Tires.....	59
Wheel pressure and braking force.....	60
Cleaning.....	60
Disinfection.....	60
Troubleshooting.....	61
Technical review.....	62
Storage	63
Disposal of electric drive and wheelchair.....	63
Warranty.....	64
Symbols	65
Contact details	67

Introduction from mdh Sp. z o.o.

This product has been manufactured by skilled and dedicated personnel. It has been designed and manufactured to the highest quality standards. Before using for the first time, please read these operating instructions. Improper use may result in personal injury or traffic accidents. If you still have questions after reading this manual, please contact your specialist dealer.

ATTENTION

Always follow the basic precautions according to the information given below regarding safety when using the product. Read the entire manual carefully before use. Keep this manual for future reference.

WARNING!

Improper use can lead to death, serious injury, unsafe conditions or product malfunction.

ATTENTION!

Improper use may lead to personal injury and/or damage to the product.

SUGGESTIONS!

Follow the instructions below to keep the product in good working order.

Indications and contraindications

Indications

The product is intended for disabled people with musculoskeletal conditions to increase their mobility.

Contraindications

Contraindications to upright standing. Injury or paresis of the upper limbs preventing the use of a power wheelchair.

Usage

The product is packaged in the assembled position. After removing from the package, check that it contains all parts according to the list of components above. Make sure that all parts have no visible mechanical damage. If in doubt, contact the manufacturer. Do not use the wheelchair until all components are installed and the battery is charged. Use the wheelchair and its drive on a level and stable surface, indoors or outdoors, under the following environmental parameters:

Temperature: 50C - 400C

Humidity: 30% - 70%

ATTENTION!

The product is intended for people with disabilities with musculoskeletal disorders. If you are not sure whether your health condition allows you to use the product consult a health care professional.

Safety instructions

- a. Before using the wheelchair, please read this manual and keep it for future reference..
- b. Apply the brakes when getting on and off the wheelchair!
- c. Before sitting down or getting out of the wheelchair, do not stand on the footrest plate, before that, fold or completely tilt the footrests aside.
- d. Patients requiring additional head stabilization should move the wheelchair ONLY under the supervision of an experienced attendant.
- e. It is unacceptable to move the wheelchair with the user.
- f. When moving the wheelchair through doorways and other passageways, make sure there is enough space on both sides to avoid injuring your hands or damaging the wheelchair.
- g. When using public roads, follow traffic regulations.
- h. As with other vehicles, do not drive the wheelchair under the influence of alcohol or other drugs or intoxicants. This also applies to indoor driving.

- i. The wheelchair is not designed to transport the patient up stairs, as a stroller or by pulling the patient on it. The structure of the wheelchair and the connecting elements are not designed for such heavy loads if it is used in this way, it risks weakening or damaging the entire load-bearing structure of the wheelchair and a possible accident!
- j. The wheelchair is designed for indoor and outdoor use only on solid and level surfaces with a slope not exceeding 8°.
- k. Use of the wheelchair on a surface that does not provide secure support, i.e. wet, icy, greasy, heavily dusty, is unacceptable due to the risk of falling.
- l. For overcoming difficult architectural barriers (steep or long ramps, high steps, etc.), the assistance of a caregiver is required.
- m. If the user does not have enough experience in using the wheelchair, ask a caregiver to help him/her overcome obstacles.
- n. It is forbidden to lean outside the wheelchair while carrying people. There is a risk of accident.
- o. Perform a visual inspection of all structural components of the product including the stability of the product and its connections before each use.
- p. Do not use the wheelchair if there is any visible or detectable damage (e.g., a crack in the body), as there is a risk of an accident.
- q. Wheelchair is used to transport one user, it is forbidden to stand on the wheelchair or stand with the whole body weight on the footrests.
- r. Do not leave the wheelchair with the user in the direct contact with an open flame source.
- s. Note that the surface temperature of the handles and frame components may increase after exposure to the sun or other heat source, while in cold temperatures there is a danger of very high cooling of the wheelchair's metal parts. This can cause burns or frostbite on body parts.
- t. The wheelchair should be moved, after unhooking the drive, by grasping the seat frame structure. Never move the wheelchair by the armrests or by the footrests.
- u. Wheelchair can not be used for use as a seat in a motor vehicle (bus, train).
- v. The wheelchair is designed for users with a maximum weight of 100 kg.
- w. The wheelchair should always be easy to fold and unfold. If you have problems with folding, read the instructions again, or contact your dealer.
- x. When moving on sloping surfaces, special care must be taken. The maximum permissible inclination angles of the ground up and down are 8°(14%).
- y. Brakes are not used to slow down the wheelchair as it moves, they should only be used to prevent unwanted movement of the wheelchair.
- z. Both brakes should be properly adjusted. Improper adjustment of one, resulting in the locking of only one rear wheel, which on steep inclines can cause the wheelchair to rotate around the locked wheel, which is difficult to control.
- aa. Tire pressure should be checked regularly. If the tire pressure is too low than the recommended pressure, the effectiveness of the brakes may be significantly reduced!
- ab. All repairs should be performed by qualified personnel of the manufacturer.
- ac. Hold the wheelchair with both hands at the same time.

- ad. When picking up objects from the floor, do not reach backwards as you risk falling backwards.
- ae. Do not reach objects from the floor in front if it requires you to get outside the wheelchair.
- af. Check that the locking mechanism used to fold and unfold the wheelchair is working properly.
- ag. Any activity such as getting on, getting off, changing the patient's position, or stopping the wheelchair requires locking the running wheel brakes.
- ah. Do not leave the patient in the wheelchair (even with the brakes applied) on exits, ramps, mezzanines, etc.
- ai. Each time after installing the wheels, check whether the wheel joint axle lock has worked, when installing should be heard the characteristic click of the locking pin, in addition, strongly pull the wheel without pressing the button in the wheel axle.
- aj. Using the parking brake while driving can risk an accident.
- ak. Always after each adjustment, check the correct functionality of the brake.
- al. If wheel tires are contaminated by oil, grease, water or other agents, braking efficiency may be reduced.
- am. The manufacturer prohibits independent modifications to the wheelchair especially the installation of additional wheels. This may result in the risk of falling and serious injuries. The manufacturer is not responsible for such actions.
- an. After long use of the wheelchair, the braking force may weaken due to the gap created between the brake pressure pin and the wheel tire. It is possible to adjust the brakes yourself. To do this, loosen the screws locking the position of the brake mechanism, move the entire mechanism towards the wheel to such a position that when the brake is applied again, the wheel does not rotate (test the correctness of the adjustment on the wheelchair with the patient). After completing the adjustment, again secure the position of the mechanism with the locking screws.
- ao. It is unacceptable to use the wheelchair if any visible or noticeable damage has been found, or the degree of wear and tear of the wheelchair's parts indicates a potential risk of an accident (e.g. wear of the running wheels, cracks in the upholstery, etc.). All repairs should be entrusted to the manufacturer's service department.
- ap. For wheelchairs equipped with inflatable wheels, do not exceed the maximum allowable wheel inflation pressure for the value listed on the wheel tire.
- aq. There is a possibility of jammed fingers during the installation of the drive to the wheelchair and in the area of unfolding and folding of the moving parts of the wheelchair frame, footrests, armrests and brakes.
- ar. Keep out of reach of children under five.
- as. When transporting (e.g. in the trunk), do not load the wheelchair and drive with any weight, it may affect its structure.
- at. The product is not intended for use by a person with vision, reading or cognitive disabilities.
- au. The wheelchair is not suitable for driving up and down stairs and escalators. Always use the elevator.

- av. When traveling with a wheelchair, you may need to use a lift/elevator to assist in transportation. It is recommended that you carefully read the instructions, specifications and safety information set forth by the elevator/lift manufacturer before using the product.
- aw. Do not modify the drive in any way. This will lead to loss of warranty.

WARNING!

Any “serious incident” related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides.

Information on electromagnetic interference - EMI

Studies of electromagnetic interference (EMI) have shown that it adversely affects the operation and control of electrically powered mobile devices.

EMI can come from a variety of sources, such as cell phones, two-way radios, radio stations, TV stations, amateur radio (HAM) transmitters, wireless computer links, microwave signals, paging transmitters and transmitters used by emergency vehicles. EMI waves can cause unintended movement of the wheelchair or damage to the controller. Every electrically powered mobile device has a certain immunity to electromagnetic interference. The higher the immunity level, the greater the protection. The intensity of interference can be measured in volts per meter (V/m). Current technology allows for an immunity level of 20 V/m, which provides protection against the most common sources of electromagnetic interference. As you get closer to the source, electromagnetic interference becomes more intense. Special attention should be paid to electromagnetic interference from hand-held radio devices. It is possible to unintentionally bring high levels of EMI close to the front drive control system, which can affect movement and braking. To prevent possible interference with the front drive control system, follow the warnings below:

WARNING!

Do not use hand-held radio transmitters such as a CB-radio or cell phone near a front drive control system.

WARNING!

Pay attention and stay away from nearby transmitters, such as radio or TV stations.

Safety in handling the battery

Positive and negative terminals of the battery must not be connected to conductive objects, such as metal, as they may risk a short circuit or electric shock.

Stay away from flammable materials.

Do not go near or be exposed to heat sources such as open flames or sparks.

There are chemicals in the battery case. It is strictly forbidden to disassemble the battery case.

It is strictly forbidden to short-circuit the battery or throw the battery into a fire. Explosion hazard.

Description and installation for models Q1-10 i Q2-16

Description of the components for model Q1-10

Part for connecting to rigid frame

Part for connecting to cross frame


Figure 3

Description of the model components Q2-16

Ilustracja 4

Ilustracja 5


Figure 6

Turning on the front light:


- hold the speed reduction button for 3 seconds.

Specifications of the drives Q1-10 i Q2-16

Model	Q1-10	Q2-16
Maximum user weight	100 kg	100 kg
Diameter of the front wheel	10" (25,4 cm)	16" (40,64cm)
Maksymalna prędkość	20 km/h	28 km/h
Motor power	500 W	400 W
Battery	12 Ah / 36 V - litowy	16 Ah / 36 V - litowy
Energy range	30±5 km	30±5 km
Weight	14 kg	22 kg

Assembly and adjustment

Folding the steering wheel

To fold the steering handles, slide the spring-lock rings on the left side to the left and on the right side to the right.

Figure 7

After assembly

Steering wheel height adjustment


1. Loosen the allen screw in the handlebar bracket as shown in Figure 8.
2. Either slide the handlebar in or out to the desired height.
3. Tighten the screw.

Ilustracja 8


Steering wheel angle adjustment

1. Loosen the allen screw in the steering wheel bracket as shown in Figure 9.
2. Either slide the steering wheel in or out to the desired height.
3. Tighten the screw.

Ilustracja 9

Assembly of a connector to the cross frame of a wheelchair

1. Install the mounting brackets on the wheelchair's lower frame at the front of the wheelchair. Figure 10.
2. Tighten the screw in the handle.
3. Install the telescopic handle in the mounting bracket.
4. Secure the connector in the holder with the ratchet lock and tighten the lock screws.
5. Position the telescopic connector so that the front part of it coincides with the center axis of the wheelchair.
6. Tighten the locking screws.
7. Install the depth connector securing it with the ratchet lock and locking screw.
8. Adjust the depth of the extension, bearing in mind the distance of the footrests from the drive module.
9. Add a connecting element.
10. Adjust the height of the connecting element so that there is free access of the drive module to the wheelchair.
11. Connect the drive module to the finished mounting system.

WARNING: Use extreme caution during installation, lack of attention may risk body parts being crushed by moving and assembling parts.

12. Tighten the locking screw in the connector.
13. Fold the support foot.

WARNING: Before driving, check that the rear wheels are clipped in correctly on the quick couplers.

The connector is located on both sides of the wheelchair frame.

Step 1: Adjust the “width adjustment lever”.

Step 2: Place the connector in the correct position, then adjust and install it.

Step 3: Connect to the module and raise the lever upwards.

Disassembly of the drive module

1. Unfold the drive support foot.
2. Unscrew the locking screw of the connector, then press down firmly and detach the drive module.
3. Disassembly of the drive module to the wheelchair should be done in the reverse order of the connector assembly procedure.

WARNING: Tighten the lock nut before use.

WARNING: Never forget to fasten seat belts before starting up.

WARNING: Do not turn violently at high speed.
This will prevent the wheelchair from tipping over.

Mounting the drive to a rigid wheelchair frame

1. After locking the brakes in the wheelchair, attach the drive module.
2. Unlock the safety pin and raise the rotation lock.
3. Adjust the angle of the connecting element to fit the front frame of the wheelchair.
4. Install the connecting elements on the front frame of the wheelchair without fully tightening the locking device (make sure that the connecting elements are at the same height, and the axis of the wheelchair and the drive module coincide).
5. Completely tighten the locks of the connecting elements.
6. Release the swivel lock, while holding the drive module by the steering wheel, tilt it forward to the right angle, so that the front wheels are not in contact with the ground (at least 5 cm).
7. Make sure that the locking mechanism of the ratchet locks the connector and tighten the locking pin.
8. Assemble the support foot.

WARNING: Before driving, check that the rear wheels are clipped in correctly on the quick-release couplers.

Step 1: "Connecting element" is placed on both sides of the wheelchair frame.


Step 2: Connect to the module and raise the lever.

Mounting the drive module to the footrest step


WARNING: The connector can only be used for a fixed and one-piece footrest step.

Adjustment by screw

1. Adjust the connector to the depth of the footrest step.

Angle adjustment
screws for the connector

Fastener mounting
screws for the drive

2. Tighten the connector to the drive module.

3. Adjust the angle of the connector.

Battery and charging

WARNING: Turn off power to the drive before removing the battery.

WARNING: Before recharging, first install the battery and then connect the charger.

When the charging cycle is complete, first unplug the charger from the power source and then unplug the battery pack from the battery jack on the drive!

Assembly and disassembly of the battery

1. Turn off the drive.
2. Turn the key to release the battery lock.
3. Tilt and pull out the battery pack in the direction from bottom to top.

To view the battery charge status, hold down the button on the case.

Battery switch

Charging port

Figure 11

Charging the battery

Note that the first recharging of the batteries should be done when they are fully depleted and requires more than 12 hours.

1. Turn off the battery switch. (Figure 11)
Remove the battery from the module.
2. Plug the charger's power cord into an outlet.
3. Uncoil the power supply cable and plug it into the outlet. It is recommended to charge the battery pack for 4-6 hours.
4. When the battery is fully charged, unplug the power adapter plug from the outlet, and then disconnect it from the charger.

Safety when charging

- a. Before using the drive for the first time, charge the battery for a minimum of 12 hours.
- b. It is forbidden to use pliers or wires to directly connect the terminals of the battery to charge it.
- c. It is forbidden to use the drive while the battery is charging. Do not use non-standard power sources (such as generators or inverters), even if their voltage and frequency appear to meet requirements.
- d. It is forbidden to trample, squeeze, bend, tug or bundle the power cord, especially at the charging port.
- e. Pay special attention that children and pets do not touch the power cord during charging. When disconnecting the power cord, hold the plug and carefully pull it out.
- f. To keep the long life of the battery, it should be completely discharged and then proceed to recharge to maintain the correct cycle.
- g. Be sure to check that the power supply is operational before use. After the second and each subsequent discharge of the battery, the recommended charging time is 4-6 hours.
- h. Be sure to charge the battery frequently. Fully charge the battery after each use of the drive. If the battery is discharged for a long time, it may be permanently damaged.
- i. If you use a battery with the wrong specification, voltage and capacity, the drive may be damaged and may not achieve proper performance.
- j. In case of a damaged or cracked battery, immediately place it in a plastic bag and contact the distributor for further disposal instructions. The distributor is obliged to provide the necessary information on the recommended recycling of the battery pack.

Description and installation - Q5 model

Description of the Q5 model components

Connecting element
for rigid frame -
connection to legrest

Connecting element
for cross frame

Connecting element
for rigid frame


Technical parameters - model Q5

Model	Q5
Maximum user weight	100 kg
Front wheel diameter	8"(20 cm)
Maximum speed	30 km/h
Motor	300 W
Battery	5,8 Ah / 36 V - litowy
Energy range	20±5 km
Weight	8 kg

Assembly and regulation

Assembly of steering wheel

Install the steering wheel according to the instructions in figure 12.

Ilustracja 12

Steering wheel height adjustment

Open the "lock", adjust the handle to the desired height and then close the "lock".
Figure 13.



Figure 13

Switching the drive on and off:

1. To turn on the drive, press and hold the multifunction button for 2.2 seconds until the indicator light turns green (after a short press, the light turns purple, to turn on press longer until the light turns green).
3. To turn off the drive, press and hold the multifunction button longer.

ATTENTION.

When starting the Q5 drive, release the accelerator lever and hold down the brake lever.

The Q5 has a 3-speed transmission, 1st gear is the slowest, 3rd gear is the fastest. When the drive is off, the indicator light is also off.

When activated, the indicator light indicates the operating status by illuminating the following colors:

1-gear - green

2-gear - blue

3-gear - yellow.

When a failure occurs, the light glows red. When the battery charge is low, the indicator light blinks red.

If cruise control is activated in the drive, the indicator light blinks in the color corresponding to the color of the gear in use.

Shifting gears: After startup, press the multifunction button which cycles "1st gear, 2nd gear, 3rd gear".

Acceleration and braking

Brake:

The Q5 drive uses an electronic brake. When you press the brake lever harder, you will get more braking power. Don't accelerate when braking.

Cruise control:

After driving at a constant speed for 6 seconds, the drive module will be in cruise control mode and move at a constant speed (no need to press the accelerator lever). When increasing speed, slowing down, or changing gear, the module will come out of cruise control mode.

Battery and charging

Warning:

Turn off the power of the drive before removing the battery.

Warning: before recharging, first install the battery and then connect the charger. When the charging cycle is complete, first unplug the charger from the power source and then unplug the battery pack from the battery jack on the drive!

Assembly and disassembly of the battery


Figure 14


Figure 15

1. Press the button, as shown in figure 14.
2. At the same time, lift the battery and then remove it from the drive (Figure 15).

Charging the battery

1. Turn off the power to the drive with the multifunction switch.
2. Remove the battery from the drive..
3. Plug the charger's power cord into a USB socket.
4. Uncoil the charger cable and plug it into a proper power outlet. It is recommended to charge the battery for 4-6 hours.
5. When the battery is fully charged, unplug it from the power outlet, and then disconnect the drive from the charger.

Please note that the battery is lithium, so use a charger designed for this type of battery. For safe, more efficient and sustainable charging, only the original charger should be used.

WARNING - To avoid the risk of electric shock, the device must be connected only to a power supply with a protective ground connection.

WARNING - Removal of ground pins may cause electrical hazards. If necessary, with access to a 2-pin socket, use a proper, approved 3-pin adapter.

WARNING - Never use an extension cord to plug in the charger. Plug the charger directly into a properly installed standard electrical outlet.

WARNING - It is forbidden, while charging the batteries, for children to play unsupervised near the wheelchair, drive and battery..

WARNING - Explosive gases may be produced when charging batteries. Keep the drive and charger away from sources that can cause ignition, such as fire, sparks, and provide adequate ventilation when charging batteries.

WARNING - Charge the drive batteries with the supplied external charger. Do not use any other type of charger.

WARNING - Before each use, check for damage to the charger, cables and connectors. If you found damage, contact your supplier.

WARNING - Do not open the battery charger housing. If the charger does not function properly, contact your authorized supplier.

WARNING - Be aware that during charging, the battery charger housing may be hot. Avoid contact with skin and other surfaces that may be damaged by heat.

Technical parameters of the manual wheelchair


Wheelchair size	16 inches	18 inches	20 inches
Maximum user weight	100 kg		
Total height of wheelchair	910 mm		
Wheelchair width	580 mm	630 mm	680 mm
Wheelchair height	1100 mm		
Seat width	400 mm	450 mm	510 mm
Seat depth	400 mm		
Seat height (from the floor)	500 mm		
Backrest height	420 mm		
Dimensions of the folded wheelchair	910 x 800 x 280 mm		
Weight of wheelchair with all components	18,5 kg	19 kg	20 kg
Weight of wheelchair without components	17 kg	17,5 kg	18,5 kg
Minimum turning radius (without drive)	105°		
Materials	Frame: steel / upholstery: nylon		
The size of rear wheels equipped with a quick-release system (quick-release)	24 inches		
Size of solid front wheels	8 inches		
Maximum safe inclination	8°		
Description of the environment in which the product will be used			
Storage and use temperature	From 5°C to 40°C		
Humidity	From 30% to 70%		
Tools needed for adjustment	allen wrench 5 mm, wrench 10 mm, 13 mm		

Care and maintenance

Tires

The tires are made of a rubber compound, which can leave hard-to-remove marks on some types of flooring (such as plastic, wooden floors, parquet, carpets and linings). The manufacturer is not responsible for damage caused by abrasion of tires or due to the following chemical processes in them.

Tire inflation - applies to wheelchairs equipped with inflatable wheels.

- Before inflating tires, identify the type of valve that is at the tire. It may be one of the two valves shown above.
- The Presta valve also has an additional small nut, which should be loosened before inflating.
- Once you have identified the valve, prepare an inflator designed for the type of valve installed on the wheelchair.
- Most new inflators are equipped with both Presta and Schrader-type valve support. Many of them allow you to reverse the tip. The smaller hole facing outward will be for Presta, the larger for Schrader. If the pump does not support Presta, you will need to equip yourself with a "Presta adapter", which is available at any bicycle store.
- At this point, if the wheel is equipped with a Presta valve, you will need to unscrew the nut located on the valve. This is a small nut that should be able to be unscrewed with your fingers. When doing this, the nut should move upward. Then the valve is considered open.
- Then remove the plastic cap from the valve and put it in a safe place. It is important not to lose it, as the valve will become dirty, leading to difficulty in pumping and faster loss of air.
- Apply the pump nozzle to the valve and check that it adheres properly. If necessary, it is advisable to press the pump to the valve so that air does not escape to the sides.
- If the pump has a lever with a clamp at the nozzle, clamp it. If the pump nozzle and the valve have a threaded end, they should be screwed together. This will ensure that the air will fall straight into the inner tube, and the nozzle will not detach from the valve.

- When proceeding to pumping, periodically check what the tire pressure is. The best way to do this is by compressing the tire with your thumbs at the end of the tread. If it is possible to make only a millimeter dent, it means that there is enough air in the tire. If necessary, continue pumping.
- When you have finished pumping, remove or unscrew the nozzle from the valve by applying pressure to its outlet with your thumb. It is important to do this smoothly enough so that not too much air escapes.

Observe the correct degree of wheel pressure in accordance with the manufacturer's figures on the tire sidewalls. Use a pump with a pressure gauge or a bicycle hand pump.

Wheel pressure and braking force

Adjust the pressure in the rear wheels to the weight of the user and then, with the assistance of the attendant, check the operation of the parking brakes (the wheelchair, loaded with the weight of the user, positioned forward in the direction of the slope, must not roll on a surface with an incline of 15%). After all, adjustments, remember to carefully check all connections.

Cleaning

- Clean the wheelchair frames with a damp cloth.
- Remove contamination on a regular basis.
- In order to maintain the operational parameters of the product, it is extremely important to keep all its elements clean and to adapt to the method of its storage.
- Storage in rooms with very high humidity (baths, laundries, etc.), or in the open air, combined with the lack of hygiene, will result in a rapid loss of functional and aesthetic parameters.
- Do not expose the wheelchair and the drive to extreme weather conditions (rain, snow, very strong sun), as some elements may change its color.

Disinfection

If the product is used by several people (e.g. in a nursing home), commercially available disinfectants should be used. Before disinfecting, clean the upholstery and handles. Only tested and approved disinfectants may be used for disinfection.

Information on recognized and approved disinfectants and disinfection methods is available from your national health authority.

WARNING!

The use of disinfectants can damage the outer surfaces, which can reduce the long-term functionality of components. Follow the instructions provided by manufacturer of the chemical product.

Not allowed to use:

- Solvents
- Toilet cleaners
- Sharp brushes and hard objects
- Agents containing chlorine
- Harsh cleaners

Troubleshooting

ATTENTION!

Every electromechanical device requires occasional troubleshooting. Many of these problems arise through batteries that are not fully charged or due to wear and tear and their inability to continue holding a charge.

What if none of the systems are responding?

1. Ensure that the batteries are fully charged.
2. Ensure that the battery is properly connected.
3. Check that the battery switch is turned on.
4. Ensure that the wiring is properly connected.
5. Remove and reinsert the key or turn the drive power off and on again.

What to do if the drive's electrical power fails?

1. Charge the module battery more often.
2. If the problem remains, take the batteries to the distributor for inspection.
3. Perform a battery load test yourself. Load meters, are available at most auto parts stores. Follow the instructions that come with the meter.

What if the battery gauge goes down and the engine doesn't run properly and has trouble starting when the accelerator lever is pushed in?

1. Fully charge the drive's battery.
2. The authorized service will carry out a load test.

For models Q1-10 and Q2-16 - the meaning of the status of errors that may appear on the display:

ERRO 1: motor failure

ERRO 2: steering wheel malfunction

ERRO 3: controller error

ERRO 4: power and brake cut-off

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ERRO 6: connection error. The device is not detected by the output of the controller

ERRO 7: connection error. The controller is not detected by the device output

Technical review

The inspection should be performed regularly. Should be preceded by a thorough cleaning of the drive and wheelchair (so as to reveal any hidden damage).

During the technical inspection you should:

Frequency	Subject of control	Remarks
Before first assembly	Inspection of the technical condition of all elements included in the package.	Check that no components have visible mechanical damage.
Every two weeks	Fastening bolts, nuts and making sure they are correct mountings.	Check that no parts have visible mechanical damage.
Every two months	Cleanliness, general condition.	In case of any mechanical damage that occurs, contact the supplier immediately.
Once a month	Wheel axles	Wheel axles should be maintained by spraying a few drops of maintenance lubricant between the axle and the hub, in addition, the wheelchair should be raised so that the wheel under test is above the ground, and the wheel should rotate freely but without clearance on the driving axle.
Once a month	Condition of threaded and soldered connections	Make attention to cracks in the frame material at the soldering points and at the locking holes of the bayonet connections).
Once a month	Condition of the tyres	Ensure that the tire tread is not worn out, and that the tires have no visible damage.

Storage

If you plan not to use the drive for a period of time then follow these rules:

- Fully charge the batteries before storage.
- Remove the batteries from the drive.
- Store the drive in a warm, dry room.
- Avoid storing in a place where the drive will be exposed to temperature changes.
- Recommended storage temperatures: from 18°C to 50°C (from 64.4°F to 122°F).

SUGGESTIONS!

Batteries that are regularly completely discharged, not regularly charged, stored at inappropriate temperatures, or stored without being charged may be permanently damaged resulting in unreliable operation and shortened battery life. To ensure proper operation, it is recommended that batteries be charged regularly, throughout the period of prolonged storage of the drive.

Disposal of electric drive and wheelchair

The electric drive must be disposed of following applicable local and national regulations. For information on the proper disposal of packaging, metal frame parts, plastic parts, electronics, batteries, neoprene, silicone, and polyurethane materials, please contact your local waste disposal company or authorized supplier.

Warranty

Wheelchair frame	12 months
Other parts (components: mechanical, electrical and electronic)	12 months
Batteries	6 months

Warranty does not include

The warranty does not cover components that may need to be replaced due to normal wear and tear during operation (tires, belts, light bulbs, upholstery, plastic bumpers, engine brushes, fuses and batteries) or have been damaged during inappropriate use or accident for which the manufacturer or supplier is not responsible. This warranty does not cover labor or service calls.

ATTENTION!

The warranty does not include deterioration of battery performance caused by being left in a fully discharged state, in cold conditions for an extended period of time, or by excessive use.

Inspection and warranty service

Warranty service can be performed by an authorized supplier. Please contact your authorized supplier for the current costs associated with a service visit.

Spare parts

The availability of replacement components and/or spare parts is subject to the capabilities of the supplier, not the manufacturer. For more information on replacement components and/or spare parts, please contact your authorized supplier.

Symbols

	Medical device
	Production date
	Manufacturer
	Country of manufacture
	Distributor
	Authorized representative in the European Community / European Union
	Conformité Européene (compliance with the European Community))
	Catalog number (to re-order)
	Serial number
	LOT code
	Unique Device Identifier
	Follow the instructions for use
	Refer to the instructions for use
	Warning

	Type BF application part
	Protect from humidity
	Keep away from sunlight
	Possibility of jamming fingers
	Permitted temperature
	For indoor/outdoor use
	Maximum user weight
	Maximum safe inclination
	The product is designed to be transported by car
	The product is not intended to be used as a seat for transport in a motor vehicle
	Do not stand on the footrest

DISPOSAL OF WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT

This symbol on the products or accompanying documents means that used electrical or electronic products should not be disposed of with the normal household waste. For proper disposal, refurbishment or recycling, please take these products to the collection points for this type of waste.

Alternatively, in some EU countries or other European countries, you may return your product to your local retailer when purchasing a similar new product. By disposing of this product correctly, you will help to conserve valuable natural resources and support the prevention of potential negative effects on the environment and human health, which could result from incorrect disposal. For further information, please contact your local municipality or collection point waste collection points. Penalties may be imposed in accordance with local regulations if this product is disposed of in an inappropriate manner. For parties within the European Union, If you wish to discard electrical or electronic equipment, please obtain the necessary information from your retailer or supplier. Elimination in countries outside the European Union. This symbol is valid in the European Union. If you wish to decommission this appliance, please obtain the correct decommissioning information from your local authorities or from your retailer.

Contact details

Manufacturer

mdh Sp. z o.o.

Address: ul. Maratońska 104, 94-007 Łódź, Polska

tel. +48 42 674 83 84, fax. +48 42 636 52 21

www.mdh.pl

www.viteacare.com

mdh Sp. z o.o.
ul. Maratońska 104, 94-007 Łódź, Polska
tel. +48 42 674 83 84, fax. +48 42 636 52 21
www.mdh.pl
www.viteacare.com