CAREBED

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
English version on next page

Users with visual impairments can find the PDF files together with further information on our website:

< www.viteacare.com >.

☞ Contact your specialist dealer when required.

Alternatively users with visual impairments can have the documentation read out by a helper.

Table of contents

Introduction	35
Indications.....	36
Contraindications	36
User safety.....	37
General safety guidelines	38
Restrictions before using the bed.....	38
Residual risk	39
Safety control	39
Safety rules regarding electrical hazards	39
Electromagnetic Compatibility (EMC) information	40
Product description.....	41
Description of bed elements - model PREMIUM	42
Description of bed elements - model SMART	43
Wheel with brake	44
Installation of bed elements.....	44
Assembly of bed tops.....	45
Assembling the bed frame.....	45
Engine installation.....	46
Connecting the wires to the control unit.....	47
Installation of safety railings.....	47
Extension arm assembly	48
Use	49
Railings - lifting and lowering	49
Description and functions of the controller	50
Manual adjustment of leg flexion angle - „rastomat“ type mechanism.....	51
Optional equipment.....	51
Bedside table.....	51
Extension arms	52
Safety railing	52
Technical specifications	53
Electrical parameters	54
Care and maintenance	54
Troubleshooting.....	54
Cleaning.....	55
Disinfection.....	55
Product life cycle.....	56
Disposal of the rehabilitation bed	56
Maintenance plan.....	56
Replacement of worn parts.....	57
Symbols	58
Contact Information	60

Introduction

The CAREBED bed is a device created in accordance with the latest technologies. Its design allows for maximum adaptation to the user's needs. It combines a robust and durable construction, modern electronic control solutions and safety systems. To maintain this condition for a long time, professional service and care of the bed is necessary. Therefore, please read these instructions for use carefully. It will be helpful both during the initial setup and in the daily use of this product. In order to maintain safety and maintain the full efficiency of the product as long as possible, it should be regularly inspected at the service point (at least once a year). In case of any doubts or damage, please contact the distributor. Make sure that you have filled the attached warranty form and that you have submitted it to the distributor, as this will enable us to ensure the best possible service.

ATTENTION

Always observe the basic precautions according to the information given below regarding safety when using the product. Read the entire manual carefully before use. Please keep these instructions for future use.

WARNING!

Improper use can lead to death, serious injury, dangerous conditions or product malfunction.

ATTENTION!

Improper use may lead to personal injury and/or damage to the product.

SUGGESTIONS!

Follow the instructions below to keep the product in good working order.

Indications

The CAREBED medical device is used for home care and is intended for use as a bed. The product is intended to be used in private homes, healthcare homes and hospitals.

1. People with mobility difficulties who need support and stability when lying down.
2. Post-surgical patients who need an adjustable bed to improve comfort and support during recovery.
3. People with neurological, rheumatic or cardiac diseases who require prolonged lying in bed.
4. Patients in advanced stages of terminal illness who need a specialized bed for palliative care. People undergoing rehabilitation at home who need a bed with height and angle adjustment to support the healing process.

The device is intended to come into contact with the undamaged skin of the user (patient) and has been tested in accordance with the relevant compliance standards, including: ISO 10993-1:2018, EN ISO 10993-5:2009 i EN ISO 10993-10:2013.

The medical bed also meets the requirements of the standards:

EN60601-2-52:2010+A1:2015, EN60601-1-2:2015+A1:2021.

Contraindications

1. Patients without medical indications for the use of a specialised rehabilitation bed.
2. Patients who may not be able to use the bed function safely, even with caregiver supervision.
3. Persons with cardiac problems for whom the Trendelenburg position is medically contraindicated.
4. Patients with spinal injuries who require specialised orthopaedic beds with more precise adjustments and supports.

Medical device in accordance with Regulation (EU) 2017/745 of the European Parliament and Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC (Text with EEA relevance.)

User safety

WARNING!

Every “serious incident” that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State where the user resides.

ATTENTION!

The product is intended for persons over 12 years of age. If you are not sure whether your health condition allows you to use the device, consult a health care professional.

WARNING!

Do not modify the product in any way; this will void the warranty.

WARNING!

Exceeding the maximum user weight will void the warranty. The manufacturer is not responsible for injury and damage if the maximum user weight is exceeded.

WARNING!

Please read these instructions for use before using the product and keep them for future reference.

General safety guidelines

1. Please do not assemble the product before reading these instructions for use.
2. Any changes to the product made without the manufacturer's knowledge are the full responsibility of the person making the changes. Changes made without authorization will invalidate the guarantee and void the declaration of conformity with the applicable directives.
3. It is forbidden to use the product in any other way than described in these instructions for use.
4. The manufacturer is not responsible for malfunctions or damage to objects or persons caused by inappropriate use of the product.
5. Operators must carefully read this instructions for use, follow the information provided, and follow the procedures for use and maintenance of this product.
6. When the patient is not under supervision, the medical bed should be placed in the lower position with the side handrails raised to reduce the risk of falling and injury.
7. Avoid improper use of the power cord. Danger may be, for example, bends, cuts and other mechanical damage of the cord.
8. In the event that there are cables of other devices inside the medical bed, precautions should be taken to avoid them being crushed between other elements of the medical bed.
9. Due to the limited space under the medical bed, it is necessary to use appropriate lifting devices.
10. The product should be used and maintained in accordance with the instructions contained in this document and only use spare parts and technical support approved by the manufacturer.
11. Do not use the product for purposes other than those for which it was designed and built.
12. Keep this instructions for use for consultation or training of personnel. In case of sale or use by new users, it should be stored together with the product.

Restrictions before using the bed

WARNING!

- Cannot be used when moving parts are not joined or assembled, or show defects.
- The device must not be used in an environment exposed to splashes of water or in critical and adverse weather conditions.
- It must not be used if it shows defects - this may cause unexpected danger situations.

Residual risk

If the requirements presented in this manual are not respected, the following risks may occur:

- Risk of crushing the limbs if the limbs are placed under the moving parts during downward movements.
- Risk of injury or fire if electrical components come into contact with liquids.
- Risk of damage to the product in the event of critical or adverse weather conditions.


Safety control

ATTENTION!

The following points should be checked before use:

- Electrical connections and wiring harnesses – protect and remove corrosion.
- Brakes – check if all the brakes in the wheels are working.
- Actuator mechanism - make sure when the bed is fully extended.
- Condition of the wheels – make sure that they are not damaged and excessively worn - if there is a problem, contact the distributor.

Maximum weight of the patient is 135 kg.

Maximum load capacity is 180 kg

- (mattress, patient, additional equipment)

Safety rules regarding electrical hazards

Power must be supplied to the product by a protected line to avoid overloading by the automatic circuit breaker. These devices must provide sufficient electrical power to properly protect the power cables.

1. If it is necessary to use extension cords, the power cables must be terminated with a non-separable plug or other devices that prevent water from entering during the bed cleaning process.
2. To avoid damage to the power cord, the rehabilitation bed must be equipped with adequate equipment to keep the power cord away from any moving part or mechanism of the bed during its operation or movement.
3. The power cables must be adequately protected against damage that may result if they come into contact with one or more moving parts or rub against sharp angles and surfaces inside the rehabilitation bed.

**WARNING!**

Before operating the moving parts of the bed, make sure that they do not damage any pipes, electrical cables, catheters or other devices, and that they do not cause discomfort and other problems for the user (patient).

**WARNING!**

The casters are only used to move the bed. In a situation where sanitary procedures, therapy, etc. are carried out, the wheels should be blocked.

Electromagnetic Compatibility (EMC) information

To ensure proper use of the bed without electromagnetic interference, it is recommended to use only dedicated cables and accessories indicated by the manufacturer. Potential electromagnetic interference from surrounding equipment should not significantly affect the operation of the bed during the expected service life, as long as the use is as intended and as described in the main instructions for use.

Use of accessories and wires: Only use accessories, transformers and wires provided by the bed manufacturer. The use of other components may lead to electromagnetic interference or reduce the electromagnetic immunity of the bed and cause it to malfunction.

Prohibition on the use of electro-surgical equipment: Do not use electro-surgical equipment on the bed as this may affect the functioning of the bed in an unpredictable way.

Avoid close proximity to other equipment: If it is necessary to use the appliance in close proximity to other equipment, monitor that all equipment is working properly. Malfunctions may result from close proximity to other devices.

Portable communication devices: Communication devices (e.g. mobile phones, radios) operating at high frequencies and their accessories (antenna cables, external antennas) should not be used at a distance of less than 30 cm from electrical elements and wires in the described bed. Improper attention to this recommendation may limit the device's functions.

Electromagnetic environment: The bed is designed to operate in a specific electromagnetic environment. The user should make sure that the bed will be used in the right conditions.

WARNING!

The use of other accessories, other transmitters and other wires than those provided by the bed manufacturer may cause increased electromagnetic interference emissions or reduced electromagnetic immunity of the bed and cause it to malfunction.

WARNING!

The use of electro-surgical equipment on the bed is prohibited, as this can lead to unexpected interference with bed function.

WARNING!

Absolutely avoid using this device in the immediate area of other devices, as this could result in malfunctions. If such use is necessary, the described device and other devices should be observed to ensure that they work properly.

WARNING!

Portable communication devices operating at high frequencies (radio devices, mobile phones, etc.) with accessories (e.g. antenna cables and external antennas) should not be used at a distance of less than 30 cm from electrical elements and wires in the described bed. Failure to comply with the above recommendation may result in a limitation of the operating functions of the device.

Product description

Typical features of home-care beds include adjustable lying surfaces, adjustable height to at least 65 cm for ergonomic care and locking castors with a minimum diameter of 10 cm and lockable castors. Multi-segmented, often electronically powered lying surfaces can be adapted to a variety of positions, such as comfortable sitting positions, shock positions, or cardiac positions. Rehabilitation beds are also often equipped with pull-up aids (trapezoidal bars) and/or bed sides (side rails) to prevent falls. Thanks to the adjustable height of the home care bed, they enable both ergonomic working height for nurses and therapists, as well as a number of appropriate positions for easy access to the patient.

Description of bed elements - model PREMIUM


Illustration 1. List of bed elements - model PREMIUM.

No.	Description
1	Extension arm
2	Handle
3	Bed base - top - head side
4	Backrest section
5	Wired controller
6	Wheels with brakes
7	Bed base - bottom - leg side
8	Side rails
9	Thigh support section
10	Pelvic support section
11	Mattress side locks
12	Leg support section leg flexion adjustment system "rastomat" type
13	Electric actuator
14	Control unit with motors
Equipment options	Table above the bed, Extension arm I, Extension arm II

Description of bed elements - model SMART


Illustration 2. List of bed elements - model SMART.

No.	Description
1	Extension arm
2	Handle
3	Bed base - top - head side
4	Backrest section
5	Wired controller
6	Wheels with brakes
7	Bed base - bottom - leg side
8	Side rails
9	Thigh support section
10	Pelvic support section
11	Mattress side locks
12	Leg support section leg flexion adjustment system "rastomat" type
13	Electric actuator
14	Control unit with motors
Equipment options	Table above the bed, Extension arm I, Extension arm II

Wheel with brake

To apply the brake, press the brake down with your foot **1** until the wheel is completely locked.

To release the brake, press the brake release lever **2**


Illustration 3. Wheel details with brake.

WARNING!

- Wear full shoes when operating the bed to avoid finger injuries.
- Make sure that the brake is applied on at least three wheels.

Installation of bed elements

WARNING!

The product must be transported in accordance with the instructions on the packaging.

WARNING!

After unpacking the product, make sure that it is not damaged. Damage caused during improper transport, unpacking or assembly is not covered by the warranty.

WARNING!

The bed should be lifted and moved slowly to avoid injuries caused by bio-mechanical overload of the spine.

- if handling the product is necessary, we recommend the use of protective gloves and the presence of at least two operators to lift the product.
- at least two operators to assemble the product.

Illustration 4. Bed components - preparing for assembly.

Assembly of bed tops

First prepare the frame to be connected to the top of the bed. A safe position for assembly is to place the top on a flat surface so that there is easy access to the bed assembly mount.

- Unlock the bracket locking mechanism by extending and then turning by 90° so that the pin can be used to assemble the frame (Illustration 5.).


Illustration 5. Assembly of the frame to the top of the bed.

Slide the two bed brackets into the holders on top of the bed.

- The locking mechanism must be turned by 90° to ensure that the pin engages in the holes of the bed hooks (Illustration 6.).
- Check the connection is stable and the bolt locks against pulling out of the bed hooks.


Illustration 6. Assembly of the frame to the top of the bed.

Assembling the bed frame

Prepare the connected tops to the frame by tilting them to the side so that the profile can be inserted, as shown in Illustration 7.

- Then lift and carefully push the profiles as far as possible against the support profile flange (Illustration 8.).


Illustration 7. Assembly of the bed frame.

Secure the folded frame of the bed by tightening the thumbscrews firmly (Illustration 8.).

- Check the connection is sufficiently stable and the screws are tightened.


Illustration 8. Butterfly-head screws.

Engine installation

After assembling the bed frame, before fixing it to the top and bottom - the base must be inverted so that its hooks fixing the module with motors are facing upwards.

- To locate the hooks for assembly, pay attention to the icons placed on the module with motors, which will facilitate correct orientation. Figures 6 and 7.
- Pull out the two covers from the slots of the module where the bed hooks are to be placed (illustration 5 **A**).
- Flip the module and align the slots with the two protruding hooks and press them in.
- Place the two covers back into the module over the hooks to secure it.

IMPORTANT! - the module should be fixed in the right place at the hooks! as shown in illustration 5 **B**.



Illustration 9. Module flaps and motors module mounting location.


Illustration 10. Back support section.


Illustration 11. Leg support section.

Connecting the wires to the control unit

The connection of the actuator wires to the motor is shown in Illustration 9.

- A** Connecting the top head side.
- B** Connecting the bottom foot side.
- C** Connecting the control panel.


Illustration 12. Control unit.

Follow the instructions below to prevent damage to the wires when starting the motors:

- Connect the bottom top (foot side) motor cable directly to the motor module.
- Connect the top (head side) motor cable directly to the motor module.
- Connect the control panel cable to the motor module.
- Connect the supply cable to the 230 V power source.
- Set the actuators of both, head and foot side tops in the highest position.
- Perform all movements of the bed lift section, bed frame tops and Trendelenburg position to check for obstacles and folds between cables and other moving parts.
- Ensure that all cables are free and there is no risk of them becoming trapped in moving parts of the bed.

WARNING!

Care must be taken when installing the bed tops as there is a risk of crushing and damage to the product if the actuators are not securely locked in the bed frame.

Installation of safety railings.

- Raise the base of the bed about 1/4 of the height to make it easier to insert the panels.
- Insert the panels from below until the stops snap into position.


Illustration 13. Connector rail.


Illustration 14. Mounting socket for railing connectors.

Illustration 16. Place the connecting rail in the rail socket.

Illustration 17. Safety rail is ready to be inserted into the bed frame.

Illustration 18. Slide the rails into the bed frame rails.

Illustration 19. Model PREMIUM - slide in and tighten the bed rails locking mechanism.

Illustration 20. In the SMART model is a button to lock the railing.


Illustration 21. Set the appropriate height of the rails.

Extension arm assembly

The installation of the extension arm should start with the location of the socket, which is located on the bed frame rails. Insert the extension arm into the seat so that the locking pin is in the position shown in illustration 10. Then mount the handle on the extension arm hook.


Illustration 15. Extension arm mounting socket.

Use

Daily use of the product is only allowed when all protections and functional components are in working order and properly maintained. The bed is intended for private use or in nursing homes, hospitals, by users over 12 years of age. The only possible use is described in this instructions for use. We recommend that you read the instructions carefully to be fully aware of the functionality of the product.

WARNING!

Store a medical device with electrical power supply in a place where the plug can be easily disconnected from the power outlet.

WARNING!

Before moving the bed, remove the plug from the power socket.

WARNING!

To reduce the risk of falling while the patient is unsupervised, the bed should be set at the lowest height and the rails raised when unsupervised, .

SUGGESTIONS!

The wired controller is equipped with a safety lock that switches on automatically when the device has not been used for more than 30-40 seconds.

ATTENTION!

The device is not equipped with a main switch, just plug it into the mains socket to access all the functions. Removing the plug from the mains socket acts as an emergency response to the power supply.

Railings - lifting and lowering

To move the side rails, press the brown button that activates the locking pin at the top and bottom, then follow these instructions:

• Raising :

- grasp the end of the panel and pull up the higher barrier to activate the locking pin - a click is an indication of the correct lock.
- repeat this operation for the top and then the bottom, or in reverse.

• Lowering :

- grasp the end of the railing, press the button on the bed frame, then lower the railing.
- repeat this operation with the other end.

Illustration 22. Railing lock release button.

Description and functions of the controller

No.	Description
1	Backrest section adjustment
2	Leg support section adjustment
3	Backrest section and leg support height adjustment
4	Adjusting the height of the entire bed frame
5	Activation of * Trendelenburg position / Anti-Trendelenburg
6	Power activation indicator light


Illustration 23. Bed functions controller.

WARNING!

If you are not sure whether your health condition allows you to use the device, consult a health care professional.

* The Trendelenburg position is a positioning method in which the patient lies on his back with his head and chest positioned below the level of the lower limbs, at an angle 15° to 30° .

The CAREBED bed allows you to obtain a position of 12° .

Locking the remote control function with the key.

- For patient safety, it is possible to lock the function of the remote control. To block the functions of the remote control, turn the triangular lock towards a closed padlock with the supplied key.
- To unlock the function of the remote control, turn the lock towards the open padlock.

Illustration 24. Remote control lock.

Manual adjustment of leg flexion angle - „rastomat” type mechanism

To adjust the bending angle of the legs, lift the lower section of the sleep surface using the side mattress locks, as shown in illustration 26. The “rastomat” type locking adjustment system is unlocked by lifting the section as high as possible to return it to the bending angle previously set by the actuator.


Illustration 25. Adjustment of leg bending angle.

Optional equipment

Bedside table

Table assembly.

The table has two parts, which have to be connected.

- Place the base with castors on level ground and lock the brakes. Loosen the height adjustment screw (Illustration 26).
- Install the top element with the table-top by sliding the profile into the profile of the base support, paying attention to the plastic insert, which must be inserted together with the top profile.
- After adjusting the height, tighten the adjustment screw and check that the connection is stable.
- The table is adjustable for the angle of the top and the angle of the base support.
- Adjustment is possible by extending the lever shown in Illustration 27. Once the tilt angle is set, check that the locking pin is in the plate hole and that the connection is stable.


Illustration 26. Table assembly.

Illustration 27. Tilt angle adjustment.

Extension arms


Illustration 28. Extension arm I-R.

Illustration 29. Extension arm II-S.

Safety railing



Illustration 30. Safety railing.

Technical specifications

Model	SMART	PREMIUM
Total weight of bed	91 kg	98 kg
Weight of the main bed frame	80 kg	86 kg
Safety rails weight (side panels)	9 KG	
Weight of extension arm	3 kg	
Maximum user weight (patient)	135 kg	
Maximum load capacity (mattress+patient +equipment)	180 kg	
Maximum extension arm load	75 kg	
Table option - maximum table load	45 kg	
Minimum patient height	1400 mm	
Wheels with brake	125 mm	
Maximum bed size	2120x1030 mm	2150x1020 mm
Permissible dimensions of the mattress	1900x830mm	
Height of the bed base from (min.-max.)	390-800 mm	
Height of top/bottom panels	460 mm	
Bed base size	2000x900 mm	
Dimensions of the head section	800x740 mm	
Dimensions of pelvic section	800x220 mm	
Dimensions of leg sections	800x300 mm	
Dimensions of foot section	800x460 mm	
Backrest section motion Angle (Max.)	70°	
Leg section motion angle (Max.)	35°	
Maximum lift height of the leg section from the bed base	210 mm	
Storage and packaging temperature (in a dry environment)	From -10° to +50°	
Operating environment temperature range	From -10° to +50°	
Maximum mattress weight	30 kg	
Trendelenburg position / Anti-Trendelenburg	12°	
Bottom slats (crossbars)	metal or wooden	
Extension arm I-R - maximum load	75 kg	
Extension arm II-S - maximum load	100 kg	

Electrical parameters

Model	SMART	PREMIUM
Supply voltage	230V AC	
Operating voltage	24V DC	
Maximum power consumption	120 W	
Maximum actuator force	6000 N	
Maximum control panel voltage	24V DC	
Motor unit (IP rating)(Ingress Protection)	IP20 - for central motor IP54 - for actuator	
Control unit (IP rating)(Ingress Protection)	IP54	
Line cord length	3500 mm	

Care and maintenance

ATTENTION!

Insufficient or neglected care and maintenance of the bed leads to limitation of the manufacturer's liability.

Troubleshooting

Malfunction	Solution
The bed does not rise or fall.	Make sure that the plug is connected correctly. Repeat the process of connecting the remote control. Disconnect the power supply from the main socket and check, clean the pins of the connection of the control controller to the motor module. Reduce the load on the bed.
The light does not work.	Contact the distributor - replace the LED lamp.
Irregular movement during product operation.	Disconnect the controller and clean dirty contacts. Contact the distributor.
The safety rails are not inserted correctly.	Lubricate the railing guides with slide guide grease. Contact the distributor.
Unusual sounds when the device is operating.	Contact the manufacturer's service center.

Cleaning

In order to ensure long-term and efficient operation of the bed, it should be subjected to basic maintenance and inspections carried out by the service point. Inspections should take place once a year. If you have any doubts or damage, contact the service center, which has all the competences to solve your problem.

WARNING!

Observe the safety instructions given by the manufacturer of cleaning and disinfecting products. Failure to comply with this requirement may cause skin irritation or premature destruction of the product.

It is not permitted to use:

- solvents,
- toilet cleaners,
- agents with chlorine content,
- sharp brushes and hard objects,
- corrosive cleaners.

Do not clean the product with pressure or steam units!

Disinfection

If the product is used by many people (e.g. in a social care home), commercial disinfectants should be used. Before proceeding with disinfection, clean the surface. Only biocidal products authorised in Poland by the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products or by another competent authority of a given EU Member State may be used for disinfection.

WARNING!

Disinfectants may damage the outer surface, reducing the long-term functionality of the components. Always observe the instructions concerning the use of the products specified by the manufacturer.

ATTENTION!

In order to ensure proper handling during cleaning and disinfection, we recommend appropriate training for users and medical personnel. During this training, it should be clarified that it is necessary to follow the following guidelines:

- **Transportation of a clean bed:** When transporting a clean bed to the patient's apartment, care must be taken to ensure that the bed is not soiled or biologically contaminated.
- **Disassembly and preparation of the bed:** Users or medical personnel should immediately wash and disinfect the bed after its disassembly. It is worth informing them about the specifics of the cleaning and disinfection process and providing guidelines for preparing the bed again.

- **Clean materials:** Use only fresh, clean materials to prepare the bed.
- **Disinfection of the hands:** After finishing preparing the bed, specialists should disinfect the hands before moving to the next place. A suitable hand sanitizer dispenser should be available.
- **Direct cleaning on site:** Direct cleaning of the bed on site has the advantage of not moving contaminated beds or bed parts together with clean ones in the same vehicle. This helps avoid the transmission of pathogens and reduces the risk of hospital infection.
- **Storage of the bed:** If the beds are not used immediately, they should be stored under cover to protect them from dust, contamination and contamination.

Product life cycle

The expected life of the medical device is 5 years if it is used as intended and if all maintenance and repair activities are performed as required. The life cycle of the medical device also depends on the frequency of use, the operating environment and care of the device. The life cycle of the device can be extended using spare parts. Spare parts are normally available for 5 years after the model is withdrawn from production. The expected life cycle of the medical device is not equivalent to an additional guarantee.

Disposal of the rehabilitation bed

The bed must be disposed of in accordance with current local and national regulations. For proper disposal of packaging, metal components, plastic components, neoprene, silicone and polyurethane materials, contact your local waste disposal company or authorised distributor.

WARNING!

Remember that safe handling of sharp-edged waste is crucial for health and safety. Observe the guidelines and follow the recommendations for the disposal of these products.

Maintenance plan

These operations can be performed by the user. The device requires regular maintenance, just like any other product, so that it is kept in a safe and perfect condition. Food, hair and any other accumulated contamination - such as dust - must be removed daily. A domestic vacuum cleaner can be used for this operation. The table below shows the recommended maintenance frequency. Constant monitoring and user attention should encourage actual unnecessary activities.

When	Description of activities
Monthly	<p>Cleaning of metal parts. Clean these parts according to their intended use and needs.</p> <p>Cleaning of plastic parts. Clean these parts according to their intended use and needs.</p> <p>Cleaning of wooden parts. Clean the parts according to their intended use and needs.</p> <p><u>Make sure that the electric motors are working properly. Lubricate the fasteners of the back, pelvis and leg support sections - see description of bed elements on pages 11-12.</u></p>
Every 6 months	<p>General inspection: Check the wheels and tighten the screws Check the paint for rust or chipping Make sure that the plastic parts are intact and that they are not worn.</p> <p><u>Make sure that the electric motors are working properly. Lubricate the fasteners of the back, pelvis and leg support sections - see description of bed elements on pages 11-12.</u></p>
In case of transfer to third parties or renewal	<p>Warning! Disinfection Disinfectants may only be used by qualified personnel only disinfection of this equipment must be carried out by specialized centers in accordance with the regulations issued by the competent health authority. For more information, contact an authorized support service.</p>

Replacement of worn parts

WARNING!


Replacement of worn parts may only be carried out by qualified personnel of the manufacturer.

If damaged or worn parts are revealed during the annual routine inspection, they must be replaced immediately as recommended by the personnel conducting the inspection. We recommend an independent check at least once every 12 months. When replacing damaged parts, use original spare parts. The manufacturer is not responsible for damage to objects or persons in the event of non-compliance with this recommendation.

Symbols

The symbols used below are used to identify warnings, mandatory actions and prohibited actions. You should read everything and understand their meaning.

	Medical device
	Date of manufacture
	Manufacturer
	Distributor
	Authorised representative in the EU
	Conformité Européene (compliance with the European Community)
	Catalog number (to re-order)
	Serial number
	LOT number
	Unique Device Identifier
	Refer to the instructions for use
	Warning

	Protect from moisture
	Possibility of jamming fingers
	Permitted temperature
	Indoor use only
	Type B applied part
	Class II equipment
IPX4	Protection against splashing water from all directions
	Maximum patient weight
	Safe working load
	The packaging material is 100% recyclable. Regulations must be observed for disposal. Keep packaging material (plastic bags, polystyrene, etc.) out of reach of children as a potential source of danger.

Selective collection of electrical and electronic equipment.

This symbol on the products or accompanying documents indicates that used electrical or electronic products and batteries must not be disposed of in municipal waste. For proper disposal, restoration or recycling, such products should be returned to collection points for this type of waste. Alternatively, in some countries of the European Union or other European countries, you can return your products to a local seller at the time of purchase of a similar new product. By correctly disposing of this product, you will help preserve valuable natural resources and support the prevention of potential negative impacts on the environment and human health, which could result from incorrect waste disposal. More information is available from the communal offices or waste collection sites. In case of incorrect disposal of this product, penalties may be imposed in accordance with local regulations. If you want to dispose of an electrical or electronic device, get the information you need from your dealer or supplier.

Liquidation in countries outside the European Union: This symbol is valid in the European Union. If you want to dispose of this device, please obtain the necessary information on the correct disposal at your local authorities or from your dealer.

Contact Information

Manufacturer

Sp. z o.o.

Address: Maratońska 104, 94-007 Łódź, Poland

tel. +48 42 674 83 84, fax. +48 42 636 52 21

www.mdh.pl

www.viteacare.com

We reserve the right to make technical changes due to technical progress.