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PRODUCT DESCRIPTION:

Back rest is made of highly durable aluminium frame, protected with non-toxic lacquer layer. Simple frame adjustment system enables for five different positions and height of the rest. The rest itself is made of easy to clean material.

INDICATIONS:

The back rest simplifies all the hygiene and care tasks (e.g. feeding) that are being carried out for laying patients. It enables for easy retaining of half-sitting position which guarantees passive orthostatic position of a patient – thus it eliminates consequences of prolonged immobilization, enhances comfort and wellbeing of a patient.

CONTRAINDICATIONS:

Contraindications against half-sitting position.

USING (ADJUSTMENTS):

Position the back rest in the upper part of a patient's bed, so that his torso is laying on the rest. Adjust the height of the back rest by positioning the rear part of the frame in an adequate latch indent, situated in the lower element of the frame.

SAFETY OF USAGE:

1. Remember that staying in a single position for a prolonged period of time may lead to irritations and skin disorders of a patient.
2. Before letting patient lay against the Rest it is required to check the height of the rest has been adjusted and secured – the rear frame element is correctly positioned in the latch system.
3. One has to pay attention to eliminate the possibility of frame elements pressing patient's body which may lead to irritations in case of patients suffering from skin trophic disorders.

MAINTENANCE AND SERVICE:

1. Frame of the rest is to be cleaned with a damp cloth.
2. Support of the rest is to be cleaned with a damp cloth and soft detergent.

For the product to maintain its proper characteristics, it is extremely important to maintain all its elements clean at all times. Lack of care in this regard, will lead to quick loss of the initial utility and aesthetical properties.

Rest width: 62 cm

Materials: Aluminium, polyester fabric

Batch No/production date: printed on the product.

Usage period: 2 years since purchase date.