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**English version** on next page

**Users with visual impairments can find the PDF files together with further information on our website:**

**<[www.viteacare.com](http://www.viteacare.com)>.**

**👉 Contact your specialist dealer when required.**

**Alternatively users with visual impairments can have the documentation read out by a helper.**

# Table of contents

<b>Introduction</b> .....	<b>48</b>
<b>Product characteristics</b> .....	<b>49</b>
Intended use / indications.....	49
Contraindications .....	49
<b>Safety instructions</b> .....	<b>50</b>
General safety rules .....	51
Modifications .....	51
Security check.....	52
Capacity .....	52
Transport the wheelchair.....	52
Inclination guidance.....	53
Driving outdoors and precautions in bad weather conditions .....	53
Physical limitations .....	54
Information on electromagnetic compatibility (EMC) .....	54
<b>Product information</b> .....	<b>55</b>
Clutch lever for disengaging the drive .....	57
Power button / status indicator.....	57
Battery charge indicator.....	57
Speed selection switch .....	58
Seat function selection button.....	58
Battery cut-circuit key.....	58
Speed control dial .....	58
Horn.....	59
Hazard lights.....	59
Directions lights .....	59
Control lever for drive and wheelchair functions - joystick.....	59
<b>Operate</b> .....	<b>60</b>
Before you start using.....	60
Getting on the wheelchair .....	60
Pre-ride inspection and adjustment.....	60
Operating a wheelchair.....	61
Getting off the wheelchair .....	61
<b>Regulation and adjustment of wheelchair components</b> .....	<b>62</b>
Headrest adjustment.....	62
Backrest angle adjustment .....	62
Seat and backrest angle adjustment.....	63
Control panel position adjustment.....	64
Armrests width and height adjustment .....	64
Footrest adjustment .....	65
Footrest with mechanical height adjustment .....	66
Removing the footrests.....	67
Footrest length adjustment.....	67

<b>Disassembly and assembly the wheelchair.....</b>	<b>67</b>
Disassembly.....	67
Disassembly procedures.....	68
Assembly.....	68
<b>Basic troubleshooting.....</b>	<b>69</b>
<b>Batteries and charging.....</b>	<b>70</b>
First charge.....	70
Connecting the charger.....	71
Operating temperatures.....	71
Battery operation and charging.....	71
Battery charging process.....	71
<b>Reporting technical defects and repairing the device.....</b>	<b>73</b>
<b>Care and maintenance.....</b>	<b>73</b>
Tires.....	73
Wheel Replacement.....	73
External surfaces.....	73
Cleaning.....	74
Disinfection.....	74
Battery terminal connections.....	74
Electric wiring harnesses.....	74
Daily checks.....	74
Weekly checks.....	75
Monthly checks.....	75
Yearly checks.....	75
ABS plastic covers.....	75
Axle bearings and motor/transaxle.....	75
Console, Charger, and Rear Electronics.....	75
Nylon Lock Nut Replacement.....	75
Storage of the wheelchair.....	76
Re-use (reusable product).....	76
Product life cycle.....	76
Replacement parts.....	76
Disposal of an electric wheelchair.....	77
<b>Warranty.....</b>	<b>77</b>
<b>Technical parameters.....</b>	<b>78</b>
<b>Product safety symbols.....</b>	<b>80</b>
Symbols.....	80
<b>Contact details.....</b>	<b>84</b>

## Introduction

The STAR wheelchair has been designed and manufactured to the highest quality standards. Before using this product for the first time, please read these operating instructions. Improper use of the wheelchair may result in personal injury or traffic accidents. If you still have questions after reading this manual, please contact your specialist dealer.

Ensure that the 'Warranty Card' document is completed and enclosed and handed over to the distributor, as this will enable claims to be processed efficiently.

# ATTENTION

Always follow the basic precautions according to the information given below regarding safety when using the product. Read the entire manual carefully before use. Keep this manual for future reference.

### WARNING!

Improper use can lead to death, serious injury, unsafe conditions or product malfunction.

### ATTENTION!

Improper use may lead to personal injury and/or damage to the wheelchair.

### SUGGESTION!

Follow the instructions below to keep your wheelchair in good working order.

If the electric wheelchair is used by another person, you must provide them with copies of these instructions for use.

These instructions for use relate to the design and specifications of the product at the time of publication. Due to design changes, some illustrations and photos in the manual may not correspond to the wheelchair you have purchased. We reserve the right to make design changes. The product shipped in its original packaging requires partial assembly - follow the information in the "Disassembly and assembly the wheelchair" section. The manufacturer shall not be responsible for damage resulting from improper or unsafe use of the medical device, if such use was against the instructions and recommendations contained in the product documentation. The manufacturer shall ensure that the product has been designed and manufactured in accordance with the applicable safety and quality standards.

## **Product characteristics**

An electrically-powered wheelchair, intended as a medical device, is a mobility aid for people with disabilities, the elderly or those with reduced mobility.

**Construction:** The Star wheelchair is available in a four-wheel version which provides good stability.

**Drive:** The vehicle is powered by an electric motor.

**Energy range:** On a single charge it can be up to several tens of kilometers, depending on battery capacity and conditions of use.

**Safety:** The vehicle has road lights, electromagnetic brakes and seat belts.

**Comfort:** The product is equipped with an adjustable seat, adjustable steering column angle, armrests, transport basket and cushioned suspension..

## **Intended use / indications**

The product is intended for disabled people with musculoskeletal conditions to increase their mobility. Designed for people with paralysis and paresis, amputees, people with stiff lower limbs and those who periodically need to use a wheelchair.

## **Contraindications**

Contraindications to upright standing. Injury or paresis of the upper limbs preventing the use of a power wheelchair.

## Safety instructions

- a. Please note that this is an electric wheelchair and is intended for outdoor use. Special care should be taken when driving around bends, obstacles and inclines.
- b. Always ensure that the electric wheelchair is fully assembled before getting on.
- c. Read and follow the information in the instructions for use. It is forbidden to use this product and the available optional equipment without first carefully reading and understanding this manual. If you do not understand the warnings, cautions or instructions, contact your supplier, otherwise personal injury or accident may result. Never attempt to use the electric wheelchair beyond its limitations as described in this manual.
- d. Place your feet on the footrest when driving. Do not stand on the footrest.
- e. To operate an electric wheelchair safely, ride on flat surfaces where all four wheels make sufficient contact with the ground.
- f. Keep metal objects away from the battery connections. Risk of electric shock.
- g. When negotiating any corner, reduce speed and keep the centre of gravity stable.
- h. When reversing your electric wheelchair, make sure there are no obstacles behind you.
- i. It is prohibited for children to play near or operate the electric wheelchair.
- j. Some parts of the power wheelchair are susceptible to temperature changes. The controller can only operate at temperatures between +5°C to +40°C.
- k. In extremely cold temperatures, the batteries may freeze and the electric wheelchair may not operate. In extremely high temperatures, the wheelchair may run at a lower speed due to the controller's heat protection function, which prevents damage to the motors and other electrical components.
- l. Store in clean and dry conditions.
- m. It is forbidden to drive at night without lighting.
- n. Do not remove the anti-tip wheels or transit wheels.
- o. Keep hands and feet away from moving parts when driving. Take care that loose clothing does not get caught in the drive wheels.
- p. Do not connect the battery charger to an extension cord.
- q. Surface temperatures of the wheelchair can increase when exposed to external sources of heat.
- r. Disassembly of the controller, motor or charger by anyone other than an authorized service is prohibited and will invalidate any applicable warranty conditions.
- s. It is forbidden to drive the wheelchair on expressways and motorways.
- t. It is forbidden to use a power wheelchair under the influence of alcohol, medication or any other substance affecting the ability to drive.
- u. Be careful when driving in busy areas and shopping centers.
- v. Under no circumstances should the wheelchair be used as a seat when being transported by motor vehicles.
- w. Do not attempt to lift the wheelchair by gripping anything other than the frame (for example, the seat, armrest or chassis components).

- x. Please note that the temperature of the components and frame parts may increase when exposed to the sun or other heat source, while at low temperatures there is a danger of very high cooling of the metal parts of the wheelchair. This can cause burns or frostbite on parts of the body.

**WARNING!**  
The product is intended for use by people with disabilities and movement disorders. If you are not sure whether your health condition allows you to use the product, consult a healthcare professional.

**WARNING!**  
Any “serious incident” related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides.

## General safety rules

Read the instructions for use carefully with understanding before using the electric wheelchair.

Users may encounter difficulties with maneuvering when passing through narrow doorways, ascending and descending ramps, driving around corners and navigating uneven terrain. Reduce speed, take your time and maneuver the electric wheelchair carefully.

## Modifications

Do not modify, remove, disable or add any parts, components or features to your electric wheelchair. The manufacturer offers a wide range of electric wheelchair accessories to meet your needs.

**WARNING!**  
Do not modify the electric wheelchair in any way - this will void the warranty.

## Security check

### ATTENTION!

Check the following before use:

- Wheelchair frame - is there any damage
- Electrical connections and wiring harnesses - check for damage
- Batteries - fully charged and not corroded.
- Brakes - check that the electromagnetic brake is functioning.
- Condition of tires - make sure they are not damaged or excessively worn - if there is a problem, contact your dealer.

## Capacity

STAR wheelchair is designed and built to transport one person with a maximum weight of 136 kg. The wheelchair is not designed to carry goods or objects other than those previously described.

### Warning!

Exceeding the maximum user weight will void the warranty. The manufacturer is not responsible for injury or damage caused by failure to comply with the maximum user weight.


## Transport the wheelchair

### Warning!

Do not sit in an electric wheelchair while it is being transported in a motor vehicle.

When transporting the electric wheelchair by vehicle, store it safely in the back of a van, truck or in the boot of a car. Adjustable parts should be disassembled or adequately secured for transport. To ensure safety and to avoid damage to the powered wheelchair during transport, tie-down straps must be hooked to the attachment points (A) and (B) on the wheelchair and in accordance with the belt system supplier's instructions.

## Inclination guidance

Remain cautious when going up and down hills. If necessary, reduce your speed when going downhill. If you are moving up a hill try to keep the wheelchair moving forward. If necessary, stop. Be careful and slowly accelerate your wheelchair forward. If you are descending, reduce the speed of the electric wheelchair to the lowest setting and then proceed with caution. If the wheelchair is descending faster than you anticipated, slowly release the drive control lever until it stops and then push the lever slightly to continue descending safely.

### WARNING!

Drive the electric wheelchair forward down or up a slope, never reverse. Irregular movements or making back and forth movements may increase the probability of tipping over. Changes in slope can reduce stability. Under no circumstances should you travel up a slope with a gradient greater than the maximum stability angle of this electric wheelchair.

## Driving outdoors and precautions in bad weather conditions

Although the electric wheelchair is designed to be used indoors to help with daily living activities, the manufacturer has designed it to perform perfectly outdoors on level surfaces.

### ATTENTION!

When driving outdoors, avoid uneven terrain, soft surfaces, tall grass, loose gravel, loose sand, fresh and salt water, stream banks, lakes and oceans. If you are unsure about a surface, avoid it.

Exposure to adverse weather conditions should always be avoided. If you find yourself in bad weather conditions while using your electric wheelchair, immediately go to the nearest shelter. Dry the electric wheelchair completely before use, before charging or before storing.

## Physical limitations

Use awareness, caution, prudence and common sense when operating the wheelchair. Always be aware of your own limitations and rules of use when operating an electric wheelchair.

**WARNING!**

It is forbidden to use the wheelchair while under the influence of alcohol or psychotropic drugs. Never smoke while sitting in the wheelchair. Damaged or worn upholstery must be replaced as it increases the risk of fire.

**ATTENTION!**

When operating an electric wheelchair, please be aware of precautions, warnings and safety issues if you are taking medications prescribed by a doctor or available without a prescription.

**WARNING!**

Drive the electric wheelchair in a straight line up or down a slope, never reverse. Irregular forward and backward movements with changes in slope can increase the risk of wheelchair rollover.

## Information on electromagnetic compatibility (EMC)

Despite complying with applicable EMC directives and standards, an electrically powered wheelchair may interfere with the operation of other electrical devices, such as mobile phones or alarm systems. In addition, this may affect the performance of the wheelchair itself. If you notice this behavior, it is a recommended idea to switch off the third-party device or power wheelchair and move it out of interference range. Remember that other devices emitting electromagnetic fields can also affect the operation of the wheelchair. An electromagnetic risk analysis was carried out.

The wheelchair maintains basic safety and functional performance in the presence of electromagnetic interference. All sources were identified and evaluated electromagnetic risks were identified and evaluated.

Level of electromagnetic immunity:

The product has been tested and passed the 20 V/m immunity level.

**WARNING!**

Pay attention and stay away from nearby transmitters, such as radio or TV stations.

## Product information

The STAR electric power wheelchair is designed for outdoor and indoor use. The following illustration shows all parts relevant to the user/driver. The drive disconnect clutch levers are located on the sides between the front and rear wheels.

Illustration 1. STAR model Orange wheelchair component list - front.

- |                  |                                     |
|------------------|-------------------------------------|
| 1. Headrest      | 8. Rear wheel                       |
| 2. Backrest      | 9. Drive disengagement clutch lever |
| 3. Armrest       | 10. Front wheel                     |
| 4. Control panel | 11. Side cover                      |
| 5. Seat          | 12. Footrest                        |
| 6. Safety belt   | 13. Footrest platform               |
| 7. Headlights    | 14. Guide handle                    |


Illustration 2. List of STAR model Orange wheelchair components - rear.

- |                                |                    |
|--------------------------------|--------------------|
| 1. Headrest adjustment knobs   | 5. Battery case    |
| 2. Headrest attachment bracket | 6. Anti-tip wheels |
| 3. Handle                      | 7. Rear lights     |
| 4. Power cut-circuit key       |                    |

1. Battery charge indicator.
2. Power on/off button. Illumination of the button indicates that the power is turned on (if it flashes red, an error has occurred).
3. Speed selection switch.
4. LEDs indicating the selected speed.
5. Sound signal button - horn.
6. Left indicator button and hazard lights.
7. Speed selection dial.
8. Seat function selector.
9. LEDs indicating the selected seat and back adjustment function.
10. Right indicator button and position light switch.
11. Drive control lever - Joystick.


Illustration 3. Control panel elements.

## Clutch lever for disengaging the drive

### ATTENTION!

It is forbidden to regularly move the wheelchair by pushing it over long distances - this can damage the internal drive components.

If you ever need to push the wheelchair for a short distance, there is an option to switch the clutch lever disengaging the drive.

1. The drive disengagement lever is located at the bottom of the wheelchair between the wheels.
2. Push up the clutch lever to disengage the drive and electromagnetic brake. You can now push the wheelchair.
3. To engage the drive and brake system, switch the clutch lever downwards. Drive mode activated.

Illustration 4. Location of the drive disconnect lever.

### WARNING!

Power off the wheelchair before disconnecting or engaging the drive clutch lever. Never sit in a wheelchair that is on an incline while switching the drive clutch lever. If the wheelchair has the drive disconnected then the brakes are also not active.

- Disconnect the motor drive only on a level ground.
- Ensure the wheelchair is powered off.
- Take a sit next to the wheelchair when engaging or disengaging the drive clutch lever.
- Never sit in the wheelchair while engaging or disengaging the drive.
- When you have finished pushing, engage the drive to lock the brakes.

## Power button / status indicator

To switch on the powered wheelchair, press the power button. The power button is the only button that can activate the wheelchair system.

- If there is no fault in the system, the status indicator (using the power button) lights up green.
- If a fault occurs when switching on the system, the status indicator will signal the fault with a series of red flashes - contact your distributor.

## Battery charge indicator

When the power button is pressed, power is switched on and the charge indicator shows the power level of the batteries. For more information about the battery charge indicator, look at chapter "Batteries and Charging."

## Speed selection switch

The STAR Orange wheelchair has a selectable maximum speed.

- Press the top of the button to select a higher maximum speed for the wheelchair.
- Press the lower part of the button to select a lower speed value.

The selected speed is displayed on the indicator to the right of the button.


Illustration 5. Speed selection.

### WARNING!

When the power is switched off, the control panel stores the set maximum speed for the interval. Be careful when the power is switched on again and set the speed according to the ambient conditions.

## Seat function selection button

STAR Orange wheelchair has a choice of seat functions.

- Press the upper part of the button to select the next seat function.
- Press the lower part of the button to select the previous function.

The selected function is indicated by illuminating the LED, as shown in the illustration opposite.


Illustration 6. Key position.

## Battery cut-circuit key

- Turn the key to the "ON" position to switch on the power to the wheelchair.
- Turn the key to the "OFF" position to switch off the battery power.

After driving / using the product, it is possible to remove the key in the "OUT" position to prevent unauthorized use of the wheelchair and to protect the batteries.


Illustration 7. Key position.

## Speed control dial

The dial allows you to set the maximum speed limit.

- The left-most position represents the slowest attitude.
- The right-most position represents achieving maximum speed.

## Horn

To switch on the hazard lights, press and quickly release the left-hand lighting button once.

LEDs behind the left-hand button and right-hand button will flash on and off for the duration the hazard lights are operating.

To switch off the hazard lights, press the left-hand or right-hand lighting button.

## Hazard lights

To signal the direction of travel, activate the indicator by pressing once the button for the particular direction you intend to turn. To deactivate the indicator, press again the button for the direction that is being signaled.

## Directions lights


To activate the left indicator, press and hold the left light button once.

The LED in the left button will flash.

To activate the right indicator, press and hold the right light button once.

The LED in the right button will flash.

To switch off either of the turn indicators, short press the left or right lighting button.

## Control lever for drive and wheelchair functions - joystick

### Direction control - driving.

The lever allows you to drive forwards or reverse up to the maximum speed set by the speed control knob.

- Using your hand, move the lever in a given direction to disengage the wheelchair's brakes and move on.
- To stop, release the lever completely and allow the trolley to brake completely.
- When the lever is fully released, it automatically returns to the center "stop" position, the brakes are automatically applied.

### Direction control - seat functions.

The lever allows control of the backrest recline and seat tilt with backrest functions. First select the seat function by pressing the button on the control panel.

- When the joystick is deflected forward from the neutral position, the seating motion of selected actuator will extend.
- When the joystick is deflected backward from the neutral position, the seating motion of selected actuator will retract.

# Operate

## Before you start using

- Have you fully charged the batteries? See section “Batteries and charging”.
- Is the manual lever of the drive disengagement clutch is in the drive engaged position.
- Never leave the lever in the disengaged position if you are not going to push the wheelchair manually.

## Getting on the wheelchair

### WARNING!

Never attempt to get on or off the wheelchair without first turning off the power. This prevents the wheelchair from being moved by accidentally switching the drive control lever.

1. Make sure the power is switched off.
2. Take a seat next to the wheelchair.
3. Raise the footrest platforms and rotate both footrests outwards on the wheelchair.
4. Assume a safe and comfortable position on the seat of the wheelchair.
5. Rotate the footrests to the position where they lock and then lower the footrest platforms.
6. Ensure legs are securely positioned on footrests.
7. Fasten seat belt.

## Pre-ride inspection and adjustment

- Is the seat in the correct and comfortable position? See the section “Regulation and adjustment of wheelchair components”.
- The backrest is in an upright position.
- The seat is in the standard horizontal position.
- All adjustment mechanisms are correctly locked.
- Is the seat securely locked?
- Is the power button pressed down and the status indicator light up green.
- Is the horn work correctly?
- Are there no people, animals or other obstacles along the route?
- Have you planned your route to avoid going through difficult terrain and, if possible, as few inclines as possible?

## Operating a wheelchair

### WARNING!

The following operations may have a negative effect on the steering and stability of the wheelchair:

- Keeping and walking an animal attached to a leash.
- Carrying passengers (including animals).
- Hanging an object on the armrests.
- Towing or pushing by another vehicle.

### WARNING!

At all times when using your wheelchair, keep both hands on the armrests and your feet on the footrests. This position provides the best control of the vehicle.

- Set the speed control dial to the desired speed.
- Move the drive control lever - joystick, in the desired direction.
- The automatic electromagnetic brake is disengaged and the wheelchair accelerates smoothly to the set speed.
- Move the joystick to the left to turn left.
- Move the joystick to the right to turn right.
- Push the joystick forward to drive straight ahead.
- To stop, slowly release the joystick. The electromagnetic brakes will automatically engage when your Powered Wheelchair comes to a stop.

### ATTENTION!

The wheelchair, with the same speed control dial setting, moves slower in reverse than in forward gear.

## Getting off the wheelchair

### WARNING!

Never attempt to get on or off the wheelchair without first turning the power off. This prevents the wheelchair from being moved by accidentally switching the drive control lever

1. Stop the wheelchair completely.
2. Press down the power button to power off.
3. Loosen the safety belt.
4. Swivel the footrests outwards
5. Carefully and safely get off the seat and take a seat next to the wheelchair.
6. Rotate both footrests inwards until they restore to the original positions.

# Regulation and adjustment of wheelchair components

## Headrest adjustment

1. Turn the adjusting screw **1** counter-clockwise.
2. Loosen the headrest height adjustment profile.
3. Pull the headrest up or down to set the desired height. Pay attention to the gap between mounting profile to avoid potential pinching or crushing of hands.
4. Tighten the adjustment screw **1** clockwise to lock the height adjustment profile.
5. Turn the adjusting screw **2 3 4** counterclockwise.
6. Adjust the headrest to a comfortable position.
7. Turn the adjusting screw **2 3 4** clockwise to lock it in place.


Illustration 8. Headrest adjustment.

## Backrest angle adjustment

1. Press the seat function selector.

2. The selected seat adjustment function will be displayed on the indicator.  

3. Pull or push the joystick forward or backward from the neutral position to activate the backrest tilt.

Illustration 9. Backrest adjustment.

**WARNING!**

During the backrest angle adjustment process, take care not to put your hands in the gap between the backrest and the seat plate.

## Seat and backrest angle adjustment

1. Press the seat function selector.

2. The selected seat adjustment function will be displayed on the indicator.

3. Pull or push the joystick forward or backward from the neutral position to activate the seat tilt including the backrest.

Illustration 10. Seat and backrest adjustment.

**WARNING!**

During the process of adjusting the angle of the seat with backrest, take care not to put your hands in the gap between the frame of the wheelchair and the seat plate.

**WARNING!**

Do not drive with the backrest or seat inclined.

Do not drive the electric wheelchair when the backrest or seat is in a tilted position. Riding in this position may cause:

- Loss of stability and tipping of the wheelchair.
- Uneven weight distribution.
- Change in the centre of gravity of the vehicle.
- Increased risk of accident.

Failure to heed this warning could lead to serious personal injury.

## Control panel position adjustment

1. Loosen the controller profile locking screw by unscrewing counterclockwise.
2. Pull the profile of the panel forwards or backwards to the desired position. Pay attention to the gap between the profile and the mounting rail to avoid potential pinching or crushing.
3. Turn the screw clockwise to tighten and lock the control panel profile.

Illustration 11. Control panel position adjustment.

## Armrests width and height adjustment

Adjusting the armrest width:

1. Unscrew the two armrest profile locking screws located underneath the seat plate.
2. Adjust the desired spacing width by extending or retracting the armrest profile, ensuring that the two screws can lock the armrest profile in place .
3. Lock by tightening both screws firmly.
4. Repeat steps 1-3 for the armrest on the opposite side.
5. Check that the armrests are at the desired width.


Illustration 12. Armrest spacing adjustment.

Adjusting the armrest height:

1. Unscrew the two screws locking the armrest profiles located on the side edge of the armrests.
2. Set the desired height by extending or retracting the armrest profile, ensuring that the two screws can lock the armrest profile in place.
3. Lock by tightening both screws firmly (do not leave the handle in the outer/horizontal position - this risks damaging the handle when moving in tight spaces).
4. Repeat steps 1-3 for the armrest on the opposite side.
5. Check that the armrests are at the desired height.

Illustration 13. Height adjustment of armrests.

## Footrest adjustment

Before using the new wheelchair for the first time, extend and adjust the footrests so that the lowered footrest platforms do not hit the front wheels when moving.

Footrest position adjustment:

1. Loosen the two profile locking screws with a 5 mm Allen key (illustration opposite).
2. Adjust by extending or retracting the footrest profile.
3. Lock by tightening both screws firmly.
4. Repeat steps 1-3 for the footrest on the opposite side.
5. Check that the footrests are in the desired position.

Illustration 14. Footrest position adjustment.

## Footrest with mechanical height adjustment

### WARNING!

When adjusting the height, do not put your hands in the adjustment mechanism. Ask an assistant to secure the footrests to prevent them from accidentally falling. Be careful that the footrest does not fall under its own weight. Risk of injury!

Adjusting the angle of the footrests:

1. Before raising/lowering the footrest, it must be relieved of weight; to do this, the accompanying person should lift it temporarily..
2. Then release the lock lever. **1**
3. Ask your companion to slowly raise/ lower the footrest to the desired height.

Illustration 15. Adjustment of the angle of the footrests.

4. After adjusting the footrest, move the lock lever back into position. **1**
5. Ensure that the footrest is properly locked.

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Height-adjustable footrest tilt:

1. To tilt the footrest to the side, press the lock lever shown in the photo opposite.
2. The lock mechanism works automatically and you will hear a click when you turn the footrest to its starting position.

Illustration 16. Tilting footrests.

## Removing the footrests


Illustration 18. Removing the footrests.

Removing the footrests:

1. Tilt the footrest platform upwards.
2. By pulling the locking lever, tilt the footrest outwards from the wheelchair frame.
3. Remove the footrest by lifting it from the wheelchair frame profile.
4. Repeat steps 1 to 3 for the other footrest.

## Footrest length adjustment

To adjust the length of the footrests:

1. Loosen by unscrewing the footrest profile locking screw (illustration opposite).
2. Adjust the length by extending or retracting the footrest profile (check that the maximum extension will allow safe movement of the wheelchair).
3. Tighten the footrest profile locking screw.
4. Repeat steps 1 to 3 for the other footrest.


Illustration 17. Footrest length adjustment.

## Disassembly and assembly the wheelchair

### Disassembly

The power wheelchair can be disassembled into several parts: seat cushion, backrest cushion, headrest, wheelchair frame, footrests, control panel. No tools are required to dismantle the power wheelchair. Always dismantle or assemble components on a level, dry surface with sufficient space to work and move 1.5 meters in all directions.

## Disassembly procedures

### WARNING!

Do not lift weights beyond your physical capabilities. When disassembling and assembling the wheelchair, if necessary, ask for assistance.

#### Identification key:

1. Seat cushion
2. Backrest cushion
3. Headrest
4. Wheelchair frame
5. Footrests
6. Control panel

Illustration 19. Removal of wheelchair components.

1. The seat cushion is attached to the seat frame with Velcro straps. Remove the seat cushion by lifting it upwards.
2. The backrest cushion is attached to the seat frame with Velcro straps. Remove the backrest cushion by pulling it forwards.
3. Remove the headrest by loosening the height adjustment screw located at the headrest profile handle on the back of the backrest.
4. Dismount the footrests according to the instructions in the section "Removing the footrests".
5. Remove the control panel by loosening the adjustment screw located on the armrest frame.

## Assembly

1. Reinstall the seat cushion.
2. Reinstall the back cushion.
3. Reinstall the headrest and fully secure the headrest adjustment screw.
4. Reinstall the footrests and fully secure the adjustment screws.
5. Reinstall the control panel and fully secure the adjustment screws.

### WARNING!

Be especially careful when unfolding and folding the wheelchair, lack of attention may cause body parts to be jammed by moving and folding parts!

## Basic troubleshooting

### WARNING!

If the battery symbol is flashing on the battery level indicator, the wheelchair has detected a fault. If this happens, follow the procedures below:

1. Release the drive control lever, then switch the wheelchair off and on again. If this does not clear the fault, proceed below to step 2.
2. Switch off the wheelchair and fully charge the battery. If this does not clear the fault, go below to step 3.
3. Count the number of flashes. Power switch will flash several times before an extended pause. Count the number of flashes between breaks and then pass this on to your dealer. The following fault codes are listed below:

Number of flashes	Fault description	Meaning and procedure
1	Battery low	Batteries are discharged • Recharge the batteries
2	Low battery fault	Batteries are discharged. Recharge the batteries. Check the batteries and their associated wires and connections.
3	High battery fault	The battery voltage is too high. This could be the result of overcharging and/or moving down the slope. When moving down the slope, reduce speed to minimize regenerative charging.
4	Current Limit Time-out or Controller too hot	Motor exceeds its maximum rated current for too long. • The wheelchair may have stopped. Turn the power off, wait a few minutes, then turn the power back on. • The motor may be defective. Check the motor and the cables and their connections.
5	Failure of the parking brake	The parking brake release switch is active or the parking brake is faulty. • Check the parking brake cables and their connections. • Ensure any associated switches are in their correct positions.

6	Drive lock	The stop function is active or a charger locked or OONAP* conditions have occurred. <ul style="list-style-type: none"> <li>• Release stop condition (e.g. seat raised)</li> <li>• Disconnect the battery charger.</li> <li>• Ensure that the drive control lever is in the neutral position when starting.</li> <li>• The drive control lever may need to be calibrated.</li> </ul>
7	Failure of the speed potentiometer	The drive control lever, speed potentiometer, SRW or wires and their connections may have failed. <ul style="list-style-type: none"> <li>• Check the drive control lever, speed potentiometer and wires and their connections.</li> </ul>
8	Motor voltage fault	Motor or wiring and connections have failed. <ul style="list-style-type: none"> <li>• Check the motor and wiring and its connections.</li> </ul>
9	other error	The controller has internal failure. <ul style="list-style-type: none"> <li>• Check all connections and wiring.</li> </ul>

\* The full name of OONAPU in item 6 is Out Of Neutral At Power Up.

### ATTENTION!

If questions or concerns arise about the performance of the product, please contact your distributor.

## Batteries and charging

Before charging the batteries, read the instructions and prepare the charger. The battery status indicator on the control panel indicates the approximate power of the batteries by means of 6 LEDs. All lit LEDs indicate fully charged, and only 1 LED lit indicates an immediate need to charge the batteries.. To determine this exactly, check the battery charge indicator while driving on dry, level ground at top speed.

### First charge

The wheelchair requires two long, full charge cycles (14 hours), sealed and maintenance-free batteries.

Batteries are charged by the charging system outside the vehicle. Before first use, first charge the wheelchair batteries for 8 to 14 hours dependent on discharge condition. Keep the batteries fully charged so that the wheelchair can run smoothly.

## Connecting the charger

The wheelchair can be charged via the charging socket on the bottom of the control panel. Connect the charger to the XLR socket and then plug into a grounded electrical outlet. When the batteries are fully charged, disconnect the charger first from the socket and then from the wheelchair.


Illustration 20. Connecting the battery charger.

## Operating temperatures

The wheelchair is designed for use, ideally in temperatures between 5°C~41°C. In extremely cold conditions, the batteries may freeze, depending on charge, use and composition. In extremely high temperatures, the wheelchair may run at a lower speed. This is due to a built-in mechanism to protect against damage to the gearbox and other electronic components when operating in hot conditions.

## Battery operation and charging

If you use your electric wheelchair on a daily basis, charge the batteries as soon as possible after you have finished using it for the day. The wheelchair will be ready every morning. When using the wheelchair on a daily basis, we recommend charging for 8 to 14 hours. We also recommend that you continue to charge for 4 hours after the charger has indicated the end of the cycle. If you use your wheelchair once a week or less frequently, charge the batteries at least once a week for 24 hours.

## Battery charging process

The LED light on the charger indicates the different charging status.

Red light - means that, it is in the process of charging.

Green light - means that, bulk charging (95% charged) is complete and the charger can be disconnected.

### ATTENTION!

Keep the batteries charged and avoid discharging them completely. To keep the batteries in good condition and to extend their life, we recommend charging a minimum for 48 hours, once a month.

### ATTENTION!

There is a charger in the wheelchair which blocks its functions. The wheelchair will not start and the battery charge indicator will not be active when charging the battery pack fitted to the wheelchair.

### WARNING!

- Removal of the ground pins may cause an electrical hazard. If necessary, when accessing a 2-pin socket, use a suitable, approved 3-pin adapter.
- Never use an extension cord to plug in the charger. Plug the charger directly into a properly installed standard grounded electrical outlet.
- It is prohibited, while charging the batteries, for children to play unsupervised near the wheelchair. We recommend that you do not charge the wheelchair batteries if the wheelchair is occupied by user.
- Explosive gases may be generated when charging the batteries. Keep the wheelchair and charger away from sources that can cause ignition such as fire, sparks and ensure proper ventilation when charging the batteries.
- Charge the wheelchair batteries with the supplied external charger. Do not use any other type of charger.
- Before each use, check the charger, cables and connectors for damage. If you find any damage, contact your supplier.
- Do not open the battery charger case. If the charger does not work correctly, contact your authorized supplier.
- Be aware that during charging, the battery charger case may be hot. Avoid contact with skin and other surfaces that may be damaged by heat.
- Do not leave the charger switched on on the seat of the pushchair. The charger generates heat during the charging process which can damage the upholstery.

**WARNING!**

If the charger has not been tested and approved for outdoor use then it should not be exposed to extreme weather conditions. If the battery charger is exposed to bad or extreme weather conditions, before using it indoors, leave it unplugged to allow it to adapt to differences in ambient conditions.

## Reporting technical defects and repairing the device

If a defect or damage to the product is discovered during use or maintenance, it must be reported immediately to the distributor or manufacturer.

**WARNING!**

Under no circumstances should an unqualified person repair or alter the product or its design.

The only person allowed to make any repairs or adjustments is an employee authorised by the manufacturer.

## Care and maintenance

Your wheelchair requires minimal care and maintenance. If you are not sure of your care and maintenance capabilities listed below, you can have the inspection and service carried out by an authorized supplier. The components listed below require periodic inspection and maintenance.

### Tires

Regularly check your electric wheelchair's tires for signs of damage or wear.

### Wheel Replacement

The wheelchair is equipped with polyurethane (PU) tires. Replace the entire wheel if the tire has become worn or damaged. For information on replacing wheels, contact your authorized electric wheelchair supplier.

**WARNING!**

The wheels on the wheelchair must be serviced/replaced by a qualified technician. During this procedure, make sure that, the power is switched off and the wheelchair does not have the clutch lever switched disengaging the drive.

### External surfaces

Bumpers, tires and rims can be occasionally cleaned with rubber and vinyl cleaner.

## Cleaning

Use a damp cloth and a mild, non-abrasive cleaner to clean the plastic and metal parts of the electric wheelchair. Avoid using products that can scratch the wheelchair's surfaces.

## Disinfection

If necessary, use approved disinfectants. Before using, make sure the disinfectant is safe for the wheelchair's surface. Only tested and approved disinfectants may be used for disinfection. Information on recognized and approved disinfectants and disinfection methods is available from your national health authority.

### WARNING!

Follow the safety instructions provided by the manufacturer of cleaning and disinfecting products. Not complying may result in skin irritation or early damage to the upholstery or electric wheelchair.

## Battery terminal connections

Make sure the connections remain tight and non-corroded.

## Electric wiring harnesses

- Regularly check wire connections.
- Regularly check the insulation of the wires, including those from the charger, for wear or damage.
- Have an authorized supplier repair or replace any damaged connectors, connections or insulation before using your electric wheelchair again.

## Daily checks

- With the power off, check the operation of the drive control lever. Make sure that the lever is not bent or damaged and that when released, it returns to the neutral position. Self-repair is prohibited. In case of problems, contact an authorized supplier.
- Visually inspect the control panel. Make sure it is not burnt, cut or wires sticking out. Contact your authorized supplier if there are problems.
- Check for worn (flat) spots on the tires. Worn tires can affect stability.
- Check armrests for loose, strained or damaged parts. If there are problems, contact your authorized supplier.
- Check brakes. Checks should be made on a level surface without obstacles within 1 meter of the wheelchair.

## **Weekly checks**

- Check the controller and charger connections for corrosion. If necessary, contact an authorized supplier.
- Check the pressure if there is a pneumatic tire on the equipment. If the tire does not hold pressure, contact an authorized supplier to replace the inner tube.

## **Monthly checks**

- Check that the anti-tip wheels do not rub against the ground when operating the electric wheelchair.
- Check the anti-tip wheels for wear and tear. Replace if necessary.
- Check for tire wear. Contact an authorized supplier for repairs.
- Keep the wheelchair clean and away from foreign objects such as mud, sand, hair, food, drinks, etc.

## **Yearly checks**

- Take your electric wheelchair to an authorized supplier for an yearly inspection, especially if you use your wheelchair daily. This will ensure the correct operation of the wheelchair and help avoid future breakdowns.

## **ABS plastic covers**

If your electric wheelchair has a shiny chassis surface then it has been treated with a clear fixative coating. You can apply a small coat of car wax to help restore the glossy surface. If your wheelchair's chassis has a matte finish, use only products designed for matte finish surfaces. Do not use wax or other products designed for glossy paintwork.

## **Axle bearings and motor/transaxle**

All these components are pre-lubricated, sealed and do not require further lubrication.

## **Console, Charger, and Rear Electronics**

Keep these areas free of moisture. Allow these areas to dry thoroughly if they have been exposed to moisture before operating your electric wheelchair again.

## **Nylon Lock Nut Replacement**

During periodic maintenance, assembly or disassembly, each securing nylon nut must be replaced with a new one. Nylon nuts must not be reused as this may damage the nylon threads and reduce the effectiveness of the fit. Replacement nylon nuts are available from local hardware stores or an authorized supplier.

## **Storage of the wheelchair**

If you plan not to use the wheelchair for a period of time then follow these rules:

- Fully charge batteries before storage.
- Switch the wheelchair battery off.
- Store the wheelchair in a warm, dry room.
- Avoid storing in a place where the wheelchair will be exposed to temperature changes.
- Recommended storage temperatures: from 5°C to +40°C.
- Recommended ambient humidity for storage: 30% to 70%.

Batteries that are regularly completely discharged, not regularly charged, stored at inappropriate temperatures, or stored without being charged may be permanently damaged resulting in unreliable operation and shortened battery life. To ensure proper operation, it is recommended that batteries be charged regularly, throughout the period of prolonged storage of the wheelchair.

## **Re-use (reusable product)**

The product is suitable for reuse. The number of times it can be reused depends on how often and in which way the product is used. Before reuse, thoroughly clean and disinfect the product and check its safety features. If necessary, take repair action. If necessary, take repair action.

## **Product life cycle**

The expected life of a medical device is 5 years, as long as it is used as intended and all maintenance and servicing is performed. The lifespan of a medical device also depends on the frequency of use, the environment in which it is used and its care. By using spare parts, the life of the product can be extended. Spare parts are generally still available for 5 years after the model is discontinued. This expected lifetime of a medical device does not imply an additional guarantee.

## **Replacement parts**

The availability of replacement components and/or spare parts is dependent on the capabilities of the supplier, not the manufacturer. For more information on replacement components and/or spare parts, please contact your authorized supplier.

## **Disposal of an electric wheelchair**

The electric wheelchair must be disposed of in accordance with applicable local and national regulations. For information on the proper disposal of packaging, metal frame parts, plastic parts, electronics, batteries, neoprene, silicone and polyurethane materials, please contact your local waste disposal company or authorized supplier.

## **Warranty**

Details of the warranty can be found in the attached “warranty card” document.

### **ATTENTION!**

The warranty does not cover deterioration of battery performance caused by being left in a fully discharged state, in cold conditions for an extended period of time, or by excessive use.

## Technical parameters

Product name	STAR Power Wheelchair
Model	Orange
Classification according to the standard EN 12184:2014	Class B
Basic UDI-DI	59017804DRVI0NNT
Maximum user weight	136 kg
Maximum permissible total weight	242 kg
Maximum speed (driving forward)	10 km/h
Maximum speed (driving backward)	5,3 km/h
Minimum braking distance at max. speed	1500 mm
Energy range	30 km
Maximum obstacle height	45 mm
Maximum safe inclination	10° (17,63%)
Minimum turning radius	900 mm
Minimum width when reversing	900 mm
Dimensions without seat (L x W x H)	1080 x 590 x 1215 mm
Carton size	1020 x 725 x 850 mm
Total mass with / without batteries	106,5 kg / 76,5 kg
Battery weight	30 kg
Effective seat depth	450 mm
Effective seat width	430-480 mm
Backrest height	520 mm
Seat surface height at front edge	75 mm
Footrest length adjustment range	130-200 mm
Seat surface angle	0°-20°
Backrest angle	90°-140°
Overall length of the wheelchair (with footrest)	1155-1820 mm
Armrest to seat distance	170 - 290 mm
Footrest to seat distance	450 - 575 mm

Motor	2x340 W, electromagnetic brake
Batteries	52 Ah / 24 V Lead-acid-gel batteries
Controller	Dynamic controller Drive+Light+Actuators
Protection class	IPX4
Battery charging	Wyjście DC 24 V , 6 A
Protection class of battery charger	IPX1
Insulation class of battery charger	Class 2
Maximum motor sound level	68 dB
EMC compatible according to	ISO 7176-21
Ignitability of upholstery according to	ISO16840-10
Diameter rear wheels	14 x 2.6 cali (356 x 66 mm) (Air)
Diameter front wheels	9 x 2.6 cali (229 x 66 mm) (Air)
Suspension	Rear
Required doorway entry depth	900 mm
Control	Lever to start/stop movement
Operating console	Rotatable switch for speed control
Storage and use temperature	from +5°C to +40°C
Humidity to store and use	30% - 70%
Safety	Electromagnetic brake, safety belt
Frame material	steel
Other materials (upholstery, cushion)	Oxford cloth, foam in compliance with ISO 16840-10:2021
<p>We reserve the right to introduce technical changes. The tolerance in the measurements is <math>\pm 15</math> mm /1.5 kg / 1.5°. The theoretical turning radius decreases if the wheelchair is often used on inclines, uneven surfaces or for driving up curbs. The maximum distance covered has been tested under ideal conditions in accordance with ISO7176-4.</p> <p>Speed and range depends on the user's weight, type of terrain, charge and battery condition. Information is current at the time of publication.</p>	

## Product safety symbols

The symbols used below on the electric wheelchair are used to identify warnings, mandatory actions and prohibited actions. Be sure to read everything and understand their meaning.

It does not comply with the requirements of the ISO 7176-19 standard for transporting people in a wheelchair located in a motor vehicle. It is prohibited to sit in a wheelchair while traveling in a motor vehicle.

Before use, charge the battery to the maximum level. When the wheelchair is not in use - switch off the power supply.

## Symbols

<b>MD</b>	Medical device
	Date of manufacture
	Manufacturer
<b>EU</b>   <b>REP</b>	Authorized representative in the European Union
<b>CE</b>	Conformité Européene (compliance with the European Community)
<b>REF</b>	Catalogue number
<b>SN</b>	Serial number
<b>LOT</b>	Lot number
<b>UDI</b>	Unique Device Identifier

	Model number
	Follow the instructions for use
	Refer to the instructions for use
	Indicates the maximum weight capacity. Do not exceed the specified weight capacity of the product
	Warning
	Indicates the danger of electric shock
	EMI-RFI - This product has been tested and passed at an immunity level of 20 V/m
	Indicates the risk of explosion
	Corrosive chemicals contained in batteries
	Indicates a mandatory action that should be performed as specified in order to avoid a hazardous condition/situation. Failure to perform mandatory actions can lead to personal injury and/or equipment damage
	Protect from humidity
	Possibility of jamming fingers
	Pinch/Crush danger points
	Permitted temperature

	For indoor/outdoor use
	Contains Lead
	Maximum user weight
	Maximum permissible total weight
	Maximum safe inclination
	Indicates attachment points on the wheelchair
	The product is designed to be transported by car
	The product is not intended to be used as a seat for transport in a motor vehicle
	Do not stand on the footrest
	Use only AGM or Gel-Cell batteries to reduce the risk of leakage or explosive conditions
	Wear eye protection
	N = Neutral (Brakes disengaged; lever disengaging the drive turned on)
	D = Drive (Brakes engaged; unit in Drive Mode)
	Indicates a prohibited action that should not be performed at any time or in any circumstances. Performing a prohibited action can cause personal injury and/or equipment damage

	Do not use batteries with different capacity (Ah). Do not mix old with new batteries. Always replace both batteries at the same time
	Do not mix AGM and Gel-Cell batteries. Always use batteries of the same type and chemistry
	Keep tools and other metal objects away from battery terminals. Contact with tools can cause electric shock
	It is forbidden for children to play, unsupervised, while charging the wheelchair's batteries
	Removal of the grounding pin may cause an electrical hazard. If necessary, properly install an approved 3-pin adapter in an electrical outlet with access to a 2-pin plug
	Do not connect an extension cord to the AC/DC converter or the battery charger
	Indicates flammable material. Do not expose to heat sources such as open flame or sparks. Do not transport batteries with flammable or combustible items
	Avoid exposure to rain, snow, ice, salt, or standing water whenever possible. Maintain and store in a clean and dry condition
	Do not remove anti-tip wheels
	Keep your hands away from the tires when driving. Be aware that loose fitting clothing can become caught in the drive tires
	Do not operate the drive disengagement clutch on an incline
	Do not operate the drive disengagement clutch on a decline
	Contact your local recycling center or authorized provider for information on proper recycling of product components

## SELECTIVE COLLECTION OF ELECTRICAL AND ELECTRONIC EQUIPMENT.

This symbol on the products or accompanying documents indicates that used electrical or electronic products and batteries must not be disposed of in municipal waste. For proper disposal, restoration or recycling, such products should be returned to collection points for this type of waste. Alternatively, in some countries of the European Union or other European countries, you can return your products to a local seller at the time of purchase of a similar new product. By correctly disposing of this product, you will help preserve valuable natural resources and support the prevention of potential negative impacts on the environment and human health, which could result from incorrect waste disposal. More information is available from the communal offices or waste collection sites. In case of incorrect disposal of this product, penalties may be imposed in accordance with local regulations. If you want to dispose of an electrical or electronic device, get the information you need from your dealer or supplier.

**Disposal in countries outside the European Union:** This symbol is valid in the European Union. If you want to dispose of this device, please obtain the necessary information on the correct disposal at your local authorities or from your dealer.

## Contact details

### Manufacturer

mdh Sp. z o.o.

Address: ul. Maratońska 104, 94-007 Łódź, Poland

tel. +48 42 674 83 84, fax. +48 42 636 52 21

[www.mdh.pl](http://www.mdh.pl)

[www.viteacare.com](http://www.viteacare.com)

We reserve the right to make technical changes due to technical progress.