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
**English version** on next page

**Users with visual impairments can find the PDF files together with further information on our website:**

**<[www.viteacare.com](http://www.viteacare.com)>.**

**👉 Contact your specialist dealer when required.**

**Alternatively users with visual impairments can have the documentation read out by a helper.**

# Table of contents

<b>Introduction</b> .....	<b>50</b>
<b>Product characteristics</b> .....	<b>51</b>
Intended use / indications.....	51
Contraindications .....	51
<b>Safety instructions</b> .....	<b>52</b>
General safety rules .....	53
Modifications .....	53
Security check.....	54
Capacity .....	54
Transport the wheelchair.....	54
Inclination guidance.....	55
Driving outdoors and precautions in bad weather conditions .....	55
Physical limitations .....	56
Information on electromagnetic compatibility (EMC) .....	56
<b>Product information</b> .....	<b>57</b>
Electrical circuit diagram.....	58
Control panel (Joystick).....	58
Clutch lever for disengaging the drive .....	59
Power button / Battery charge indicator .....	59
Speed selection buttons .....	60
Horn.....	60
Control lever for drive and wheelchair functions - joystick.....	60
Charging Port.....	60
<b>Installation and adjustment</b> .....	<b>61</b>
Mounting and removing control panel .....	61
Installation and removal of the battery.....	62
Fastening and unfastening seat belts .....	62
Anti-tip wheels .....	63
Armrests - reclining.....	63
Unfolding the wheelchair .....	64
Folding the wheelchair.....	65
Installation 24-inch rear wheels .....	66
Removal of 24-inch rear wheels .....	66
<b>Operate</b> .....	<b>67</b>
Before you start using .....	67
Operating a wheelchair - electric drive mode.....	67
Using the wheelchair - manual wheelchair mode.....	68
<b>Basic troubleshooting</b> .....	<b>74</b>
<b>Batteries and charging</b> .....	<b>76</b>
Precautions for Li-ion Battery.....	76
First charge.....	76

Connecting the charger .....	77
Operating temperatures .....	77
Battery operation and charging .....	77
The process of charging batteries outside the wheelchair .....	78
<b>Reporting technical defects and repairing the device .....</b>	<b>79</b>
<b>Care and maintenance .....</b>	<b>79</b>
Tires.....	79
Wheel Replacement.....	79
External surfaces .....	79
Cleaning.....	80
Disinfection.....	80
Battery terminal connections.....	80
Electric wiring harnesses .....	80
Daily checks .....	80
Weekly checks.....	81
Monthly checks .....	81
Yearly checks .....	81
Plastic covers .....	81
Console, Charger, and Rear Electronics .....	81
Storage of the wheelchair .....	82
Re-use (reusable product) .....	82
Product life cycle.....	82
Replacement parts .....	82
Disposal of an electric wheelchair.....	82
<b>Warranty.....</b>	<b>83</b>
<b>Technical parameters .....</b>	<b>84</b>
<b>Symbols .....</b>	<b>86</b>
<b>Contact details .....</b>	<b>88</b>

## Introduction

The TWIST wheelchair has been designed and manufactured to the highest quality standards. Before using this product for the first time, please read these operating instructions. Improper use of the wheelchair may result in personal injury or traffic accidents. If you still have questions after reading this manual, please contact your specialist dealer.

Ensure that the 'Warranty Card' document is completed and enclosed and handed over to the distributor, as this will enable claims to be processed efficiently.

# ATTENTION

Always follow the basic precautions according to the information given below regarding safety when using the product. Read the entire manual carefully before use. Keep this manual for future reference.

### WARNING!

Improper use can lead to death, serious injury, unsafe conditions or product malfunction.

### ATTENTION!

Improper use may lead to personal injury and/or damage to the wheelchair.

### SUGGESTION!

Follow the instructions below to keep your wheelchair in good working order.

If the electric wheelchair is used by another person, you must provide them with copies of these instructions for use.

These instructions for use relate to the design and specifications of the product at the time of publication. Due to design changes, some illustrations and photos in the manual may not correspond to the wheelchair you have purchased. We reserve the right to make design changes. The product shipped in its original packaging requires partial assembly - follow the information in the "Disassembly and assembly the wheelchair" section. The manufacturer shall not be responsible for damage resulting from improper or unsafe use of the medical device, if such use was against the instructions and recommendations contained in the product documentation. The manufacturer shall ensure that the product has been designed and manufactured in accordance with the applicable safety and quality standards.

## **Product characteristics**

TWIST is an electric wheelchair designed to transport disabled people both indoors and outdoors. The wheelchair can be used in both electric and manual mode when fitted with the appropriate 24-inch wheels.

**Construction:** The wheelchair has a folding mechanism for easy storage and transport. When folded, it can be placed upright to save space. A unique feature is the ability to switch between electric and manual modes:

- Electric mode with joystick and motor drive.
- Manual mode with 24-inch self-propelled wheels.

**Drive:** The vehicle can be powered by an electric motor

**Energy range:** On a single charge it can be up to several kilometers, depending on battery capacity and conditions of use.

**Safety:**

- Anti-tipping wheels to prevent the wheelchair from tipping over.
- Seat belts with "click" fastening system.
- Raised armrests for easy entry and exit.

## **Intended use / indications**

The product is intended for disabled people with musculoskeletal conditions to increase their mobility.

## **Contraindications**

Contraindications to upright standing. Injury or paresis of the upper limbs preventing the use of a power wheelchair.

## Safety instructions

1. Please note that the device may be either a powered (electric) wheelchair or a manual wheelchair and is intended for use both indoors and outdoors. Special care should be taken when driving around bends, obstacles and inclines.
2. Always ensure that the electric wheelchair is fully assembled before getting on.
3. Read and follow the information in the instructions for use. It is forbidden to use this product and the available optional equipment without first carefully reading and understanding this manual. If you do not understand the warnings, cautions or instructions, contact your supplier, otherwise personal injury or accident may result. Never attempt to use the electric wheelchair beyond its limitations as described in this manual.
4. It should not be used by people who clearly suffer from physical or mental limitations that prevent them from safely driving an electric wheelchair in traffic.
5. Place your feet on the footrest when driving. Do not stand on the footrest.
6. To operate an electric wheelchair safely, ride on flat surfaces where all four wheels make sufficient contact with the ground.
7. Keep metal objects away from the battery connections. Risk of electric shock.
8. When negotiating any corner, reduce speed and keep the centre of gravity stable.
9. When reversing your electric wheelchair, make sure there are no obstacles behind you.
10. It is prohibited for children to play near or operate the electric wheelchair.
11. Some parts of the power wheelchair are susceptible to temperature changes. The controller can only operate at temperatures between  $-5^{\circ}\text{C}$  to  $+40^{\circ}\text{C}$ .
12. In extremely cold temperatures, the batteries may freeze and the electric wheelchair may not operate. In extremely high temperatures, the wheelchair may run at a lower speed due to the controller's heat protection function, which prevents damage to the motors and other electrical components.
13. Store device in clean and dry conditions.
14. It is forbidden to drive at night without lighting.
15. Do not remove the anti-tip wheels or transit wheels.
16. Keep hands and feet away from moving parts when driving. Take care that loose clothing does not get caught in the drive wheels.
17. Do not connect the battery charger to an extension cord.
18. Surface temperatures of the wheelchair can increase when exposed to external sources of heat.
19. Disassembly of the controller, motor or charger by anyone other than an authorised technician is prohibited and will invalidate any applicable warranty conditions.
20. It is forbidden to drive the wheelchair on expressways and motorways.
21. It is forbidden to use an power wheelchair under the influence of alcohol, medication or any other substance affecting the ability to drive.
22. Be careful when driving in busy areas and shopping centers.
23. Under no circumstances should the wheelchair be used as a seat when being transported by motor vehicles.

24. Do not attempt to lift the wheelchair by gripping anything other than the frame (for example, the seat, armrest or chassis components).
25. Please note that the temperature of the components and frame parts may increase when exposed to the sun or other heat source, while at low temperatures there is a danger of very high cooling of the metal parts of the wheelchair. This can cause burns or frostbite on parts of the body.
26. To avoid the risk of electric shock, the device must only be connected to a mains supply with protective grounding.

**WARNING!**

The product is intended for use by people with disabilities and movement disorders. If you are not sure whether your health condition allows you to use the product, consult a healthcare professional.

**WARNING!**

Any “serious incident” related to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides.

## **General safety rules**

Read the instructions for use carefully with understanding before using the electric wheelchair.

Users may encounter difficulties with maneuvering when passing through narrow doorways, ascending and descending ramps, driving around corners and navigating uneven terrain. Reduce speed, take your time and maneuver the electric wheelchair carefully.

## **Modifications**

Do not modify, remove, disable or add any parts, components or features to your electric wheelchair. The manufacturer offers a wide range of electric wheelchair accessories to meet your needs.

**WARNING!**

Do not modify the electric wheelchair in any way - this will void the warranty.

## Security check

### ATTENTION!

Check the following before use:

- Wheelchair frame - is there any damage
- Electrical connections and wiring harnesses - check for damage
- Batteries - that they are fully charged and the connections are free of corrosion.
- Brakes - check that the electromagnetic brake is functioning.
- Condition of tires - make sure they are not damaged or excessively worn - if there is a problem, contact your dealer.

## Capacity

TWIST wheelchair is designed and built to transport one person with a maximum weight of 125 kg. The wheelchair is not designed to carry goods or objects other than those previously described.

### Warning!

Exceeding the maximum user weight will void the warranty. The manufacturer is not responsible for injury or damage caused by failure to comply with the maximum user weight.

## Transport the wheelchair

### Warning!

Do not sit in an electric wheelchair while it is being transported in a motor vehicle.

### Warning!

Disconnect the power supply to the battery for transport by setting the button on its cover to the "OFF" position.

When transporting the electric wheelchair by vehicle, store it safely in the back of a van, truck or in the boot of a car. Adjustable parts should be disassembled or adequately secured for transport. To ensure safety and avoid damage to the wheelchair and its components during transport, it should be carefully secured to prevent movement in the motor vehicle.

## Inclination guidance

Remain cautious when going up and down hills. If necessary, reduce your speed when going downhill. If you are moving up a hill try to keep the wheelchair moving forward. If necessary, stop. Be careful and slowly accelerate your wheelchair forward. If you are descending, reduce the speed of the electric wheelchair to the lowest setting and then proceed with caution. If the wheelchair is descending faster than you anticipated, slowly release the drive control lever until it stops and then push the lever slightly to continue descending safely.

### WARNING!

Drive the electric wheelchair forward down or up a slope, never reverse. Irregular movements or making back and forth movements may increase the probability of tipping over. Changes in slope can reduce stability. Under no circumstances should you travel up a slope with a gradient greater than the maximum stability angle of this electric wheelchair.

## Driving outdoors and precautions in bad weather conditions

Although the electric wheelchair is designed to be used indoors to help with daily living activities, the manufacturer has designed it to perform perfectly outdoors on level surfaces.

### ATTENTION!

When driving outdoors, avoid uneven terrain, soft surfaces, tall grass, loose gravel, loose sand, fresh and salt water, stream banks, lakes and oceans. If you are unsure about a surface, avoid it.

Exposure to adverse weather conditions should always be avoided. If you find yourself in bad weather conditions while using your electric wheelchair, immediately go to the nearest shelter. Dry the electric wheelchair completely before use, before charging or before storing.

## Physical limitations

Use awareness, caution, prudence and common sense when operating the wheelchair. Always be aware of your own limitations and rules of use when operating an electric wheelchair.

### WARNING!

It is forbidden to use the wheelchair while under the influence of alcohol or psychotropic drugs. Never smoke while sitting in the wheelchair. Damaged or worn upholstery must be replaced as it increases the risk of fire.

### ATTENTION!

When operating an electric wheelchair, please be aware of precautions, warnings and safety issues if you are taking medications prescribed by a doctor or available without a prescription.

## Information on electromagnetic compatibility (EMC)

### Electromagnetic Interference - Warning and Safety Recommendations

Despite compliance of the device with applicable directives and standards regarding electromagnetic compatibility (EMC), the electric-powered wheelchair may cause interference with the operation of other electrical or electronic devices, such as mobile phones, radios, alarm systems, and medical devices. At the same time, the wheelchair itself may be susceptible to the influence of strong electromagnetic fields emitted by other devices.

### Safety Recommendations:

- If you notice unstable or improper operation of the wheelchair or other devices, immediately stop the wheelchair and switch off both the wheelchair and any devices that may be causing interference.
- Move the wheelchair to another location, away from potential sources of electromagnetic interference. After relocating, restart the device and ensure it operates correctly before continuing use.
- Avoid using the wheelchair near strong sources of electromagnetic radiation, such as mobile phone base stations, large electric motors, radio stations, transmitters, or repeaters.
- If you use other medical or electronic devices simultaneously with the wheelchair, consult with a doctor or the manufacturers of those devices to confirm their safe co-use.

### Information Regarding Electromagnetic Risk Analysis:

A detailed electromagnetic risk analysis was conducted in accordance with applicable standards and guidelines. The results confirm that the electric-powered wheelchair maintains basic safety parameters and functionality under the electromagnetic compatibility conditions specified by standards. All potential sources of electromagnetic hazards have been identified, assessed, and considered during the device's design process.

### Additional Recommendations:

- Perform regular technical inspections and check of the wheelchair in accordance with manufacturer’s recommendations to ensure appropriate resistance to electro-magnetic interference.
- Users should carefully read the entire user manual and be aware of all warnings and safety instructions.
- In case of doubts regarding the operation of the wheelchair or upon detecting interference, immediately contact an authorized distribution center or the manufacturer.
- In environments with exceptionally strong electromagnetic fields (e.g., hospitals, radio stations), exercise particular caution and consult specialists responsible for medical equipment safety.

### Product information

The TWIST electric power wheelchair is designed for outdoor and indoor use. The following illustration shows all parts relevant to the user/driver. The drive disconnect clutch levers are located on the sides between the front and rear wheels.


Illustration 1. TWIST wheelchair component list.

- |                |                              |
|----------------|------------------------------|
| 1. Backrest    | 9. Rear wheel                |
| 2. Armrest     | 10. Front Wheel              |
| 3. Joystick    | 11. Battery assembly         |
| 4. Handbrake   | 12. Mudguard                 |
| 5. Seat        | 13. Motor                    |
| 6. Footrest    | 14. Anti-tip wheels          |
| 7. Push handle | 15. 24" wheel (Left + Right) |
| 8. Brake lever | 16. 24" wheel adapter        |

## Electrical circuit diagram


Illustration 2. Draft of the electrical diagram.

## Control panel (Joystick)


Illustration 3. Control panel components.

1. Battery charge indicator
2. Power button
3. Audible horn
4. Speed indicator
5. Button for decreasing speed
6. Button for increasing speed
7. Joystick
8. Charger socket connection

## Clutch lever for disengaging the drive

### ATTENTION!

It is forbidden to regularly move the wheelchair by pushing it over long distances - this can damage the internal drive components.

If you ever need to push the wheelchair for a short distance, there is an option to switch the clutch lever disengaging the drive.

1. The drive disconnection levers are located at the wheel motors.
2. Push the clutch lever down to disengage the drive and electromagnetic brake. You can now push the wheelchair.
3. To engage the drive and brake system, switch the clutch lever upwards. Drive mode activated.


Illustration 4. Location of the drive disconnect lever.

### WARNING!

Power off the wheelchair before disconnecting or engaging the drive clutch lever. Never sit in a wheelchair that is on an incline while switching the drive clutch lever. If the wheelchair has the drive disconnected then the brakes are also not active.

- Disconnect the motor drive only on a level ground.
- Ensure the wheelchair is powered off.
- Take a sit next to the wheelchair when engaging or disengaging the drive clutch lever.
- Never sit in the wheelchair while engaging or disengaging the drive.
- When you have finished pushing, engage the drive to lock the brakes.

## Power button / Battery charge indicator

Power button: the control system to activate and deactivate the wheelchair. Press the button to power on and off the wheelchair. 5 segments are displayed, 2 in green, 2 in yellow and 1 in red.

2 green segments displays: the right time to charge the battery to prolong its lifespan. 1 red light flashes, which indicate the battery is in a critical state requiring urgent recharge the battery. To continue use the wheelchair is not allowed because it will be shortened the battery lifespan once deep discharge happens.

## Speed selection buttons

The TWIST wheelchair offers a choice of maximum speed.

- Press the plus “+” symbol button to select a higher maximum speed.
- Press the minus button “-” to select a lower speed value.

The selected speed is indicated by LEDs.


Illustration 5. Speed selection.

### WARNING!

When the power is switched off, the control panel stores the set maximum speed for the interval. Be careful when the power is switched on again and set the speed according to the ambient conditions.

## Horn

This button activates the warning signal. The wheelchair must be switched on for the horn to be active. Do not hesitate to use the horn if it can prevent an accident or injury.

## Control lever for drive and wheelchair functions - joystick

### Direction control - driving.

The lever allows you to drive forwards or reverse up to the maximum speed set by the speed control knob.

- Using your hand, move the lever in a given direction to disengage the wheelchair's brakes and move on.
- To stop, release the lever completely and allow the trolley to brake completely.
- When the lever is fully released, it automatically returns to the center “stop” position, the brakes are automatically applied.


Illustration 6. Using the joystick.

## Charging Port

Positioned under the joystick (See section “Batteries and charging”).

# Installation and adjustment

## Mounting and removing control panel

1. Insert the joystick bracket into the armrest and tighten the round handle screw.
2. Once the joystick is mounted, connect the two connectors and tighten them.


Illustration 7. Joystick assembly.

## Installation and removal of the battery

1. **Removal:** Grip the handle on battery and lift it up to remove it.
2. **Installation:** Align the slot of battery with the one on the battery rack, then push downwards to secure the battery in place.

### ATTENTION!

1. When inserting the battery box, the button on the battery box should be in the "ON" state, i.e., the power is on.
2. Hold the battery with both hands to prevent dropping when to install or remove the battery.


Illustration 8. Battery assembly.

## Fastening and unfastening seat belts

1. Hold the buckle on the seat belt and press together to unbuckle the seat belt.
2. Align the buckle with the insert, press down until you hear the "click" sound.

### WARNING!

Fasten your seatbelt when sitting in a wheelchair.


Illustration 9. Use of seat belts.

## Anti-tip wheels

The wheelchair is equipped with wheels that prevent it from tipping over in an emergency. The anti-tip wheels cannot be removed.

**WARNING!**  
Do not drive the wheelchair after the anti-tip wheels are damaged.

**WARNING!**  
The anti-tip wheels are only suitable for driving on hard surfaces. If the road is muddy, the wheel may get stuck.


Illustration 10. Anti-tip wheels.

## Armrests - reclining

1. Press the red button on the side.
2. The armrest can be lifted up or down, making it easy to get in or out of the wheelchair.


Illustration 11. Reclining armrests during transfer to/from a wheelchair.

## Unfolding the wheelchair

1. Take out the wheelchair from the cardboard box and open the footrest.
2. Locate the release pin under the seat, pull it outwards.
3. Press down on both side crossbars while stabilising the frame until you hear an audible click, until the wheelchair is fully extended.
4. Once the seat is unfolded, lower both armrests.
5. Place the seat cushion on the seat.


Illustration 12. Unfolding the wheelchair.

## Folding the wheelchair

1. Remove the seat cushion and set it aside.
2. Raise both armrests.
3. Firmly pull the release lock mechanism under the seat.
4. Grasp the seat and backrest frame while pressing inwards.
5. The folding stage of the frame continues until you hear an audible click.
6. Lift the footrest platform into the folded position.
7. The wheelchair is fully folded.

### ATTENTION!

1. Folded wheelchair can be stored upright to save space.
2. To move it in the folded position, hold the handles and push steadily on the rear and support wheels without tilting it.


Illustration 13. Folding the wheelchair.

### WARNING!

1. Before opening or folding the wheelchair, make sure that there are no obstacles around the wheelchair.
2. Always ensure that your fingers are away from the mechanical joints of the folded wheelchair to prevent your hands from getting caught.
3. Do not sit on the seat or guide the wheelchair before fully opening the folded wheelchair.
4. Ensure the wheelchair is fully open before use.

## Installation 24-inch rear wheels

Press the quick release button in the middle of the 24" wheel, then insert it into the inner hole of the wheel adaptor until you hear a click and the wheel is correctly fitted.


Illustration 14. Mounting of 24-inch wheels.

## Removal of 24-inch rear wheels

Hold the wheelchair steady with one hand while pressing the quick release button with the other, then pull the 24-inch wheel outwards to remove it.


Illustration 15. Removal of 24-inch wheels.

# Operate

## Before you start using

- Have you fully charged the batteries? See section “Batteries and charging”.
- Is the manual lever of the drive disengagement clutch is in the drive engaged position.
- Never leave the lever in the disengaged position if you are not going to push the wheelchair manually.
- Is the power button pressed down and the status indicator light up green.
- Is the horn work correctly?
- Are there no people, animals or other obstacles along the route?
- Have you planned your route to avoid going through difficult terrain and, if possible, as few inclines as possible?

## Operating a wheelchair - electric drive mode

### WARNING!

The following operations may have a negative effect on the steering and stability of the wheelchair:

- Keeping and walking an animal attached to a leash.
- Carrying passengers (including animals).
- Hanging an object on the armrests.
- Towing or pushing by another vehicle.

### WARNING!

At all times when using your wheelchair, keep both hands on the armrests and your feet on the footrests. This position provides the best control of the vehicle.

## Getting on the wheelchair - electric drive mode

### WARNING!

Never attempt to get on or off the wheelchair without first turning off the power. This prevents the wheelchair from being moved by accidentally switching the drive control lever.

1. Ensure that the power supply is switched off in electric drive mode. In manual wheelchair mode, apply the brakes.
2. Take a seat next to the wheelchair.
3. Raise the footrest platform.
4. Assume a safe and comfortable position on the seat of the wheelchair.
5. Lower the footrest platform.
6. Ensure legs are securely positioned on footrests.
7. Fasten seat belt.

- Set desired speed.
- Move the drive control lever - joystick, in the desired direction.
- The automatic electromagnetic brake is disengaged and the wheelchair accelerates smoothly to the set speed.
- Move the joystick to the left to turn left.
- Move the joystick to the right to turn right.
- Push the joystick forward to drive straight ahead.
- To stop, slowly release the joystick. The electromagnetic brakes will automatically engage when your Powered Wheelchair comes to a stop.

**ATTENTION!**

The wheelchair, with the same speed control dial setting, moves slower in reverse than in forward gear.

**Getting off the wheelchair - electric drive mode**

**WARNING!**

Never attempt to get on or off the wheelchair without first turning the power off. This prevents the wheelchair from being moved by accidentally switching the drive control lever

1. Stop the wheelchair completely.
2. Press down the power button to power off.
3. Loosen the safety belt.
4. Fold the footrest platform.
5. Carefully and safely get off the seat and take a seat next to the wheelchair.

**Using the wheelchair - manual wheelchair mode**

**Getting in and out of the wheelchair**

- Place the wheelchair as close to the place to which you want to move as possible.
- Engage the parking brakes.
- Remove the armrests or pull them upwards.
- Fold the footrests platform.
- Place your feet on the ground.
- Hold on to the wheelchair and, where necessary, to a nearby object.
- Slowly move to or from the wheelchair.

Illustration 16. Getting in and out of the wheelchair.

**WARNING!**

Tipping risk.

When getting in or out of the wheelchair, there is a high risk of tipping. Only get in or out of the wheelchair when you are physically able to do so. When getting in or out, sit as far back in the seat as possible. This will prevent damage to the upholstery and stop the wheelchair from tipping forward.

**WARNING!**

If the brakes are released or damaged, the wheelchair may roll down in an uncontrolled fashion – do not support yourself on the brakes when getting in or out of the wheelchair.

**WARNING!**

Side guards and armrests may be damaged – never sit on the side guards or armrests when getting in or out of the wheelchair.

### Guiding and controlling the wheelchair

- Before moving off with the wheelchair, always make sure that no part of the patient's body or clothes may get caught on the moving parts of the wheelchair.
- If the patient cannot maintain a seated position on their own, the wheelchair should be additionally provided with stabilising belts (preventing the patient from sliding out of the seat).
- Guide the wheelchair with both hands.
- The wheelchair is guided and steered by the user using the handrims. Before moving off without a caregiver, identify the critical tipping point of the wheelchair.

Illustration 17. Identifying the critical tipping point.

**WARNING!**

Tipping risk.

The wheelchair may tip backwards if it is not provided with tipping protection (anti-tip wheels). To prevent the wheelchair from tipping, install tipping protection (anti-tip wheels).

**WARNING!**

Tipping risk.

The wheelchair may tip forward – when configuring the wheelchair, check how it behaves, when it is tipped forward and adjust your driving style accordingly.

**ATTENTION!**

A heavy load suspended from the backrest may change the centre of gravity of the wheelchair – take caution when guiding the wheelchair.

**Going down thresholds or curbs**

**Forward without assistance from a caregiver:**

- An experienced wheelchair user can negotiate thresholds or curbs without assistance.
- You can get off a low curb by moving forward. Make sure that the footrest platforms do not touch the ground.
- Keep your balance on the rear wheels to reduce the pressure on the front wheels.
- Then, negotiate the curb.

**Backwards without assistance from a caregiver:**

- An experienced wheelchair user can negotiate thresholds or curbs without assistance.
- Turn the wheelchair so that both rear wheels face the curb.
- Lean forward to shift the centre of gravity forward.
- Move the wheelchair close to the curb edge.
- Then, using the handrims, move the wheelchair off the curb in a controlled fashion.

**Forward with assistance from a caregiver:**

- Taller curbs can be crossed with assistance.
- The caregiver should tilt the wheelchair back slightly.
- Negotiate the curb moving on the rear wheels.
- Then, the wheelchair must be supported on all four wheels.


Illustration 18. Guiding the wheelchair without assistance from a caregiver.


Illustration 19. Guiding the wheelchair with assistance.

## Moving on slopes

### WARNING!

Risk of losing control over the wheelchair When moving on inclines or slopes, the wheelchair may tip backwards, forward or to the sides.

- To prevent tipping when moving up an incline, lean slightly to the front to bring the centre of gravity forward.
- When negotiating long inclines, a caregiver should always walk behind the wheelchair.
- Avoid side slopes.
- Avoid slopes over 9°.
- Avoid sharp movements when changing direction on an incline.
- If the wheelchair is standing on a slope, use parking brakes.

### Going down slopes:

- When going down slopes, make sure to control your direction and speed.
- Lean backwards and slow down the handrims with your hands.
- You should be able to stop the wheelchair at any time by grabbing the handrims.


Illustration 20. Driving down slopes without assistance from a caregiver.

### WARNING!

Risk of hand burns. Prolonged braking generates a lot of heat and friction on the handrims (particularly non-slip handrims) – wear suitable gloves.

## Going up thresholds or curbs

### Forward without assistance from a caregiver:

- An experienced wheelchair user can negotiate thresholds or curbs without assistance.
- Move up to the edge of the threshold or curb with your front wheels.
- Make sure that the footrests do not touch the threshold or curb.
- You can now move the front wheels up the curb.
- Then, move the rear wheels up the curb.


Illustration 21. Guiding the wheelchair without assistance from a caregiver.

### Forward with assistance from a caregiver:

- Make sure that the footrests do not touch the curb.
- The caregiver should tilt the wheelchair back until the front wheels are lifted above the curb.
- The person sitting in the wheelchair should lean backwards to move the centre of gravity over the rear wheels.
- Move the front wheels onto the curb and then move the rear wheels onto the curb.
- Taller curbs can be crossed by driving backwards with assistance from a caregiver.
- Turn the wheelchair in such a way that the rear wheels face the curb.
- The person sitting in the wheelchair should lean backwards to move the centre of gravity over the rear wheels.
- The caregiver should pull the wheelchair onto the curb.
- The person sitting in the wheelchair should return to the standard position in the wheelchair.

Illustration 22. Guiding the wheelchair with assistance from a caregiver.

#### **WARNING!**

Tipping risk.

When going up or down stairs, you may lose your balance, causing the wheelchair to tip over – when going up or down stairs with more than one step, always use the assistance of two caregivers.

### **Moving forward or backwards up or down the stairs:**

- You can move up or down the stairs one step at a time, as shown in illustration opposite. The first caregiver should stand behind the wheelchair, holding the push handles. The second caregiver should hold the solid part of the front frame to stabilise the wheelchair from the front.


Illustration 23. Moving in the wheelchair with the assistance of two caregivers.

### **Stability and balance in the seated position**

#### **Leaning forward.**

Some everyday activities and actions require you to lean forward, backwards or to the side while sitting in the wheelchair. This has a big impact on the stability of the wheelchair. To maintain your balance at all times, follow the instructions shown in illustration opposite.


Illustration 24. Leaning forward.

#### **WARNING!**

Risk of falling out.

If you lean forward, you may fall out of the wheelchair. Never lean too far forward or move to the front of the seat to reach any objects. When picking up things from the floor, do not lean forward between your knees.

#### **Leaning backwards.**

1. Turn the front wheels forward. (To do this, move the wheelchair slightly forward and then reverse).
2. When reaching behind you, do not reach so far that you need to change the seated position.


Illustration 25. Reaching backwards.

**WARNING!**

Risk of falling out.

If you lean too far backwards, the wheelchair may tip over. Do not lean out over the backrest. Use tipping protection.

## Basic troubleshooting

Users can try to solve the problem using the following methods.

If the problem persists, please contact your distributor or the manufacturer.

**WARNING!**

If the battery indicator light or speed indicator light appears flashing on the control panel - this indicates that the wheelchair has detected a fault.

Below is a checklist for self-diagnosis:

Number of Flashes	Fault description	Meaning and actions
	The power indicator emits one short beep followed by seven rapid beeps with intervals.	<ul style="list-style-type: none"><li>• The joystick is not in the central position.</li><li>• First check that the joystick is in good condition/disconnect the red interface/rotate the joystick 360 degrees/connect and switch on the power to the wheelchair to check that everything is OK.</li><li>• If the problem persists, contact your distributor.</li></ul>
	The power light intermittently and continuously emits two beeps	<ul style="list-style-type: none"><li>• Poor wire connection of left motor.</li><li>• Check left motor wire has good contact/is loose/unplugged/disconnected.</li><li>• Swap left and right motor leads - after swapping, 4 alarms indicate left motor failure and 2 alarms indicate controller failure.</li><li>• If the problem persists, contact your distributor.</li></ul>

	<p>The power light illuminates intermittently and emits 3 beeps continuously</p>	<ul style="list-style-type: none"> <li>• Poor left brake connection.</li> <li>• Check left brake lever is in engaged drive mode/is brake cable in good position/loose/unplugged.</li> <li>• Swap left motor cable with right motor cable - 5 alarms after swapping indicates left brake not working and 3 alarms indicate controller failure.</li> <li>• If the problem persists, contact your dealer.</li> </ul>
	<p>The power and speed indicator emits intermittent 3 beeps and light signals.</p>	<ul style="list-style-type: none"> <li>• Brake fault</li> <li>• Check that the left and right brake levers are operating in drive-on mode.</li> <li>• Restore factory settings (electromagnetic brake).</li> <li>• If the problem persists, contact your distributor.</li> </ul>
	<p>Power and speed indicator emits 4 intermittent sound and light signals</p>	<ul style="list-style-type: none"> <li>• Poor right motor cable connection.</li> <li>• Check that the right motor cable has good contact/is not loose/disconnected/disconnected.</li> <li>• Swap left and right motor wires. Four alarms after swapping indicate a failure of the right motor and two alarms indicate a failure of the controller.</li> <li>• If this is the case, the unit will need to be sent back to distributor for repair.</li> </ul>
	<p>Power and speed indicator emits 5 intermittent sound and light signals</p>	<ul style="list-style-type: none"> <li>• Poor connection of the right brake.</li> <li>• Check that the right brake lever is in engaged drive mode / brake cable is well connected/ loose/unplugged/disconnected.</li> <li>• Swap left and right motor cables, 5 alarms after swapping indicate a failure of the right brake and 3 alarms indicate a failure of the controller.</li> <li>• In this case it will be necessary to send the unit back to the distributor centre for repair.</li> </ul>

	<p>Power and speed indicator emits intermittent 6 beeps</p>	<ul style="list-style-type: none"> <li>• Current overload protection.</li> <li>• If there are no six beeps then no attention is required.</li> </ul>
	<p>Power and speed indicator emits intermittent 7 beeps</p>	<ul style="list-style-type: none"> <li>• Incorrect joystick parameter.</li> <li>• Contact your distributor.</li> </ul>
	<p>The first red power indicator light flashes continuously.</p>	<ul style="list-style-type: none"> <li>• Is the battery voltage correct?</li> <li>• The battery needs charging.</li> <li>• Control panel fault</li> <li>• Contact your distributor.</li> </ul>

## Batteries and charging

Before charging the batteries, read the instructions and prepare the charger. The battery status indicator on the control panel indicates the approximate power of the batteries by means of 5 LEDs. All lit LEDs indicate fully charged, and only 1 LED lit indicates an immediate need to charge the batteries.. To determine this exactly, check the battery charge indicator while driving on dry, level ground at top speed.

### Precautions for Li-ion Battery

1. Do not charge the battery when it is full.
2. The battery should not be deeply discharged. The optimum level is 30%~50%, but not over 80%.
3. If the battery is not used often, it should not be charged repeatedly - it should be charged every two months.
4. The battery should be charged after use and not stored with a low charge.
5. The battery must not be near a fire, immersion in water, overheating (45°C), subject to violent vibrations, collisions and short circuits is prohibited.
6. Do not charge the battery for long periods of time, up to a maximum of 10 hours.
7. It is recommended to contact your distributor for the disposal of the battery.

### First charge

There is no need to charge for the first time. Before charging, use the device until the power indicator turns red. When the charger indicator turns green, this indicates that charging is complete. Keep the batteries fully charged so that the wheelchair can move efficiently.

## Connecting the charger

The wheelchair can be charged via the charging socket on the bottom of the control panel. Connect the charger to the XLR socket and then plug into a grounded electrical outlet. When the batteries are fully charged, disconnect the charger first from the socket and then from the wheelchair.


Illustration 26. Connecting the battery charger.

**ATTENTION!**

The button on the battery should be in position "ON" - during charging.

## Operating temperatures

The wheelchair is designed for use, ideally in temperatures between 5°C~41°C. In extremely cold conditions, the batteries may freeze, depending on charge, use and composition. In extremely high temperatures, the wheelchair may run at a lower speed. This is due to a built-in mechanism to protect against damage to the gearbox and other electronic components when operating in hot conditions.

## Battery operation and charging

- Turn the power button to the "OFF" position.
- Connect the charger plug to the charging port, which is located on the handle of the control panel.

**ATTENTION!**

Keep the batteries charged and avoid discharging them completely. To keep the batteries in good condition and to extend their life, we recommend charging a minimum once a month.

## The process of charging batteries outside the wheelchair

1. The wheelchair has a controller charging port and a battery charging port.
2. Charging is done via the charging port on the controller handle (charging method 1). You can also charge by removing the battery box (charging method 2).
3. Connect the charger to a domestic grounded power source (electrical outlet).
4. A red light will come on when charging begins. Charging time depends on the condition and temperature of the battery and is usually between 8 and 10 hours. When charging is complete, the charger will automatically switch off and the charger indicator will turn green.
5. When the charger indicator turns green, this indicates that charging is complete. At this point, first unplug the charger from the electrical outlet and then unplug the charger from the wheelchair.

### WARNING!

- Removal of the ground pins may cause an electrical hazard. If necessary, when accessing a 2-pin socket, use a suitable, approved 3-pin adapter.
- Never use an extension cord to plug in the charger. Plug the charger directly into a properly installed standard grounded electrical outlet.
- It is prohibited, while charging the batteries, for children to play unsupervised near the wheelchair. We recommend that you do not charge the wheelchair batteries if the wheelchair is occupied by user.
- Dangerous gases may be generated when charging the batteries. Keep the wheelchair and charger away from sources that can cause ignition such as fire, sparks and ensure proper ventilation when charging the batteries.
- Charge the wheelchair batteries with the supplied external charger. Do not use any other type of charger.
- Before each use, check the charger, cables and connectors for damage. If you find any damage, contact your supplier.
- Do not open the battery charger case. If the charger does not work correctly, contact your authorized supplier.
- Be aware that during charging, the battery charger case may be hot. Avoid contact with skin and other surfaces that may be damaged by heat.
- Do not leave the charger switched on on the seat of the pushchair. The charger generates heat during the charging process which can damage the upholstery.
- Please check the charger indicator to see whether the battery is full or not. Red color means “charging in progress” and green color means “fully charged”. Please do not determine whether the battery is full or not by using the 4-mode power indicator, but when the charger indicator turns green.

**WARNING!**

If the charger has not been tested and approved for outdoor use then it should not be exposed to extreme weather conditions. If the battery charger is exposed to bad or extreme weather conditions, before using it indoors, leave it unplugged to allow it to adapt to differences in ambient conditions.

## Reporting technical defects and repairing the device

**ATTENTION!**

If questions or concerns arise about the performance of the product, please contact your distributor.

If a defect or damage to the product is discovered during use or maintenance, it must be reported immediately to the distributor or manufacturer.

**WARNING!**

Under no circumstances should an unqualified person repair or alter the product or its design.

The only person allowed to make any repairs or adjustments is an employee authorised by the manufacturer.

## Care and maintenance

Your wheelchair requires minimal care and maintenance. If you are not sure of your care and maintenance capabilities listed below, you can have the inspection and check carried out by an authorized supplier. The components listed below require periodic inspection and maintenance.

### Tires

Regularly check your electric wheelchair's tires for signs of damage or wear.

### Wheel Replacement

The wheelchair is equipped with polyurethane (PU) tires. Replace the entire wheel if the tire has become worn or damaged. For information on replacing wheels, contact your authorized wheelchair supplier.

**WARNING!**

The wheels on the wheelchair must be repaired/replaced by a qualified technician. During this procedure, make sure that, the power is switched off and the wheelchair does not have the clutch lever switched disengaging the drive.

### External surfaces

Bumpers, tires and rims can be occasionally cleaned with rubber and vinyl cleaner.

## Cleaning

Use a damp cloth and a mild, non-abrasive cleaner to clean the plastic and metal parts of the electric wheelchair. Avoid using products that can scratch the wheelchair's surfaces.

## Disinfection

If necessary, use approved disinfectants. Before using, make sure the disinfectant is safe for the wheelchair's surface. Only tested and approved disinfectants may be used for disinfection. Information on recognized and approved disinfectants and disinfection methods is available from your national health authority.

### WARNING!

Follow the safety instructions provided by the manufacturer of cleaning and disinfecting products. Not complying may result in skin irritation or early damage to the upholstery or electric wheelchair.

## Battery terminal connections

Make sure the connections remain tight and non-corroded.

## Electric wiring harnesses

- Regularly check wire connections.
- Regularly check the insulation of the wires, including those from the charger, for wear or damage.
- Have an authorized supplier repair or replace any damaged connectors, connections or insulation before using your electric wheelchair again.

## Daily checks

- With the power off, check the operation of the drive control lever. Make sure that the lever is not bent or damaged and that when released, it returns to the neutral position. Self-repair is prohibited. In case of problems, contact an authorized supplier.
- Visually inspect the control panel. Make sure it is not burnt, cut or wires sticking out. Contact your authorized supplier if there are problems.
- Check for worn (flat) spots on the tires. Worn tires can affect stability.
- Check armrests for loose, strained or damaged parts. If there are problems, contact your authorized supplier.
- Check brakes. Checks should be made on a level surface without obstacles within 1 meter of the wheelchair.

## **Weekly checks**

- Check the controller and charger connections for corrosion. If necessary, contact an authorized supplier.
- Check the pressure if there is a pneumatic tire on the equipment. If the tire does not hold pressure, contact an authorized supplier for repair.

## **Monthly checks**

- Check that the anti-tip wheels do not rub against the ground when operating the electric wheelchair.
- Check the anti-tip wheels for wear and tear. Replace if necessary.
- Check for tire wear. Contact an authorized supplier for repairs.
- Keep the wheelchair clean and away from foreign objects such as mud, sand, hair, food, drinks, etc.

## **Yearly checks**

- Take your electric wheelchair to an authorized supplier for an yearly inspection, especially if you use your wheelchair daily. This will ensure the correct operation of the wheelchair and help avoid future breakdowns.

## **Plastic covers**

If your electric wheelchair has a shiny chassis surface then it has been treated with a clear fixative coating. You can apply a small coat of car wax to help restore the glossy surface. If your wheelchair's chassis has a matte finish, use only products designed for matte finish surfaces. Do not use wax or other products designed for glossy paintwork.

## **Console, Charger, and Rear Electronics**

Keep these areas free of moisture. Allow these areas to dry thoroughly if they have been exposed to moisture before operating your electric wheelchair again.

## **Storage of the wheelchair**

If you plan not to use the wheelchair for a period of time then follow these rules:

- Fully charge batteries before storage.
- Switch the wheelchair battery off.
- Store the wheelchair in a warm, dry room.
- Avoid storing in a place where the wheelchair will be exposed to temperature changes.
- Recommended storage temperatures: from -25°C do +50°C.
- Recommended ambient humidity for storage: ≤85%.

Batteries that are regularly completely discharged, not regularly charged, stored at inappropriate temperatures, or stored without being charged may be permanently damaged resulting in unreliable operation and shortened battery life. To ensure proper operation, it is recommended that batteries be charged regularly, throughout the period of prolonged storage of the wheelchair.

## **Re-use (reusable product)**

The product is suitable for reuse. The number of times it can be reused depends on how often and in which way the product is used. Before reuse, thoroughly clean and disinfect the product and check its safety features. If necessary, take repair action. If necessary, take repair action.

## **Product life cycle**

The expected life of a medical device is 5 years, as long as it is used as intended and all maintenance and check is performed. The lifespan of a medical device also depends on the frequency of use, the environment in which it is used and its care. By using spare parts, the life of the product can be extended. Spare parts are generally still available for 5 years after the model is discontinued. This expected lifetime of a medical device does not imply an additional guarantee.

## **Replacement parts**

The availability of replacement components and/or spare parts is dependent on the capabilities of the supplier, not the manufacturer. For more information on replacement components and/or spare parts, please contact your authorized supplier.

## **Disposal of an electric wheelchair**

The electric wheelchair must be disposed of in accordance with applicable local and national regulations. For information on the proper disposal of packaging, metal frame parts, plastic parts, electronics, batteries, neoprene, silicone and polyurethane materials, please contact your local waste disposal company or authorized supplier.

## Warranty

Details of the warranty can be found in the attached “warranty card” document.

### ATTENTION!

The warranty does not cover deterioration of battery performance caused by being left in a fully discharged state, in cold conditions for an extended period of time, or by excessive use.


## Technical parameters

Product name	TWIST Power Wheelchair
Seat width	460 mm
Maximum user weight (weight capacity)	125 kg
Maximum load capacity	125 kg
Maximum speed (forward)	6 km/h
Maximum speed (backward)	2.5 km/h
Seat surface height from floor	495 mm
Seat width between sides	440 mm
Seat depth	430 mm
Backrest height	420 mm
Overall length of the wheelchair	1080 mm
Overall height of the wheelchair	865 mm
Overall width of the wheelchair	625 mm
Folded Width	990x604x380 mm
Shank length (min/max)	295 mm
Armrest height from the ground	715 mm
Thickness of the seat cushion	40 mm
Footrest	243x298 mm
Total length of wheelchair with/without footrests	1080/785 mm
Weight of wheelchair with/ without batteries	24.2 kg / 21.7 kg
Height of the wheelchair for transport	380 mm
Front wheels	8" PU
Rear wheels	12" PU (electric drive) 24" PU (manual)
Maximum obstacle to overcome safely	40 mm
Maximum safe inclination (degree)	9° (15,83%)
Turning radius (degree)	780 mm
Energy range* (travel range)	13-15 km

Battery	24V 13Ah Li-Ion
Motor power	150 W (brushless motor)
Armrest height from the ground	Flip-Up
Controller type	Intelligent controller brushless
Minimum braking distance at max. speed	1,5 m
Battery charger parameters	4V 2A
Suspension	tak
Safety	Seat belt
Frame material	Aluminum
Storage and use temperature	from -5°C to +40°C
Humidity to store and use	≤80%
<p>We reserve the right to introduce technical changes. The tolerance in the measurements is <math>\pm 15</math> mm /1.5 kg / 1.5°. The theoretical turning radius decreases if the wheelchair is often used on inclines, uneven surfaces or for driving up curbs. The maximum distance covered has been tested under ideal conditions in accordance with ISO7176-4.</p> <p>Speed and range depends on the user's weight, type of terrain, charge and battery condition. Information is current at the time of publication.</p>	

## Symbols

The symbols used below on the wheelchair are used to identify warnings, mandatory actions and prohibited actions. Be sure to read everything and understand their meaning.

	Medical device
	Date of manufacture
	Manufacturer
	Authorized representative in the European Union
	Conformité Européene (compliance with the European Community)
	Catalogue number
	Serial number
	Lot number
	Unique Device Identifier
	Model number
	Follow the instructions for use
	Refer to the instructions for use
	Warning
	Protect from humidity

	Keep away from sunlight
	Possibility of jamming fingers
	Permitted temperature
	For indoor/outdoor use
	Maximum user weight
	Maximum permissible total weight
	Maximum safe inclination
	The product is not intended to be used as a seat for transport in a motor vehicle
	Do not stand on the footrest
	N = Neutral (Brakes disengaged; lever disengaging the drive turned on)
	D = Drive (Brakes engaged; unit in Drive Mode)

**SELECTIVE COLLECTION OF ELECTRICAL AND ELECTRONIC EQUIPMENT.** This symbol on the products or accompanying documents indicates that used electrical or electronic products and batteries must not be disposed of in municipal waste. For proper disposal, restoration or recycling, such products should be returned to collection points for this type of waste. Alternatively, in some countries of the European Union or other European countries, you can return your products to a local seller at the time of purchase of a similar new product. By correctly disposing of this product, you will help preserve valuable natural resources and support the prevention of potential negative impacts on the environment and human health, which could result from incorrect waste disposal. More information is available from the communal offices or waste collection sites. In case of incorrect disposal of this product, penalties may be imposed in accordance with local regulations. If you want to dispose of an electrical or electronic device, get the information you need from your dealer or supplier.

**Disposal in countries outside the European Union:** This symbol is valid in the European Union. If you want to dispose of this device, please obtain the necessary information on the correct disposal at your local authorities or from your dealer.

## Contact details

### Manufacturer

mdh Sp. z o.o.

Address: ul. Maratońska 104, 94-007 Łódź, Poland

tel. +48 42 674 83 84, fax. +48 42 636 52 21

[www.mdh.pl](http://www.mdh.pl)

[www.viteacare.com](http://www.viteacare.com)

We reserve the right to make technical changes due to technical progress.