

ANTI BEDSORE MATTRESS

*alternating pressure mattress
of tubular structure*

catalogue number **VCM502**

PRODUCT CHARACTERISTICS:

- alternating pressure, anti-bedsore mattress, with the highest number of chambers, made of PVC
- smooth manual pressure adjustment, individually adapted according to the patient's weight
- operation cycle: 10 minutes
- quiet electric pump
- can be placed on standard bed mattresses
- 1 extra chamber
- cover
- color: dark blue

USE:

Bedsore prevention and treatment.



 **VITEA[®]
CARE**

Length (cm)	190
Width (cm)	85
Height (cm)	11
Maximum capacity (kg)	140
Weight (kg)	3,4
Weight of pressure pump (kg)	1,3