

TUBE MASTER

anti bedsore mattress

catalogue number **VCM502D**

PRODUCT CHARACTERISTICS:

- pneumatic, alternating pressure, anti-bedsore mattress
- tube structure with 18 chambers in two rows arrangement alternately filled with air, with a cover
- made of nylon, to facilitate cleaning
- smooth, manual pressure adjustment, individually adapted according to the patient's weight
- thanks to the mattress structure based on separate air chambers, patient's weight is evenly distributed, the mattress provides correct ventilation (oxygenation and drainage of moisture), great insulation properties guarantee temperature comfort
- number of chambers: 18
- operation cycle: 12 minutes
- modern, quiet electric pump
- can be placed on standard bed mattresses
- color: dark blue

USE:

Bedsore prevention and treatment.



Length (cm)	200	Work cycle	12 minut
Width (cm)	90	Amount of air pumped	6-8 l/min
Height (cm)	11,5	Weight (kg)	4
Maximum capacity (kg)	145	Weight of pressure pump (kg)	1,2
Power supply	220-240 V/50-60 HZ		