

FELXI LIGHT

aluminium
wheelchair



catalogue number: VCWK9AFL



full version

PRODUCT CHARACTERISTICS:

- foldable, lightweight aluminum cross frame
- upholstery made of sponge trimmed with black nylon
- side panel with wheel arch
- removable and swing away footrests
- shank length adjustment
- seat height adjustment
- backrest height adjustment
- solid front wheels, size 6"
- pneumatic rear wheels equipped with a quick release system, size 24x1"
- possibility of shifting the center of gravity by adjusting the rear wheels
- aluminum handrim
- pressure brakes, foldable brake lever
- frame color: blue or black
- available version with additional equipment: adjustable seat depth, cushion, anti-tip wheels, armrests

INTENDED USE:

The wheelchair is used for the individual transport of disabled persons after a spinal cord injury or persons with cerebral palsy and various injuries to the lower limbs, neuromuscular diseases, craniocerebral injuries and geriatric patients. It is also intended for people who are looking for a comfortable way to actively participate in everyday life.

Size	15"	16"	17"	18"
Seat width (cm)	38	40	43	46
Seat depth (cm)	42	42	42	42
Wheelchair height (cm)	81 - 97	81 - 97	81 - 97	81 - 97
Wheelchair width (cm)	55	58	60,5	63
Wheelchair length (cm)	93	93	93	93
Width of folded wheelchair (cm)	30	30	30	30
Wheelchair weight (kg)	~12	~12	~12	~12
Maximum capacity (kg)	114	114	114	114