

CROSS SMART

W459 Power Wheelchair



catalogue number: **DRVKOP** model: **CROSS SMART**



PRODUCT CHARACTERISTICS:

- foldable frame of cross structure (easy to transport)
- foldable backrest
- seat width: 430 mm
- shock absorbing system front wheels
- 50 mm thick foam cushion trimmrd with soft, breathable fabric
- lifted armrest
- control unit mounted on the right or left side upon request
- removable footrest with lower leg support (adjustable tilt angle)
- removable and tiltable footrest
- footrest adjustable length
- anti-fall belt
- 8" full front wheels, 12" pneumatic rear wheels
- versions available:
 - standard
 - remote controller
 - caregiver control unit

INTENDED USE:

A powered wheelchair electric is a quality of life device designed to ensure mobility.

Maximum speed* (km/h)	6	Front wheels	8" PU
Seat width (cm)	43	Rear wheels	12" pneumatic
Height of backrest (cm)	44	Range* (km)	20
Maximum wheelchair height with backrest (cm)	94	Maximum capacity (kg)	120
Seat height from ground (cm)	54	Batteries (Ah)	20 Ah (24V)
Wheelchair width (cm)	66	Motor (W)	2x250 W
Wheelchair length (mm)	110	Other features	electromagnetic brake

* Speed, maximum distance of the vehicle depends on the user weight, tire pressure, battery load, battery capacity and the terrain on which the vehicle is used. The estimated range refers to laboratory conditions.