

CLEVER TRIP

power wheelchair



catalogue number: **DRVIOP**



PRODUCT CHARACTERISTICS:

- modular design – the frame can be dismantled for easy transport
- modern design
- detachable battery pack with off-board recharging capability
- swivel seat, horizontally adjustable
- backrest angle adjustment
- smoothly adjustable steering wheel angle
- set of rear-view mirrors
- shopping basket
- headrest
- full road lighting
- electromagnetic brakes
- coupling for disengaging the electric drive
- front and rear bumper
- front wheels: 10"
- rear wheels: 10"
- shock-absorbing frame

INTENDED USE:

The product is intended for disabled people with musculoskeletal conditions to increase their mobility.

MANUFACTURER:

mdh sp. z o.o. , ul. Maratońska 104, 94-107 Łódź, tel. +48 42 676 32 54,
mdh.rehabilitacja@meyragroup.com

Max speed (km/h)	10	Weight of wheelchair with/without batteries (kg)	90/66
Seat width (cm)	45	Height of wheelchair for transport (cm)	64
Seat width between armrests (cm)	49-66	Front wheel size	10"
Seat depth (cm)	45,5	Rear wheel size	10"
Backrest height (cm)	50	Ground clearance (cm)	11
Total height of the wheelchair (cm)	112	Turning radius (cm)	140
Height of armrest from seat upholstery (cm)	19	Range* (km)	37
Height of the upper surface of the seat from the ground (cm)	69	Max user weight (kg)	150
Total width of the wheelchair (cm)	59	Batteries (Ah)	2 x 45
Total length of the wheelchair (cm)	127	Motor (W)	320

* Speed, maximum distance of the vehicle depends on the user weight, tire pressure, battery load, battery capacity and the terrain on which the vehicle is used. The estimated range refers to laboratory conditions.