

ORANGE

STAR power wheelchair



catalogue number: **DRVION**



PRODUCT CHARACTERISTICS:

- seat width: 16", 18"
- shock-absorbing frame
- adjustable seat depth and width (430-480 mm)
- electrically adjustable seat angle 0-20°
- Electrically adjustable backrest angle 90°-140°
- safety belt
- soft and breathable upholstery
- adjustable height and width of armrests
- removable footrests with calf support (Flexion angle adjustment)
- adjustable footrest length
- full illumination
- disconnection of electric drive
- control unit can be switched from right to left
- horizontal adjustment and tilting of the control unit
- anti-tip wheels
- aluminium rims
- adjustable headrest
- handles for the caregiver

INTENDED USE:

The product is intended for disabled people with musculoskeletal conditions to increase their mobility. Designed for people with paralysis and paresis, amputees, people with stiff lower limbs and those who periodically need to use a wheelchair.

MANUFACTURER: mdh sp. z o.o., ul. Maratońska 104,
94-107 Łódź, tel. +48 42 674 83 84, mdh.rehabilitacja@meyragroup.com

Max speed (km/h)	10	Weight of wheelchair with/without batteries (kg)	106,5/76,5
Seat width (cm)	43-48	Height of wheelchair for transport (cm)	85
Seat depth (cm)	43-48	Front wheel size	9"
Backrest height (cm)	52	Rear wheel size	14"
Total height of the wheelchair (cm)	121,5	Max obstacle to be overcome safely (cm)	4,5
Height of armrest from seat (cm)	17-29	Turning radius (cm)	90
Seat cushion thickness (cm)	7,5	Range* (km)	30
Height of seat with cushion from ground (cm)	54	Max User weight (kg)	up to 136
Total width of the wheelchair (cm)	58	Batteries (Ah)	52 (Gel)
Overall length of wheelchair with footrests (cm)	182	Motor (W)	2 x 340

* Speed, maximum distance of the vehicle depends on the user weight, tire pressure, battery load, battery capacity and the terrain on which the vehicle is used. The estimated range refers to laboratory conditions.