

# MULTIPLUS

## Power wheelchair



catalogue no: **DRVIOM**



### PRODUCT CHARACTERISTICS::

- combines several functions: electric wheelchair, electric rollator, manual rollator, transport rollator
- aluminium frame with cross-brace
- two control modules (with display) depending on the function used
- soft seat
- swing away footrests equipped with heel pads
- possibility of step adjustment of shank length
- lifted armrests
- height adjustable handles
- back support while seated (with the possibility of flipping)
- safety belt
- manual and electromagnetic brakes
- 4 wheels with a diameter of 8" (swivel front wheels)
- shopping bag
- two colors available: champagne, graphite

### INTENDED USE:

The product is intended for disabled people with motor conditions to increase their mobility, when: • limited or no muscle strength • limited range of movement • limited trunk and body stability • other neurological disorders requiring support

### CONTRAINDICATIONS:

Contraindications to upright standing. Injuries of upper limbs which make it impossible supporting yourself in a wheelchair.

MANUFACTURER: Maratońska 104, 94-007 Łódź, Poland, tel. +48 42 674 83 84

Maximum speed (km/h)	6	Front wheel	8"
Seat width (mm)	450	Front wheel	8"
Effective seat depth (mm)	300	Distance range* (km)	do 15
Seat height from the ground (mm)	515	Maximum user weight (kg)	110
Wheelchair total height (mm)	785 - 868	Maximum load of bag under the seat (kg)	5
Wheelchair total width (mm)	650	Wheelchair weight (kg)	21
Wheelchair total length with legrests (mm)	1050	Batteries (Ah)	10
Dimensions when folded (LxWxH) (mm)	717x320x785	Motor (W)	2 x 250

\* Speed, maximum distance of the vehicle depends on the user weight, tire pressure, battery load, battery capacity and the terrain on which the vehicle is used. The estimated range refers to laboratory conditions.