

MINI III

power wheelchair



catalogue number: **DRVI00**



PRODUCT CHARACTERISTICS:

- modular design – the frame can be dismantled for easy transport
- easy operation, intuitive control panel
- „delta” steering wheel
- easy folding of the wheelchair for transport without the use of tools
- removable battery pack with the option of charging them outside the wheelchair
- modern design
- swivel seat, with folding backrest
- adjustable steering wheel angle
- set of rear-view mirrors
- shopping basket
- LED road lighting
- electromagnetic brakes
- a clutch disconnecting electric drive
- front and rear bumper
- front wheel: 8”
- rear wheel: 8”
- shock-absorbing frame

INTENDED USE:

The product is intended for disabled people with musculoskeletal conditions to increase their mobility.

CONTRAINDICATIONS:

Contraindications to upright standing. Injury or paresis of the upper limbs preventing the use of a power wheelchair.

MANUFACTURER:

mdh sp. z o.o. , ul. Maratońska 104, 94-107 Łódź, tel. +48 42 676 32 54,
 mdh.rehabilitacja@meyragroup.com

Max speed (km/h)	8	Weight of wheelchair with/without batteries (kg)	53/44
Seat width (cm)	43,5	Height of wheelchair for transport (cm)	35
Seat width between armrests (cm)	42-62	Front wheels	8”
Seat depth (cm)	41	Rear wheels	8”
Backrest height (cm)	37	Clearance (cm)	7
Overall height (cm)	92	Turning radius (cm)	130
Height of armrest from seat upholstery (cm)	20-26	Range* (km)	15
Seat height (cm)	57	Maximum user weight (kg)	136
Overall width (cm)	54-70	Batteries (Ah)	15
Overall length (cm)	113	Motor (W)	270

* Speed, maximum distance of the vehicle depends on the user weight, tire pressure, battery load, battery capacity and the terrain on which the vehicle is used. The estimated range refers to laboratory conditions.