

ROSE 4 - wheel rollator

catalogue number **VCBK222T with bag or basket**

PRODUCT CHARACTERISTICS:

- aluminum folding frame (easy to transport)
- soft seat covered with eco-leather
- lockable brakes
- four maintenance-free 20 cm diameter wheels,
- swivel front wheels
- adjustable height
- ergonomic handles
- shopping basket
- possibility of changing the direction of moving while the rollator is pushed by third person
- footrest

USE:

Walking aid which improves stability and helps to transport the luggage.



handles with brake levers



handles that helps to stand up (bag version)



Height (cm)



Weight (kg)



Maximum capacity (kg)



Width (cm)



Wheel diameter (cm)

	Height (cm)	Weight (kg)	Maximum capacity (kg)	Width (cm)	Wheel diameter (cm)
basket	93	8	110	62	20
bag	81-96	8	120	62	20