

WAVE

walker H4



catalogue no: **VCBPO032**



PRODUCT CHARACTERISTICS:

- foldable aluminum frame
- reinforced joints
- construction ensures easy standing up
- legs finished with non-slip rubber tips
- ergonomic handles
- gradual height adjustments every 2.5 cm

INTENDED USE:

Walking aid designed for disabled people for those with locomotor disorders. Ideal for rehabilitation process and physiotherapy (verticalisation, lower limb relief).

CONTRAINDICATIONS:

Contraindications to the upright position, injuries and dysfunctions of upper limbs preventing the patient from supporting themselves.

Width (cm)	62
Height (cm)	81,5 - 91,5
Weight (kg)	2,5
Maximum capacity (kg)	100