

CARTER

wheeler
walker



catalogue no: **VCBPO032**



PRODUCT CHARACTERISTICS:

- foldable, aluminium frame
- reinforced joints
- ergonomic handles
- four wheels with a diameter of 8 cm
- lockable rear wheels
- gradual height adjustments every 2.5 cm

INTENDED USE:

Walking aid designed for disabled people for those with locomotor disorders. Ideal for rehabilitation process and physiotherapy (verticalisation, lower limb relief).

CONTRAINDICATIONS:

Contraindications to the upright position, injuries and dysfunctions of upper limbs preventing the patient from supporting themselves.

Width (cm)	55
Height (cm)	81,5 - 91,5
Weight (kg)	3,7
Maximum capacity (kg)	110