

AMBONA

4-wheel walker



catalogue no: **VCBK2A_01**



PRODUCT CHARACTERISTICS:

- steel frame
- soft armrest height adjustment
- handle height and angle adjustment
- rear wheels with brakes
- 3" rubber wheels upon request
- silver with rounded triangular adjustment handles

INTENDED USE:

Rehabilitation supports are an ancillary equipment to facilitate tilting and moving of patients with different types of motor disorders. They can be used as an aid for patients after surgery (lower limbs) who are unable to use crutches, showing no coordination disorders, patients with mild dysfunction of the lower limbs or weakness. They are used in patients with in patients with impaired balance preventing them from independent vertical position and locomotion, as well as various types of circulatory dysfunction (e.g. intermittent claudication).

Supports can also be used as an aid for geriatric patients.

CONTRAINDICATIONS:

Contraindications to the upright position, injuries and dysfunctions of upper limbs preventing the patient from supporting themselves, impaired coordination and perception.

MANUFACTURER: Maratońska 104, 94-007 Łódź, Poland, tel. +48 42 674 83 84

VCBK2A_01

Max width of the walker base (mm)	53
Inner width of the walker base (mm)	46
Max length with brakes (max extended handles) (mm)	75 - 97
Max height with handles (mm)	97,5 - 137
Cushion measurements W x L x H (mm)	58 x 49 x 5
Max width of the walker (mm)	58
Cushion cut-out measurements (mm)	36 x 28
Handles (mm)	11
Handles spacing (mm)	29
Walker weight (mm)	11
Maximum capacity (kg)	125