

DUAL

2-functional walker



catalogue no: **VCBPOO31**



PRODUCT CHARACTERISTICS:

- foldable aluminum frame
- stable and durable construction
- supports both standing and walking
- walking mode can be easily blocked
- legs finished with non-slip rubber tips
- ergonomic handles
- gradual height adjustments every 2.5 cm

INTENDED USE:

Walking aid designed for disabled people for those with locomotor disorders. Ideal for rehabilitation process and physiotherapy (verticalisation, lower limb relief).

CONTRAINDICATIONS:

Contraindications to the upright position, injuries and dysfunctions of upper limbs preventing the patient from supporting themselves.

Width (cm)	54
Height (cm)	82 - 99
Weight (kg)	2,5
Maximum capacity (kg)	110