

BILLY

3-wheel rollator



catalogue no: **VCBK5221**



PRODUCT CHARACTERISTICS:

- foldable, aluminum frame (easy to transport)
- lockable brakes
- three 19 cm diameter wheels
- adjustable height
- adjustable rear wheels position
- ergonomic handles
- zipped shopping bag
- front shopping basket with tray

INTENDED USE:

Rehabilitation supports are an ancillary equipment to facilitate tilting and moving of patients with different types of motor disorders. They can be used as an aid for patients after surgery (lower limbs) who are unable to use crutches, showing no coordination disorders, patients with mild dysfunction of the lower limbs or weakness. They are used in patients with impaired balance preventing them from independent vertical position and locomotion, as well as various types of circulatory dysfunction (e.g. intermittent claudication). Supports can also be used as an aid for geriatric patients.

CONTRAINDICATIONS:

Contraindications to the upright position, injuries and dysfunctions of upper limbs preventing the patient from supporting themselves, impaired coordination and perception.

MANUFACTURER: mdh sp. z o.o. , ul. Maratońska 104,
94-107 Łódź, tel. +48 42 674 83 84, mdh.rehabilitacja@meyragroup.com

Width (cm)	68
Height (cm)	80 - 93
Weight (kg)	5,3
Maximum capacity (kg)	120