

TRIPLA

3-functional walker



catalogue no: **VCBPO031**



PRODUCT CHARACTERISTICS:

- foldable aluminum frame
- stable and durable construction
- can be used for walking, standing or as a wheeled walker
- walking mode can be easily blocked
- additional removable front wheels with 12 cm diameter
- ergonomic handles
- gradual height adjustments every 2.5 cm

INTENDED USE:

Walking aid designed for disabled people for those with locomotor disorders. Ideal for rehabilitation process and physiotherapy (verticalisation, lower limb relief).

CONTRAINDICATIONS:

Contraindications to the upright position, injuries and dysfunctions of upper limbs preventing the patient from supporting themselves.

Width (cm)	60
Height (cm)	82 - 100
Weight (with wheels) (kg)	2,5
Maximum capacity (kg)	110