

UMBRELLA

children rehabilitative stroller

Catalogue number: **DRVG0C**

PRODUCT CHARACTERISTIC:

- foldable, aluminum frame with reinforced steel elements
- upholstery made of high quality fabric attested by certificates of quality (OEKO-Tex, SGS)
- smoothly adjustable backrest tilt angle
- front wheels and rear wheels made of PU
- lockable front wheels
- front wheels attached via quick release system with independent brakes
- abduction wedge, front safety barrier covered with fabric and 5-point belt
- side supports and lumbar roller
- plastic footrest platform with gradual height adjustment (every 2.5 cm)
- footrest adjustable tilt angle with feet stabilizing straps
- headrest with height adjustment
- parent window in the canopy
- swivel handles with adjustable height
- standard equipment: canopy, shopping basket, cup holder, rain cover, sleeping cover
- optional equipment: umbrella, tray
- available in wide range of colors
- available in 3 sizes: 1, 2 and 3

USE:

Intended for children with Cerebral Palsy, Spina Bifida, muscle disorders and lower limb dysfunction.



Black and Pink



Black and Green



Black and Red



Black with Color dots



Black and White checkered



Black and Orange