



CE VITEA<sup>®</sup>  
CARE

	size 1	size 2	size 3
Seat width(cm)	30	29 - 34 cm	40
Seat depth(cm)	25	30	40
Backrest height (cm)	50	50	75
Max stroller width (cm)	55	58	63
Max stroller length (cm)	88	98	105
Max stroller height (cm)	114	116	126
Dimensions when folded (lengthx width x height)	122 x 37 x 44	118 x 47 x 44	128 x 44 x 50
Front wheel diameter (cm)	18	18	17
Rear wheel diameter (cm)	22	22	25
Stroller weight (kg)	13,5	15	19
Max child weight (kg)	30	45	75
Max child height (cm)	120	140	170